MDCStoreUtils API Reference Guide

MDCStore™ Database 2.0 or later

Microsoft SQL Server and Oracle

5000957 C August 2012



This document is provided to customers who have purchased Molecular Devices, LLC ("Molecular Devices") equipment, software, reagents, and consumables to use in the operation of such Molecular Devices equipment, software, reagents, and consumables. This document is copyright protected and any reproduction of this document, in whole or any part, is strictly prohibited, except as Molecular Devices may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such usage is intended only to designate those manufacturers' products as supplied by Molecular Devices for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

Molecular Devices makes no warranties or representations as to the fitness of this equipment for any particular purpose and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use to which the purchaser may put the equipment described herein, or for any adverse circumstances arising therefrom.

CE



For research use only. Not for use in diagnostic procedures.

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

The trademarks mentioned herein are the property of Molecular Devices, LLC or their respective owners. These trademarks may not be used in any type of promotion or advertising without the prior written permission of Molecular Devices, LLC.

Patents: http://www.moleculardevices.com/productpatents

Product manufactured by Molecular Devices, LLC.

1311 Orleans Drive, Sunnyvale, California, United States of America 94089.

Molecular Devices, LLC is ISO 9001 registered.

© 2012 Molecular Devices, LLC.

All rights reserved.

Contents

Preface
Who This Manual Is For
About this Manual
Conventions
Chapter 1: Introduction to the MDCStoreUtils API 17
What is the MDCStoreUtils API?
Main Functionality
How to Use the MDCStoreUtils API
Supported database servers
Compiling and linking with MDCStoreUtils API 18
Definitions
Functions Available in MDCStoreUtils API
Chapter 2: Database Connection Functions
MDCS_CONNECTION_GetAttributes
MDCS CONNECTION GetDetails
MDCS_CONNECTION_TestCredentials
MDCS_CONNECTION_CreateConnectionString
MDCS_CONNECTION_GetDBHandleFromString 24
MDCS_CONNECTION_GetDBHandle
MDCS_CONNECTION_GetNewDBHandle25
MDCS_CONNECTION_DestroyDBHandle
MDCS_CONNECTION_Disconnect
MDCS_CONNECTION_Reconnect
MDCS_CONNECTION_GetDetailsFromString 27
MDCS_CONNECTION_BeginTransaction
MDCS_CONNECTION_CommitTransaction 28
MDCS_CONNECTION_RollbackTransaction
MDCS_CONNECTION_GetDatabaseType
MDCS_CONNECTION_GetInfo
MDCS_CONNECTION_SetAsyncMode
MDCS_CONNECTION_CancelQueryExecution(
MDCS_CONNECTION_SetSilentMode
MDCS_CONNECTION_CheckIfDead
Chapter 3: Results Set Functions
MDCS_ImportDataToDB33
MDCS_GetResultsetData
MDCS_GetAllResultsets
MDCS_GetResultsetInfo35

Chapter 4: Result Set Data Types Functions	37
MDCS_GetResultsetDataTypes	37
MDCS_GetAllDataTypes	38
MDCS_GetUniqueDataTypes	38
MDCS_GetAllDataTypesOfAssays	
MDCS_GetDataTypesNotInAssaysByDBName	39
Chapter 5: Error Handling Functions	41
MDCS GetLastError	41
MDCS_GetLastErrorMsg	41
MDCS_GetDefaultErrorMsg	
Chapter 6: BLOB Functions	43
MDCS_BLOB_GetInfo	44
MDCS_BLOB_GetLocationID	
MDCS_BLOB_GetLocation	
MDCS BLOB Get	
MDCS_BLOB_Save	46
MDCS_BLOB_RemoveBlobData	47
MDCS_GetNewDatabaseID	48
MDCS_BLOB_UpdateBlobDescAndName	48
MDCS_BLOB_UpdateBlobDescAndNameEx	49
MDCS_BLOB_SaveBlobEx	
MDCS_BLOB_GetBLOBInfoEx	
MDCS_BLOB_GetInfoByReferenceID	
MDCS_BLOB_GetBlobInfoByReferenceIDEx	
MDCS_BLOB_GetNumBlobOfObject	
MDCS_BLOB_UpdatePlateObjectImageID	
MDCS_BLOB_GetUniqueLocationNamebyFilter	
MDCS_BLOB_GetAllBlobInfoByLocationFilterAndLocationInfo	
MDCS_BLOB_GetAllBlobInfoByPlateID	
MDCS_BLOB_DeleteUnUsedLocations	
MDCS_BLOB_DeleteFileLocation	
MDCS_BLOB_LocationIsUsed	
MDCS_BLOB_RemoveBlobData	
MDCS_BLOB_RemoveFile	
MDCS_BLOB_Attacti	
MDCS_BLOB_OpdateBioDiffice	
MDCS_BLOB_MapLocations	
MDCS_BLOB_cariwfiteroLocation	
Chapter 7: Measurement Sets and Results Functions	
MDCS_ASSAY_AppendMeasurementSet	
MDCS_ASSAY_ImportMeasurementSet	
MDCS_ASSAY_ImportMeasurementSet	
MDCS_ASSAY_InsertDValue	68

		_InsertNumericValue	
MDCS_	_ASSAY_	_InsertSValue	69
MDCS_	_ASSAY_	_InsertStringValue	70
MDCS_	_ASSAY_	_InsertShapeLines	70
		_UpdateShapeLines	
MDCS	ASSAY	 GetShapeLineBlobDesc	72
		 _GetShapeLineBlobDescBySite	
		 GetShapeLineBlobDescBySiteAndSeries	
MDCS	ASSAY		75
		 _GetShapeLineBlobDescByPlate	
		GetMeasurementAttributes	
MDCS	ASSAY		79
		_ _GetMeasurementsBySiteID	
MDCS	ASSAY		82
MDCS	ASSAY	Create	83
MDCS	ASSAY	_GetAssaySiblingFolders	84
MDCS_	ASSAY_	_AddMeasurement	85
MDCS_	_ASSAY_	_GetMeasurementByID	86
MDCS_	_ASSAY_	_GetMeasurementByName	87
MDCS_	_ASSAY_	_Delete	87
MDCS_	_ASSAY_	_DeleteAllForPlate	88
		_DeleteDataForPlate	
		_DeleteFolder	
		_CreateOutlinesTable	
		_ManageFolderSecurity	
		_CanModifyFolder	
MDCS_	_ASSAY_	_CreateFolder	91
MDCS_	_ASSAY_	_ModifyFolder	91
_		_DoesSubFolderExist	
		_GetAllInFolder	
		_GetFolderPath	
		_Reindex	
		_ReindexShapes	
		_OptimizeAll	
		_Copy	
		_Merge	
		_GetSiteCount	
MDCS_	_ASSAY_	GetSiteInfoImageByID	96
		_GetSiteInfoImageBySiteAndSeriesID	
		_MeasurementGetSiteCount	
MDCS_	_ASSAY_	_DeleteMeasurement	97

MDCS_	_ASSAY_	_CellOutlinesGetSiteCount	. 98
MDCS_	_ASSAY_	_CreateNewName	. 98
MDCS_	_ASSAY_	_CreateRun	. 99
MDCS_	_ASSAY_	_GetLatestAssayRunID	. 99
		_GetByID	
MDCS_	_ASSAY_	_CreateProfile	100
MDCS_	_ASSAY_	_AssociateWithPlate	101
MDCS_	_ASSAY_	_GetSpotID	102
MDCS_	_ASSAY_	_GetAllForPlate	102
MDCS_	_ASSAY_	_GetProfiles	103
MDCS_	_ASSAY_	_GetProfile	104
MDCS_	_ASSAY_	_GetProfileInfo	104
MDCS_	_ASSAY_	_UpdateProfile	105
MDCS_	_ASSAY_	_DeleteProfile	105
MDCS_	_ASSAY_	_GetRecord	106
MDCS_	_ASSAY_	_GetAssayByRunID.	106
MDCS_	_ASSAY_	_UsedInDatasets	107
MDCS_	_ASSAY_	_HavePermissionsToModify	108
		_GetImageSourceRecords	
		_GetStatisticalValuesForPlateAndAssay	
		_GetStatisticalValuesForPlateAndMeasurement	
		_GetValuesForPlateAndAssay	
_		_GetValuesForPlate	
		_GetAllSiteMeasurements	
		_GetMeasurementInfoByAssayAndPlate	
		_GetAssayInfoByNamePlateAndSettings	
		_GetAssayIDsOfPlate	
		_GetAssaySettingsOfPlate	
_		_GetMeasurementStatistic	
		_GetShapeLinesBySite	
		_GetShapeLinesBySiteAndSeries	
		_GetMarkedAssaysWithCallback	
		_Restore	
		_CreateForPlate	
		_GetMeasurementByName	
		_GetMeasurementByDBName	
		_GetScopeAttributeByName	
		_GetScopeAttributeByID	
		_DeleteMeasurementInAssays	
		_UpdateDataType	
		_UpdateDataTypeByAssay	
MDCS_	_ASSAY_	_CanModifyAssay	123
		_CalculateStatisticResults	
		_CalculateStatisticEx	
		_GetZPrime	
		_GetZPrimeScopeAttribute	
MDCS_	_ASSAY_	_GetUniqueMeasurementValues	175

MDCS_ASSAY_UpdateMeasurementData 12' MDCS_ASSAY_CreateAttribute 13 MDCS_ASSAY_GetAttributeInfoByDisplayName 13 MDCS_ASSAY_GetAttributeInfoByDBName 13 MDCS_ASSAY_GetAttributeValueByDBName 13 MDCS_ASSAY_GetAttributeValueByDisplayName 13 MDCS_ASSAY_GetAttributeValueByDisplayName 13 MDCS_ASSAY_AssignAttributeValueString 13 MDCS_ASSAY_AssignAttributeValueIng 13 MDCS_ASSAY_AssignAttributeValueFloat 13 MDCS_ASSAY_AssignAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_UpdateMeasurementSetName 13 MDCS_ASSAY_UpdateMeasurementSetName 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_GetUniqueMeasurementValuesByPlate 13 Chapter 8: Common Database Functions 13 MDCS_DATABASE_Optimize 14 MDCS_DATABASE_Compact 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetVersion 14 <	MDCS_ASSAY_FindMeasurementValues	
MDCS_ASSAY_CetattributeInfoByDisplayName. 13 MDCS_ASSAY_GetAttributeInfoByDBName 13 MDCS_ASSAY_GetAttributeValueByDBName. 13 MDCS_ASSAY_GetAttributeValueByDBName. 13 MDCS_ASSAY_GetAttributeValueByDisplayName 13 MDCS_ASSAY_GetAttributeValueByDisplayName 13 MDCS_ASSAY_AssignAttributeValueLong 13 MDCS_ASSAY_AssignAttributeValueIng 13 MDCS_ASSAY_AssignAttributeValueFloat 13 MDCS_ASSAY_AssignAttributeValueFloat 13 MDCS_ASSAY_RenameAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_UpdateMeasurementValueSbyPlate 13 MDCS_ASSAY_UpdateMeasurementValueSbyPlate 13 Chapter 8: Common Database Functions 13 MDCS_ASSAY_GetUniqueMeasurementValueSbyPlate 13 Chapter 8: Common Database Functions 14 MDCS_DATABASE_Optimize 14 MDCS_DATABASE_CountActiveConnections 14 MDCS_DATABASE_CountActiveConnections 14 MDCS_DATABASE_CountActiveConnections 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedDataEx 14 MDCS_UTILS_RemoveMarkedDataEx 14 MDCS_UTILS_DoesObjectExist 14 MDCS_UTILS_CreateFable 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UPD_UpdateHistoryRecord 14 MDCS_UPD_UpdateDatabaseVersion 15 Chapter 9: Plate Functions 15 MDCS_PLATE_GetAllOrderedByAttributes 15 MDCS_PLATE_GetAllOrderedByAttributes 15 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetInfoBasedOnAssay 156	MDCS_ASSAY_UpdateMeasurementData	. 129
MDCS_ASSAY_CetattributeInfoByDisplayName. 13 MDCS_ASSAY_GetAttributeInfoByDBName 13 MDCS_ASSAY_GetAttributeValueByDBName. 13 MDCS_ASSAY_GetAttributeValueByDBName. 13 MDCS_ASSAY_GetAttributeValueByDisplayName 13 MDCS_ASSAY_GetAttributeValueByDisplayName 13 MDCS_ASSAY_AssignAttributeValueLong 13 MDCS_ASSAY_AssignAttributeValueIng 13 MDCS_ASSAY_AssignAttributeValueFloat 13 MDCS_ASSAY_AssignAttributeValueFloat 13 MDCS_ASSAY_RenameAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_UpdateMeasurementValueSbyPlate 13 MDCS_ASSAY_UpdateMeasurementValueSbyPlate 13 Chapter 8: Common Database Functions 13 MDCS_ASSAY_GetUniqueMeasurementValueSbyPlate 13 Chapter 8: Common Database Functions 14 MDCS_DATABASE_Optimize 14 MDCS_DATABASE_CountActiveConnections 14 MDCS_DATABASE_CountActiveConnections 14 MDCS_DATABASE_CountActiveConnections 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedDataEx 14 MDCS_UTILS_RemoveMarkedDataEx 14 MDCS_UTILS_DoesObjectExist 14 MDCS_UTILS_CreateFable 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UPD_UpdateHistoryRecord 14 MDCS_UPD_UpdateDatabaseVersion 15 Chapter 9: Plate Functions 15 MDCS_PLATE_GetAllOrderedByAttributes 15 MDCS_PLATE_GetAllOrderedByAttributes 15 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetInfoBasedOnAssay 156	MDCS_ASSAY_GetAllMSetUniqueAnnotation	. 130
MDCS_ASSAY_GetAttributeInfoByDIsplayName. 13 MDCS_ASSAY_GetAttributeInfoByDBName. 13 MDCS_ASSAY_GetAttributeValueByDBName. 13 MDCS_ASSAY_GetAttributeValueByDIsplayName 13 MDCS_ASSAY_AssignAttributeValueString 13 MDCS_ASSAY_AssignAttributeValueIong 13 MDCS_ASSAY_AssignAttributeValueIong 13 MDCS_ASSAY_AssignAttributeValueFloat 13 MDCS_ASSAY_RenameAttribute 13 MDCS_ASSAY_RenameAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_UpdateMeasurementSetName 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_UpdateMeasurementValuesByPlate 13 Chapter 8: Common Database Functions 13 MDCS_DATABASE_CouthActiveConnections 14 MDCS_DATABASE_Contpact 14 MDCS_DATABASE_CouthActiveConnections 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetVersion 14 MDCS_DATABASE_GetVersion 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_Execute 14 MDCS_UTILS_Execute 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTD_CreateHistoryRecord 14 MDCS_UTD_DeftHistoryRecord 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllOrderedByAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155	MDCS_ASSAY_CreateAttribute	. 131
MDCS_ASSAY_GetAttributeValueByDBName		
MDCS_ASSAY_GetAttributeValueByDisplayName 13. MDCS_ASSAY_GetAttributeValueByDisplayName 13. MDCS_ASSAY_AssignAttributeValueString 13. MDCS_ASSAY_AssignAttributeValueLong 13. MDCS_ASSAY_AssignAttributeValueFloat 13. MDCS_ASSAY_AssignAttribute 13. MDCS_ASSAY_RenameAttribute 13. MDCS_ASSAY_DeleteAttribute 13. MDCS_ASSAY_GetHeaderAndFileInfo 13. MDCS_ASSAY_UpdateMeasurementSetName 13. MDCS_ASSAY_UpdateMeasurementSetDescription 13. MDCS_ASSAY_UpdateMeasurementValuesByPlate 13. MDCS_ASSAY_GetUniqueMeasurementValuesByPlate 13. MDCS_ASSAY_GetUniqueMeasurementValuesByPlate 13. MDCS_DATABASE_Optimize 14. MDCS_DATABASE_Optimize 14. MDCS_DATABASE_Compact 14. MDCS_DATABASE_CountActiveConnections 14. MDCS_DATABASE_GetSize 14. MDCS_DATABASE_GetVersion 14. MDCS_DATABASE_GetVersion 14. MDCS_DATABASE_GetVersion 14. MDCS_UTILS_RemoveMarkedData 14. MDCS_UTILS_RemoveMarkedData 14. MDCS_UTILS_RemoveMarkedData 14. MDCS_UTILS_RemoveMarkedData 14. MDCS_UTILS_DropTable 14. MDCS_UTILS_DropTable 14. MDCS_UTILS_Execute 14. MDCS_UTILS_Execute 14. MDCS_UTILS_Execute 14. MDCS_UTILS_CreateTable 14. MDCS_UTILS_CreateTable 14. MDCS_UTILS_CreateTable 14. MDCS_UTILS_CreateForeignKey 14. MDCS_UTILS_CreateForeignKey 14. MDCS_UTILS_CreateForeignKey 14. MDCS_UTILS_CreateForeignKey 14. MDCS_UTILS_CreateForeignKey 14. MDCS_UPD_UpdateHistoryRecord 14. MDCS_UPD_GetHistoryRecord 14. MDCS_UPD_GetHistoryRecord 14. MDCS_UPD_UpdateDatabaseVersion 15. Chapter 9: Plate Functions 15. MDCS_PLATE_GetInfo 15. MDCS_PLATE_GetInfo 15. MDCS_PLATE_GetInfoBasedOnAssay 15.		
MDCS_ASSAY_AssignAttributeValueString 13: MDCS_ASSAY_AssignAttributeValueLong 13: MDCS_ASSAY_AssignAttributeValueFloat 13: MDCS_ASSAY_AssignAttributeValueFloat 13: MDCS_ASSAY_RenameAttribute 13: MDCS_ASSAY_RenameAttribute 13: MDCS_ASSAY_DeleteAttribute 13: MDCS_ASSAY_DeleteAttribute 13: MDCS_ASSAY_UpdateMeasurementSetName 13: MDCS_ASSAY_UpdateMeasurementSetDescription 13: MDCS_ASSAY_UpdateMeasurementSetDescription 13: MDCS_ASSAY_UpdateMeasurementValuesByPlate 13: Chapter 8: Common Database Functions 13: MDCS_DATABASE_Optimize 14: MDCS_DATABASE_Compact 14: MDCS_DATABASE_CountActiveConnections 14: MDCS_DATABASE_GetSize 14: MDCS_DATABASE_GetSize 14: MDCS_DATABASE_GetVersion 14: MDCS_DATABASE_GetVersion 14: MDCS_UTILS_RemoveMarkedData 14: MDCS_UTILS_PropTable 14: MDCS_UTILS_DropTable 14: MDCS_UTILS_DoesObjectExist 14: MDCS_UTILS_DesObjectExist 14: MDCS_UTILS_CreateTable 15: MDCS_PLATE_GetInfo 15: MDCS_PLATE_GetAllorderedByAttributes 15: MDCS_PLATE_GetAllorderedByAttributes 15: MDCS_PLATE_GetInfoBasedOnAssay 15:		
MDCS_ASSAY_AssignAttributeValueString 13: MDCS_ASSAY_AssignAttributeValueLong 13: MDCS_ASSAY_AssignAttributeValueFloat 13: MDCS_ASSAY_RenameAttribute 13: MDCS_ASSAY_DeleteAttribute 13: MDCS_ASSAY_DeleteAttribute 13: MDCS_ASSAY_GetHeaderAndFileInfo 13: MDCS_ASSAY_UpdateMeasurementSetName 13: MDCS_ASSAY_UpdateMeasurementSetDescription 13: MDCS_ASSAY_UpdateMeasurementSetDescription 13: MDCS_ASSAY_UpdateMeasurementValuesByPlate 13: Chapter 8: Common Database Functions 13: MDCS_DATABASE_Optimize 14: MDCS_DATABASE_Compact 14: MDCS_DATABASE_CountActiveConnections 14: MDCS_DATABASE_CountActiveConnections 14: MDCS_DATABASE_GetSize 14: MDCS_DATABASE_GetSize 14: MDCS_DATABASE_GetVersion 14: MDCS_UTILS_RemoveMarkedData 14: MDCS_UTILS_RemoveMarkedData 14: MDCS_UTILS_DropTable 14: MDCS_UTILS_DoesObjectExist 14: MDCS_UTILS_DoesObjectExist 14: MDCS_UTILS_Execute 14: MDCS_UTILS_CreateTable 14: MDCS_UTILS_CreateForeignKey 14: MDCS_UPD_GetHistoryRecord 14: MDCS_UPD_GetHistoryRecord 14: MDCS_UPD_GetHistoryRecord 14: MDCS_UPD_GetHistoryRecord 14: MDCS_UPD_UpdateDatabaseVersion 15: Chapter 9: Plate Functions 15: MDCS_PLATE_GetAllPropertyAttributes 15: MDCS_PLATE_GetAllPropertyAttributes 15: MDCS_PLATE_GetAllPropertyAttributes 15: MDCS_PLATE_GetInfoBasedOnAssay 15:		
MDCS_ASSAY_AssignAttributeValueLong MDCS_ASSAY_AssignAttributeValueFloat MDCS_ASSAY_AssignAttributeValueFloat MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_DeleteAttribute 13 MDCS_ASSAY_UpdateMeasurementSetName 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_UpdateMeasurementValuesByPlate 13 MDCS_ASSAY_GetUniqueMeasurementValuesByPlate 13 Chapter 8: Common Database Functions 13 MDCS_DATABASE_Optimize 14 MDCS_DATABASE_Compact 14 MDCS_DATABASE_Compact 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetVersion 14 MDCS_DATABASE_GetVersion 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_DropTable 14 MDCS_UTILS_DropTable 14 MDCS_UTILS_DoesObjectExist 14 MDCS_UTILS_CreateTable 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UPD_GetHistoryRecord 14 MDCS_UPD_GetHistoryRecord 15 MDCS_UPD_FindFinishedUpdates 16 MDCS_UPD_FindFinishedUpdates 17 MDCS_UPD_FindFinishedUpdates 18 MDCS_PLATE_GetInfo 19 MDCS_PLATE_GetInfo 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetInfoBasedOnAssay 15		
MDCS_ASSAY_AssignAttributeValueFloat		
MDCS_ASSAY_RenameAttribute	MDCS ASSAY AssignAttributeValueFloat	. 135
MDCS_ASSAY_GetHeaderAndFileInfo. 13 MDCS_ASSAY_UpdateMeasurementSetName 13 MDCS_ASSAY_UpdateMeasurementSetDescription 13 MDCS_ASSAY_GetUniqueMeasurementValuesByPlate 13 MDCS_ASSAY_GetUniqueMeasurementValuesByPlate 13 Chapter 8: Common Database Functions 13 MDCS_DATABASE_Optimize 14 MDCS_DATABASE_Compact 14 MDCS_DATABASE_CountActiveConnections 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetAvailableDatabases 14 MDCS_DATABASE_GetVersion 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedDataEx 14 MDCS_UTILS_DropTable 14 MDCS_UTILS_DropTable 14 MDCS_UTILS_DoesObjectExist 14 MDCS_UTILS_Execute 14 MDCS_UTILS_CreateTable 14 MDCS_UTILS_CreateTable 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UPD_CreateHistoryRecord 14 MDCS_UPD_GetHistoryRecord 14 MDCS_UPD_GetHistoryRecord 14 MDCS_UPD_GetHistoryRecord 14 MDCS_UPD_GetHistoryRecord 15 MDCS_UPD_UpdateDatabaseVersion 15 Chapter 9: Plate Functions 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetInfoBasedOnAssay 156		
MDCS_ASSAY_UpdateMeasurementSetName 13: MDCS_ASSAY_UpdateMeasurementSetDescription 13: MDCS_ASSAY_GetUniqueMeasurementValuesByPlate 13: Chapter 8: Common Database Functions 13: MDCS_DATABASE_Optimize 14: MDCS_DATABASE_Compact 14: MDCS_DATABASE_CountActiveConnections 14: MDCS_DATABASE_GetSize 14: MDCS_DATABASE_GetSize 14: MDCS_DATABASE_GetAvailableDatabases 14: MDCS_DATABASE_GetVersion 14: MDCS_UTILS_RemoveMarkedData 14: MDCS_UTILS_RemoveMarkedData 14: MDCS_UTILS_RemoveMarkedDataEx 14: MDCS_UTILS_DropTable 14: MDCS_UTILS_DropTable 14: MDCS_UTILS_Execute 14: MDCS_UTILS_Execute 14: MDCS_UTILS_CreateTable 14: MDCS_UTILS_CreateForeignKey 14: MDCS_UTILS_CreateForeignKey 14: MDCS_UTILS_CreateForeignKey 14: MDCS_UTILS_CreateHistoryRecord 14: MDCS_UPD_CreateHistoryRecord 14: MDCS_UPD_GetHistoryRecord 14: MDCS_UPD_GetHistoryRecord 14: MDCS_UPD_GetHistoryRecord 15: MDCS_UPD_UpdateDatabaseVersion 15: Chapter 9: Plate Functions 15: MDCS_PLATE_GreatePlate 15: MDCS_PLATE_GetAllPropertyAttributes 15: MDCS_PLATE_GetAllPropertyAttributes 15: MDCS_PLATE_GetAllOrderedByAttributes 15: MDCS_PLATE_GetInfoBasedOnAssay 15:	MDCS_ASSAY_DeleteAttribute	. 136
MDCS_ASSAY_UpdateMeasurementSetDescription MDCS_ASSAY_GetUniqueMeasurementValuesByPlate Chapter 8: Common Database Functions MDCS_DATABASE_Optimize MDCS_DATABASE_Compact MDCS_DATABASE_CountActiveConnections 14: MDCS_DATABASE_GetSize MDCS_DATABASE_GetSize MDCS_DATABASE_GetVersion MDCS_UTILS_RemoveMarkedData MDCS_UTILS_RemoveMarkedData MDCS_UTILS_DropTable MDCS_UTILS_DoesObjectExist MDCS_UTILS_Execute MDCS_UTILS_CreateTable MDCS_UTILS_CreateForeignKey MDCS_UTILS_CreateForeignKey MDCS_UTILS_CreateHistoryRecord MDCS_UPD_CreateHistoryRecord MDCS_UPD_GetHistoryRecord MDCS_UPD_GetHistoryRecord MDCS_UPD_UpdateDatabaseVersion Chapter 9: Plate Functions MDCS_PLATE_GetAllPropertyAttributes MDCS_PLATE_GetAllOrderedByAttributes	MDCS_ASSAY_GetHeaderAndFileInfo	. 136
MDCS_ASSAY_GetUniqueMeasurementValuesByPlate 138 Chapter 8: Common Database Functions 139 MDCS_DATABASE_Optimize 144 MDCS_DATABASE_Compact 144 MDCS_DATABASE_CountActiveConnections 144 MDCS_DATABASE_GetSize 144 MDCS_DATABASE_GetSize 144 MDCS_DATABASE_GetVersion 144 MDCS_DATABASE_GetVersion 144 MDCS_UTILS_RemoveMarkedData 144 MDCS_UTILS_DropTable 144 MDCS_UTILS_DropTable 144 MDCS_UTILS_DoesObjectExist 144 MDCS_UTILS_CreateTable 144 MDCS_UTILS_CreateTable 144 MDCS_UTILS_CreateForeignKey 144 MDCS_UTILS_CreateForeignKey 144 MDCS_UTILS_CreateForeignKey 144 MDCS_UTILS_CreateHistoryRecord 144 MDCS_UPD_CreateHistoryRecord 144 MDCS_UPD_GetHistoryRecord 144 MDCS_UPD_GetHistoryRecord 144 MDCS_UPD_GetHistoryRecord 144 MDCS_UPD_GetHistoryRecord 145 MDCS_UPD_UpdateDatabaseVersion 155 Chapter 9: Plate Functions 155 MDCS_PLATE_GetAllPropertyAttributes 155 MDCS_PLATE_GetAllPropertyAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetAllOrderedByAttributes 155 MDCS_PLATE_GetInfoBasedOnAssay 156	MDCS_ASSAY_UpdateMeasurementSetName	. 137
Chapter 8: Common Database Functions139MDCS_DATABASE_Optimize144MDCS_DATABASE_Compact144MDCS_DATABASE_CountActiveConnections144MDCS_DATABASE_GetSize144MDCS_DATABASE_GetAvailableDatabases145MDCS_DATABASE_GetVersion144MDCS_UTILS_RemoveMarkedData144MDCS_UTILS_RemoveMarkedData144MDCS_UTILS_RemoveMarkedDataEx144MDCS_UTILS_DropTable144MDCS_UTILS_DoesObjectExist144MDCS_UTILS_Execute144MDCS_UTILS_CreateTable144MDCS_UTILS_AddColumnToTable144MDCS_UTILS_CreateForeignKey144MDCS_UTILS_CreateHistoryRecord144MDCS_UPD_CreateHistoryRecord144MDCS_UPD_CreateHistoryRecord144MDCS_UPD_GetHistoryRecord144MDCS_UPD_GetHistoryRecord144MDCS_UPD_FindFinishedUpdates144MDCS_UPD_FindFinishedUpdates144MDCS_UPD_TeateFunctions155Chapter 9: Plate Functions155MDCS_PLATE_CreatePlate155MDCS_PLATE_GetAllPropertyAttributes156MDCS_PLATE_GetAllPropertyAttributes156MDCS_PLATE_GetAllPropertyAttributes156MDCS_PLATE_GetInfoBasedOnAssay156	MDCS_ASSAY_UpdateMeasurementSetDescription	. 137
MDCS_DATABASE_Optimize	MDCS_ASSAY_GetUniqueMeasurementValuesByPlate	. 138
MDCS_DATABASE_Optimize	Chautau Ca Causana Batabana Functiona	120
MDCS_DATABASE_Compact MDCS_DATABASE_CountActiveConnections MDCS_DATABASE_GetSize MDCS_DATABASE_GetAvailableDatabases MDCS_DATABASE_GetAvailableDatabases MDCS_DATABASE_GetVersion MDCS_UTILS_RemoveMarkedData MDCS_UTILS_RemoveMarkedData MDCS_UTILS_DropTable MDCS_UTILS_DropTable MDCS_UTILS_DropTable MDCS_UTILS_Execute MDCS_UTILS_Execute MDCS_UTILS_CreateTable MDCS_UTILS_CreateTable MDCS_UTILS_CreateForeignKey MDCS_UTILS_CreateForeignKey MDCS_UTILS_CreateForeignKey MDCS_UTILS_CreateHistoryRecord MDCS_UPD_CreateHistoryRecord MDCS_UPD_UpdateHistoryRecord MDCS_UPD_GetHistoryRecord MDCS_UPD_GetHistoryRecord MDCS_UPD_FindFinishedUpdates MDCS_UPD_UpdateDatabaseVersion Chapter 9: Plate Functions MDCS_PLATE_CreatePlate MDCS_PLATE_GetInfo MDCS_PLATE_GetAllPropertyAttributes MDCS_PLATE_GetAllPropertyAttributes MDCS_PLATE_GetAllOrderedByAttributes MDCS_PLATE_GetInfoBasedOnAssay 156		
MDCS_DATABASE_CountActiveConnections 14 MDCS_DATABASE_GetSize 14 MDCS_DATABASE_GetAvailableDatabases 14 MDCS_DATABASE_GetVersion 14 MDCS_UTILS_RemoveMarkedData 14 MDCS_UTILS_RemoveMarkedDataEx 14 MDCS_UTILS_DropTable 14 MDCS_UTILS_DoesObjectExist 14 MDCS_UTILS_Execute 14 MDCS_UTILS_CreateTable 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateForeignKey 14 MDCS_UTILS_CreateHistoryRecord 14 MDCS_UPD_CreateHistoryRecord 14 MDCS_UPD_GetHistoryRecord 14 MDCS_UPD_GetHistoryRecord 14 MDCS_UPD_GetHistoryRecord 15 MDCS_UPD_UpdateDatabaseVersion 15 Chapter 9: Plate Functions 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllPropertyAttributes 15 MDCS_PLATE_GetAllOrderedByAttributes 15 MDCS_PLATE_GetInfoBasedOnAssay 156		
MDCS_DATABASE_GetSize		
MDCS_DATABASE_GetAvailableDatabases 14: MDCS_DATABASE_GetVersion 14: MDCS_UTILS_RemoveMarkedData 14: MDCS_UTILS_RemoveMarkedDataEx 14: MDCS_UTILS_DropTable 14: MDCS_UTILS_DoesObjectExist 14: MDCS_UTILS_Execute 14: MDCS_UTILS_CreateTable 14: MDCS_UTILS_CreateTable 14: MDCS_UTILS_CreateForeignKey 14: MDCS_UTILS_CreateForeignKey 14: MDCS_UTILS_CreateStoredProc 14: MDCS_UPD_CreateHistoryRecord 14: MDCS_UPD_UpdateHistoryRecord 14: MDCS_UPD_GetHistoryRecord 14: MDCS_UPD_GetHistoryRecord 14: MDCS_UPD_UpdateDatabaseVersion 15: Chapter 9: Plate Functions 15: MDCS_PLATE_CreatePlate 15: MDCS_PLATE_GetAllPropertyAttributes 15: MDCS_PLATE_GetAllPropertyAttributes 15: MDCS_PLATE_GetInfoBasedOnAssay 15:		
MDCS_DATABASE_GetVersion		
MDCS_UTILS_RemoveMarkedData14MDCS_UTILS_RemoveMarkedDataEx14MDCS_UTILS_DropTable14MDCS_UTILS_DoesObjectExist14MDCS_UTILS_Execute14MDCS_UTILS_CreateTable14MDCS_UTILS_CreateTable14MDCS_UTILS_CreateForeignKey14MDCS_UTILS_CreateStoredProc14MDCS_UTILS_CreateHistoryRecord14MDCS_UPD_CreateHistoryRecord14MDCS_UPD_UpdateHistoryRecord14MDCS_UPD_GetHistoryRecord14MDCS_UPD_FindFinishedUpdates14MDCS_UPD_UpdateDatabaseVersion15Chapter 9: Plate Functions15MDCS_PLATE_CreatePlate15MDCS_PLATE_GetInfo15MDCS_PLATE_GetAllPropertyAttributes15MDCS_PLATE_GetAllOrderedByAttributes15MDCS_PLATE_GetInfoBasedOnAssay15		
MDCS_UTILS_RemoveMarkedDataEx14MDCS_UTILS_DropTable14MDCS_UTILS_DoesObjectExist14MDCS_UTILS_Execute14MDCS_UTILS_CreateTable14MDCS_UTILS_AddColumnToTable14MDCS_UTILS_CreateForeignKey14MDCS_UTILS_CreateStoredProc14MDCS_UPD_CreateHistoryRecord14MDCS_UPD_UpdateHistoryRecord14MDCS_UPD_GetHistoryRecord14MDCS_UPD_FindFinishedUpdates14MDCS_UPD_UpdateDatabaseVersion15Chapter 9: Plate FunctionsMDCS_PLATE_CreatePlate15MDCS_PLATE_GetInfo15MDCS_PLATE_GetAllPropertyAttributes15MDCS_PLATE_GetAllOrderedByAttributes15MDCS_PLATE_GetInfoBasedOnAssay15		
MDCS_UTILS_DropTable		
MDCS_UTILS_DoesObjectExist		
MDCS_UTILS_Execute14MDCS_UTILS_CreateTable14MDCS_UTILS_AddColumnToTable14MDCS_UTILS_CreateForeignKey14MDCS_UTILS_CreateStoredProc14MDCS_UPD_CreateHistoryRecord14MDCS_UPD_UpdateHistoryRecordStatus14MDCS_UPD_GetHistoryRecord14MDCS_UPD_FindFinishedUpdates14MDCS_UPD_UpdateDatabaseVersion15Chapter 9: Plate Functions15MDCS_PLATE_CreatePlate15MDCS_PLATE_GetInfo15MDCS_PLATE_GetAllPropertyAttributes15MDCS_PLATE_GetAllOrderedByAttributes15MDCS_PLATE_GetInfoBasedOnAssay15		
MDCS_UTILS_CreateTable14!MDCS_UTILS_AddColumnToTable14!MDCS_UTILS_CreateForeignKey14!MDCS_UTILS_CreateStoredProc14!MDCS_UPD_CreateHistoryRecord14!MDCS_UPD_UpdateHistoryRecordStatus14!MDCS_UPD_GetHistoryRecord14!MDCS_UPD_FindFinishedUpdates14!MDCS_UPD_UpdateDatabaseVersion15!Chapter 9: Plate FunctionsMDCS_PLATE_CreatePlate15:MDCS_PLATE_GetInfo15:MDCS_PLATE_GetAllPropertyAttributes15:MDCS_PLATE_GetAllOrderedByAttributes15:MDCS_PLATE_GetInfoBasedOnAssay15:		
MDCS_UTILS_AddColumnToTable		
MDCS_UTILS_CreateForeignKey140MDCS_UTILS_CreateStoredProc141MDCS_UPD_CreateHistoryRecord142MDCS_UPD_UpdateHistoryRecordStatus143MDCS_UPD_GetHistoryRecord143MDCS_UPD_FindFinishedUpdates143MDCS_UPD_UpdateDatabaseVersion150Chapter 9: Plate FunctionsMDCS_PLATE_CreatePlate153MDCS_PLATE_GetInfo154MDCS_PLATE_GetAllPropertyAttributes154MDCS_PLATE_GetAllOrderedByAttributes154MDCS_PLATE_GetInfoBasedOnAssay156		
MDCS_UTILS_CreateStoredProc14MDCS_UPD_CreateHistoryRecord14MDCS_UPD_UpdateHistoryRecordStatus14MDCS_UPD_GetHistoryRecord14MDCS_UPD_FindFinishedUpdates14MDCS_UPD_UpdateDatabaseVersion15Chapter 9: Plate FunctionsMDCS_PLATE_CreatePlate15MDCS_PLATE_GetInfo15MDCS_PLATE_GetAllPropertyAttributes15MDCS_PLATE_GetAllOrderedByAttributes15MDCS_PLATE_GetInfoBasedOnAssay15		
MDCS_UPD_CreateHistoryRecord14MDCS_UPD_UpdateHistoryRecordStatus148MDCS_UPD_GetHistoryRecord148MDCS_UPD_FindFinishedUpdates149MDCS_UPD_UpdateDatabaseVersion150Chapter 9: Plate FunctionsMDCS_PLATE_CreatePlate153MDCS_PLATE_GetInfo154MDCS_PLATE_GetAllPropertyAttributes154MDCS_PLATE_GetAllOrderedByAttributes155MDCS_PLATE_GetInfoBasedOnAssay156		
MDCS_UPD_UpdateHistoryRecordStatus148MDCS_UPD_GetHistoryRecord148MDCS_UPD_FindFinishedUpdates149MDCS_UPD_UpdateDatabaseVersion150Chapter 9: Plate FunctionsMDCS_PLATE_CreatePlate153MDCS_PLATE_GetInfo154MDCS_PLATE_GetAllPropertyAttributes154MDCS_PLATE_GetAllOrderedByAttributes155MDCS_PLATE_GetInfoBasedOnAssay156		
MDCS_UPD_GetHistoryRecord148MDCS_UPD_FindFinishedUpdates149MDCS_UPD_UpdateDatabaseVersion150Chapter 9: Plate FunctionsMDCS_PLATE_CreatePlate150MDCS_PLATE_GetInfo150MDCS_PLATE_GetAllPropertyAttributes150MDCS_PLATE_GetAllOrderedByAttributes150MDCS_PLATE_GetInfoBasedOnAssay150	MDCS_UPD_CreateHistoryRecord	. 14/
MDCS_UPD_FindFinishedUpdates149MDCS_UPD_UpdateDatabaseVersion150Chapter 9: Plate Functions153MDCS_PLATE_CreatePlate153MDCS_PLATE_GetInfo154MDCS_PLATE_GetAllPropertyAttributes154MDCS_PLATE_GetAllOrderedByAttributes154MDCS_PLATE_GetInfoBasedOnAssay154		
MDCS_UPD_UpdateDatabaseVersion150Chapter 9: Plate Functions151MDCS_PLATE_CreatePlate152MDCS_PLATE_GetInfo154MDCS_PLATE_GetAllPropertyAttributes154MDCS_PLATE_GetAllOrderedByAttributes154MDCS_PLATE_GetInfoBasedOnAssay156		
Chapter 9: Plate Functions		
MDCS_PLATE_CreatePlate153MDCS_PLATE_GetInfo154MDCS_PLATE_GetAllPropertyAttributes154MDCS_PLATE_GetAllOrderedByAttributes154MDCS_PLATE_GetInfoBasedOnAssay156	MDCS_UPD_updateDatabaseversion	. 150
MDCS_PLATE_CreatePlate153MDCS_PLATE_GetInfo154MDCS_PLATE_GetAllPropertyAttributes154MDCS_PLATE_GetAllOrderedByAttributes154MDCS_PLATE_GetInfoBasedOnAssay156	Chapter 9: Plate Functions	151
MDCS_PLATE_GetInfo	-	
MDCS_PLATE_GetAllPropertyAttributes		
MDCS_PLATE_GetAllOrderedByAttributes		
MDCS_PLATE_GetInfoBasedOnAssay		

MDCS_PLATE_	_DeleteImages	157
MDCS_PLATE_	_GetUniqueAttributeValues	157
MDCS_PLATE_	_GetAllByAttributes	158
MDCS_PLATE_	_ManageSecurity	158
MDCS_PLATE_	_FolderManageSecurity	159
MDCS_PLATE_	_ShareSiteImagesByTimeAndZIndex	160
MDCS_PLATE_	_CreateSite	160
MDCS_PLATE_	_CreateSeries	161
	_CreateImageSource	
	_CreateImageRecord	
	_ImportLayoutData	
MDCS_PLATE_	_ApplyLayoutToAssay	163
	_DeleteLayout	
	_GetAllRecords	
	_GetRecord	
	_GetSiteRecord	
	_GetSitesByPlate	
	_GetImageRecord	
	_GetImageRecordPerPlate	
	_UpdateImageObjectID	
	_GetImageObjectIDForImage	
	_GetSeriesRecord	
	_GetSiteLocationsForPlate	
	_GetImageSourcesOfPlate	
	_GetMaxTimePointForPlate	
	_GetPlatesByDate	
	_GetFlatesbyDate	
	_opuateAcquisition	
	_GetAttributeInfoByDisplayName	
	_GetAttributeInfoByDBName	
	_AssignAttributeValueString	
	_AssignAttributeValueLong	
	AssignAttributeValueFloat	
	_UpdateInfo	
	_GetInfoBasedOnAssaySet	
MDCS PLATE	CreateLayout	179
MDCS PLATE	_CreateLayout	180
MDCS PLATE		180
MDCS PLATE	GetTemplateByName	181
	 _GetTemplates	
	 _RenameAttribute	
MDCS_PLATE	CountPlateDatasets	183
	_ _GetAcqSiteCount	
	 _GetAcqWellCount	
	 _GetAcqSeriesCount	

MDCS_PLATE_GetCompoundCount MDCS_PLATE_GetControlsCount . MDCS_PLATE_GetControlStatistic . MDCS_PLATE_GetCompleteImageInfo . MDCS_PLATE_CanModify . MDCS_PLATE_GetImageIDs . MDCS_PLATE_GetThumbImageIDs . MDCS_PLATE_ChangeStatus .	185 186 186 187 187 188
Chapter 10: Acquisitions Functions	191
MDCS_ACQUISITION_Delete MDCS_ACQUISITION_Create MDCS_ACQUISITION_CreateProfile MDCS_ACQUISITION_GetProfileRecords MDCS_ACQUISITION_GetProfile MDCS_ACQUISITION_GetInstanceRecord MDCS_ACQUISITION_GetBatchRecords MDCS_ACQUISITION_GetBatchRecord	192 193 193 194 195
Chapter 11: Security Functions	197
MDCS_SECURITY_ManageUsersDlg	198 199 199 200 200 201
Chapter 12: Datasets Functions	203
MDCS_DATASET_InsertAnalysisAttributes MDCS_DATASET_GetAnalysisAttributesByDataset MDCS_DATASET_GetAnalysisAttributes MDCS_DATASET_GetAnalysisDescriptions MDCS_DATASET_FindAnalysisAttributesRecord MDCS_DATASET_DeleteAnalysisAttributes MDCS_DATASET_UpdateAnalysisAttributes MDCS_DATASET_ManageFolderSecurity MDCS_DATASET_DeleteFolder MDCS_DATASET_Copy MDCS_DATASET_Copy MDCS_DATASET_Update MDCS_DATASET_Get MDCS_DATASET_Get MDCS_DATASET_Create MDCS_DATASET_AddAssay MDCS_DATASET_HavePermissionsToModify MDCS_DATASET_GetAnalysisDescription MDCS_DATASET_GetAssayAndFilters	206 207 208 209 209 210 210 211 211 212 212 213 213 214 214

MDCS_DATASET_GetItemFolder	. 215
MDCS_DATASET_Gettterni older	
MDCS_DATASET_GetAssaySiblingFolders	. 217
MDCS_DATASET_DoesSubFolderExist	
MDCS_DATASET_CreateFolder	
MDCS_DATASET_ModifyFolder	
MDCS_DATASET_GetSiblingFolders	
MDCS_DATASET_GetAllMeasurements	
MDCS_DATASET_AddScriptletAssay	
MDCS_DATASET_GetScriptletAssays	
MDCS_DATASET_GetScriptletAssayIDByName	
MDCS_DATASET_GetAllForPlate	
MDCS_DATASET_GetAllForAssay	
MDCS_DATASET_GetAllMSetParamValues	
MDCS_DATASET_GetResultsInfoByConfig	
MDCS_DATASET_GetResultInfo	
MDCS_DATASET_CallbackToAnalysisInfo	
MDCS_DATASET_DoesNameExist	
MDCS_DATASET_UpdateFolderItem	
MDCS_DATASET_Delete	
MDCS_DATASET_GetAnalysisCount	
—	
Chapter 13: Assay Images and Assay Normalization	
Functions	
MDCS_ASSAYIMAGES_GetAllPropertyAttributes	230
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231 . 232 . 233
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231 . 232 . 233
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231 . 232 . 233 . 234
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231 . 232 . 233 . 234 . 234
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231 . 232 . 233 . 234 . 234 . 235
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231 . 232 . 233 . 234 . 234 . 235
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231 . 232 . 233 . 234 . 234 . 235
MDCS_ASSAYIMAGES_GetAllByAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions	. 231 . 232 . 233 . 234 . 234 . 235 . 235
MDCS_ASSAYIMAGES_GetAllByAttributes	. 231 . 232 . 233 . 234 . 234 . 235 . 235
MDCS_ASSAYIMAGES_GetAllByAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions	. 231 . 232 . 233 . 234 . 234 . 235 . 235
MDCS_ASSAYIMAGES_GetAllByAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions	. 231 . 232 . 233 . 234 . 234 . 235 . 235
MDCS_ASSAYIMAGES_GetAllByAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions MDCS_APP_ManageLocationOptions	. 231 . 232 . 233 . 234 . 234 . 235 . 237 . 237
MDCS_ASSAYIMAGES_GetUniqueAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions MDCS_APP_ManageLocationOptions MDCS_APP_CreateLocationOption	. 231 . 232 . 233 . 234 . 234 . 235 . 237 . 237
MDCS_ASSAYIMAGES_GetUniqueAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions MDCS_APP_ManageLocationOptions MDCS_APP_CreateLocationOption MDCS_APP_CreateJobRecord	. 231 . 232 . 233 . 234 . 235 . 235 . 237 . 237 . 240 . 241 . 241
MDCS_ASSAYIMAGES_GetAllByAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions MDCS_APP_ManageLocationOptions MDCS_APP_CreateLocationOption MDCS_APP_CreateJobRecord MDCS_APP_CreateParameter	. 231 . 232 . 233 . 234 . 235 . 235 . 237 . 237 . 241 . 241 . 242
MDCS_ASSAYIMAGES_GetUniqueAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions. MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions. MDCS_APP_ManageLocationOptions MDCS_APP_CreateLocationOption MDCS_APP_CreateJobRecord MDCS_APP_CreateParameter MDCS_APP_GetParameterByName	. 231 . 232 . 233 . 234 . 235 . 235 . 237 . 237 . 241 . 241 . 242 . 242
MDCS_ASSAYIMAGES_GetAllByAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig. MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions. MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions. MDCS_APP_ManageLocationOptions MDCS_APP_CreateLocationOption MDCS_APP_CreateJobRecord MDCS_APP_CreateParameter MDCS_APP_GetParameterByName MDCS_APP_GetParameterByName MDCS_APP_GetParameterByName	. 231 . 232 . 233 . 234 . 235 . 235 . 237 . 237 . 240 . 241 . 242 . 242
MDCS_ASSAYIMAGES_GetUniqueAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions. MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions. MDCS_APP_ManageLocationOptions MDCS_APP_CreateLocationOption MDCS_APP_CreateJobRecord MDCS_APP_CreateParameter MDCS_APP_GetParameterByName MDCS_APP_GetParameterByName MDCS_APP_GetParameterByID	. 231 . 232 . 233 . 234 . 235 . 235 . 235 . 237 . 240 . 241 . 242 . 243 . 243
MDCS_ASSAYIMAGES_GetAllByAttributes MDCS_ASSAYIMAGES_GetUniqueAttributeValues MDCS_ASSAYNORM_CreateConfig MDCS_ASSAYNORM_UpdateConfig. MDCS_ASSAYNORM_GetConfig MDCS_ASSAYNORM_GetConfigByName MDCS_ASSAYNORM_GetConfigAll MDCS_ASSAYNORM_GetConfigForAssay Chapter 14: Quicklist Functions. MDCS_QUICKLIST_ManageFolderSecurity Chapter 15: Application Data Table Functions. MDCS_APP_ManageLocationOptions MDCS_APP_CreateLocationOption MDCS_APP_CreateJobRecord MDCS_APP_CreateParameter MDCS_APP_GetParameterByName MDCS_APP_GetParameterByName MDCS_APP_GetParameterByName	. 231 . 232 . 233 . 234 . 235 . 235 . 235 . 237 . 241 . 241 . 242 . 243 . 243

MDCS_APP_GetLocationOptionRecord	245
MDCS_APP_GetLocationRecordsByLabel	
MDCS_APP_DeleteLocationOption	246
MDCS_APP_DeleteUserLocationOptions	246
MDCS_APP_SetUserLocationOptions	
MDCS_APP_GetUserLocationOption	
MDCS_APP_GetJobQueueRecords	
MDCS_APP_GetJobQueueRecord	
MDCS_APP_GetJobQueue	
MDCS_APP_UpdateJobStatus	
MDCS_APP_UpdateJobProgress	
MDCS_APP_UpdateJobQueueRecord	
MDCS_APP_CancelJobProgress	
MDCS_APP_ClaimJob	
MDCS_APP_ResetJob	
MDCS_APP_ResetCrashedJobs	
MDCS_APP_RefreshAllJobs	
MDCS_APP_GetLocationOptionByID	255
Chapter 16: File Data and Macro Functions	257
MDCS_FILEDATA_GetRecordByID	
MDCS_FILEDATA_GetRecordByGUID	
MDCS_FILEDATA_UpdateRecord	
MDCS MACROS GetRecord	
MDCS_MACROS_GetAllByGroup	
MDCS_MACROS_UpdateRecord	
MDCS MACROS DeleteRecord	
MDCS_MACROS_Create	
MDCS_MACROS_ConvertQueryResultsToMacrosStructure	
Chapter 17: Structures, ENUMS, and Definitions	
•	
Structures	
ENUMS Types	
Definitions	
Error Code Definitions	
Chapter 18: Virtual Callback Classes	
Class MDCS_GroupInfoCallback	
Overview	
Public Methods	
Private Methods	
Function signature of MDCS_GroupInfoCallback	
Class MDCS_SaveBlobCallback	
Overview	
Public Methods	
Private Methods	
Function signatures of MDCS_SaveBlobCallback	271

Class MDCS_GetBlobCallback	273
Overview	273
Public Methods	273
Private Methods	273
Function signature of MDCS_GetBlobCallback	273
Class MDCS_ProgressCallback	
Overview	275
Public Methods	275
Private Methods	
Function signature of MDCS_ProgressCallback	275
Class MDCS_GetDBResultsCCallback	277
Overview	277
Public Methods	
Private Methods	277
Function signature of MDCS_GetDBResultsCCallback	277
Class MDCS_QueryResults	
Overview	
Public Methods	
Private Methods	
Function signature of MDCS_QueryResults	
Class MDCS_ImportMeasurementSet	
Overview	
Public Methods	
Private Methods	
Function signatures of MDCS_ImportMeasurementSet	
Class MDCS_ImportPlateLayout	
Overview	
Public Methods	
Private Methods	
Function signature of MDCS_ImportPlateLayout	
Class MDCS_DBHandleSmartPtr	
Overview	
Data Member	
Public Methods	
Private Methods	
Function signature of MDCS_DBHandleSmartPtr	
Class MDCS_CL_ImportDS	
Public Methods	
Function signature of MDCS_CL_ImportDS	290
Class MDCS_GetProgressStatus	
Public Methods	
Private Methods	
Function signature of MDCS_GetProgressStatus	
Class MDCS_CL_BlobLocationCB	
Public Methods	
Private Methods	
Return	292

class MDCS_CL_BlobInfoCB	
Public Methods	92
Private Methods	92
Return	92
Chapter 19: Usage Examples)3
How to Connect to the Database and Get a Database Handle 29	
How to Use Error Handling	96
How to Get Plate Information	97
How to Get All Sites Per Plate	98
Deriving the MDCS_GetDBResultCCallback Class 29	99
Using the Callback CGetSitesByPlate Callback	01
Deriving the MDCS_GetBlobCallback Class	ງ2
Using the CBlobCallbackEx Callback Class 30	
Index)7

The preface contains the following sections:

- Who This Manual Is For on page 15
- About this Manual on page 15
- Conventions on page 16

Who This Manual Is For

This manual is written for those who wish to build an application that accesses and exchanges data with an MDCStore $^{\text{TM}}$ database. Use of the MDCStore API is only available under the MDCStore API limited licensing agreement.

About this Manual

This manual contains the following chapters:

Title	Overview
Chapter 1 Introduction to the MDCStoreUtils API	Provides an introduction the MDCStore Utilities API
Chapter 2 Database Connection Functions	Contains functions that can be used with database connections
Chapter 3 Results Set Functions	Contains functions that can be used with results sets
Chapter 4 Result Set Data Types Functions	Contains functions that can be used with results set data types
Chapter 5 Error Handling Functions	Contains error handing functions
Chapter 6 BLOB Functions	Contains functions that can be used with BLOBs
Chapter 7 Measurement Sets and Results Functions	Contains functions that can be used with measurement sets and results
Chapter 8 Common Database Functions	Contains common database functions
Chapter 9 Plate Functions	Contains functions that can be used with plates
Chapter 10 Acquisitions Functions	Contains functions that can be used with acquisitions
Chapter 11 Security Functions	Contains security functions
Chapter 12 Datasets Functions	Contains functions that can be used with datasets
Chapter 13 Assay Images and Assay Normalization Functions	Contains functions that can be used with assay images and assay normalization
Chapter 14 Quicklist Functions	Contains quicklist functions
Chapter 15 Application Data Table Functions	Contains functions that can be used with acquisition data table functions

5000957 C 15

Title	Overview
Chapter 16 File Data and Macro Functions	Contains functions that can be used with file data and macros
Chapter 17 Structures, ENUMS, and Definitions	Describes data structures enums, definitions and error codes
Chapter 18 Virtual Callback Classes	Describes virtual call back classes
Chapter 19 Usage Examples	Contains examples of how to use the MDCStore Utilities API

Conventions

Within the scope of this manual, the following typographical conventions are used.



WARNING! A warning indicates an operation that may cause personal injury if precautions are not followed.

CAUTION! Indicates an operation that may cause damage to the instrument, device, or data, if the precautions are not followed.



Note: Provides essential information for the completion of a procedure.



Tip! Provides useful information that helps apply the techniques and procedures in the text to your specific needs, and provides shortcuts, but is not essential to the completion of a procedure.

16 5000957 C

Introduction to the MDCStoreUtils API

What is the MDCStoreUtils API?

The MDCStoreUtils API is an open application programming interface (API) created using ODBC classes to work with the MDCStore $^{\text{TM}}$ database. It provides a set of functions to allow MDC customers to build applications that access and exchange data with the MDCStore database.

The API provides users the following benefits:

- Access to the MDCStore database.
- The ability to use the highest level of security and still be able to conduct business unrestricted.
- Access to the database error-handling mechanism.
- Access to all MDCStore privileges to retrieve data, manipulate data and save data to the database.
- Ease of use and robustness.
- C++ can be used to access the API.

Main Functionality

The MDCStoreUtils API provides a set of definitions, structures and enums and virtual pure callback classes. It also provides functions that work with the following:

- Database connections
- Results sets
- Results set data types
- Database errors
- BLOB data
- Measurement sets and results of measurement sets
- Plates
- Data acquisition
- Security features
- Datasets
- Images in measurement sets
- Quicklists
- Common database functions
- Application data tables
- File data

How to Use the MDCStoreUtils API

Supported database servers

MDCStoreUtils API supports two common database servers:

- Microsoft SQL Server 2005 and Microsoft SQL Server Express 2005
- Oracle 9.2 and higher

Compiling and linking with MDCStoreUtils API

Normally, the MDCStoreUtils API distribution comes with precompiled DLLs and all libraries needed. However, users may want to recompile the source files of this API distribution.

To compile and link application with MDCStoreUtils API:

 At compile time, include the follow header file anywhere to use the functions of MDCStoreUtils API:

```
# include "MDCStoreUtilsAPI.h"
# include "MDCSStructures.h"
# include "MDCSDefinitions.h"
# include "MDCSEnums.h"
# include "AxStringCollection.h"
```

At link time, users need other libraries to link with MDCStoreUtils:

```
AxFileClient.lib
AxMFCDBUtils.lib
AxStringCollection.lib
MDCStoreUtils.lib
```

At runtime, users need the following DLLs to run the program:

```
AxFileClient.dll
AxMFCDBUtils.dll
AxStringCollection.dll
MDCStoreUtils.dll
```

18 5000957 C

Definitions

Dataset

A database query that returns a subset of objects in the database.

Instance

A sequence number that is used to denote multiple similar measurements at the same location.

Measurement Set

All the data extracted from one plate of images by a MetaXpress® application module with specific settings. Running the same application module again with different settings produces a separate measurement set.

Plate

All the images acquired from one physical microplate.

Result Set

A set of rows of data as a result of executing a SQL query.

Series

The location of an image or data value within a series in the z-dimension or the time dimension.

Site

A region within a microplate well from which an image has been acquired. Images can be acquired from more than one site per well.

Functions Available in MDCStoreUtils API

MDCStoreUtils API consists of functions divided into 15 categories. In this section, we provide all the function names, descriptions of what each function does based on these categories, and the function signatures.

There are two special functions you should always call before calling any function in MDCStoreUtils API:

To initialize the MDCStore interface:

```
BOOL MDCStoreUtils Init()
```

To detach the MDCStore interface before exiting the application:

void WINAPI MDCStoreUtils Finished()

Database Connection Functions

This chapter contains functions that you can use to work with connections to databases.

Table 2-1: Database Connection Functions

Function Name	Description
MDCS_CONNECTION_GetAttributes	To get information about connection attributes
MDCS_CONNECTION_GetDetails	To obtain login information, then return a structure that could be used in all other function calls
MDCS_CONNECTION_TestCredentials	To test user credentials to see if they can create a database connection
MDCS_CONNECTION_CreateConnectionString	To generate ODBC connection string
MDCS_CONNECTION_GetDBHandleFromString	To get the database handle from a connection string
MDCS_CONNECTION_GetDBHandle	To get database handle to establish a connection to a database
MDCS_CONNECTION_GetNewDBHandle	To create a new database handle from an existing one
MDCS_CONNECTION_DestroyDBHandle	To destroy a database handle
MDCS_CONNECTION_Disconnect	To disconnect from the database
MDCS_CONNECTION_Reconnect	To reconnect to the database using an existing connection
MDCS_CONNECTION_GetDetailsFromString	To populate an MDCS_ST_UserLogin structure from an ODBC connection string
MDCS_CONNECTION_BeginTransaction	To start a transaction
MDCS_CONNECTION_CommitTransaction	To commit a transaction
MDCS_CONNECTION_RollbackTransaction	To rollback a transaction
MDCS_CONNECTION_GetDatabaseType	To get the type of database
MDCS_CONNECTION_GetInfo	To get information about a connection: Server name, type, version, ODBC name, driver version, etc.
MDCS_CONNECTION_SetAsyncMode	To set a database connection to async mode, this affects all statements
MDCS_CONNECTION_CancelQueryExecution(To cancel the current execution statement
MDCS_CONNECTION_SetSilentMode	To set silent mode on a connection, no messages will be displayed
MDCS_CONNECTION_CheckIfDead	To check if a connection is active

MDCS_CONNECTION_GetAttributes

```
BOOL MDCS_CONNECTION_GetAttributes(
HDBHANDLE hHandle,
MDCS_ST_ConnectionAttr& stConnectionAttr
);
```

Purpose

To get information about connection attributes; DSN (Data Source) name, User name and database name.

Parameters

hHandle - database connection handle

Output

stConnectionAttr - structure that contains connection attributes.

Return

False - if error occurred

MDCS_CONNECTION_GetDetails

```
BOOL MDCS_CONNECTION_GetDetails(

MDCS_ST_UserLogin * pstLoginInfo,

BOOL* pbSavePasswordChecked = NULL,

HWND hWnd = NULL,

LPCSTR pszTitle = NULL,

const MDCS_ST_UserLogin* pstLoginInfoIn = NULL,

BOOL bDefaultCheckSavePassword = FALSE

BOOL bReloginMode = FALSE

BOOL bValidateVersion = TRUE

);
```

Purpose

To get connection details by calling a dialog to obtain login information and return a structure that could be used in all other function calls. If called for Oracle, then the structure contains the selected database name.

Parameters

```
pstLoginInfoIn - default settings
bCheckSavePassword - default settings for password
pszTitle - dialog title
hWnd - application window
bDefaultCheckSavePassword - to check password by default
bReloginMode - bring up connection dialog in relogin mode
```

22 5000957 C

Output

pstLoginInfo - selected settings
pbSavePasswordChecked - if password save option is checked, if null
the option will not appear in the dialog

Return

FALSE - if dialog is cancelled

MDCS_CONNECTION_TestCredentials

```
BOOL MDCS_CONNECTION_TestCredentials(
const MDCS_ST_UserLogin * pstLoginInfo,
LPSTR pDatabaseName,
int nSize
BOOL bConnectAsAppUser = TRUE
);
```

Purpose

To test if the user can create a database connection with the supplied credentials.

Parameters

pstLoginInfo - pointer to the structure that contains login information
pDataBaseName - pointer to a database name
nSize - size of pDatabaseName
bConectAsAppUser - flag indicates that the connection should be made
using application user permissions

Output

pDatabaseName - name of the database the connection is made to

Return

FALSE - if error occurred

MDCS_CONNECTION_CreateConnectionString

```
BOOL MDCS_CONNECTION_CreateConnectionString(
const MDCS_ST_UserLogin * pstLoginInfo,

LPSTR pszConnectionString,
int nSize
);
```

Purpose

To generate an ODBC connection string from MDCS_ST_UserLogin structure.

Parameters

PstLoginInfo - pointer to the structure that contains login information nSize - size of the connection string

Output

pszConnectionString - generated connection string

Return

FALSE - if error occurred

MDCS_CONNECTION_GetDBHandleFromString

```
HDBHANDLE MDCS_CONNECTION_GetDBHandleFromString(
LPCTSTER pszConnectionString
BOOL SetSilentMode = FALSE
);
```

Purpose

Function to generate a database connection handle from the ODBC connection string.

Parameters

SetSilentMode - flag that indicates if silent mode should be set on connection pszConnectionString - ODBC connection string

Return

Database connection handle

24 5000957 C

MDCS_CONNECTION_GetDBHandle

```
HDBHANDLE MDCS_CONNECTION_GetDBHandle(
const MDCS_ST_UserLogin& stLoginInfo
BOOL bSetSilentMode = FALSE
BOOL dConnectAsAppUser = TRUE
);
```

Purpose

Function to get a database handle.

Parameters

stLoginInfo - pointer to the structure that contains login information bSetSilentMode - flag indicating if silent mode should be set on connection

dConnectAsAppUser - flag indicating that connection should be made using application user permissions

Return

Connection handle

MDCS_CONNECTION_GetNewDBHandle

```
HDBHANDLE MDCS_CONNECTION_GetNewDBHandle(
HDBHANDLE hDBHandle
BOOL bSetSilentMode = FALSE
);
```

Purpose

Function to create a new database handle from an existing one.

Parameters

hDBHandle - existing database connection handle bSetSilentMode - flag indicating if silent mode should be set on connection

Return

Database connection handle

MDCS_CONNECTION_DestroyDBHandle

```
void MDCS_CONNECTION_DestroyDBHandle(
HDBHANDLE hHandle
);
```

Purpose

Function to destroy the database handle.

Parameters

hHandle - database connection handle

MDCS_CONNECTION_Disconnect

```
void MDCS_CONNECTION_Disconnect(
HDBHANDLE hHandle
);
```

Purpose

Function to disconnect from the database.

Parameters

hHandle - database connection handle

MDCS_CONNECTION_Reconnect

```
BOOL MDCS_CONNECTION_Reconnect(
HDBHANDLE hHandle
);
```

Purpose

Function to reconnect to the database using the existing connection; can be used to make sure that the object is connected to the database.

Parameter

hHandle - database connection handle

Return

FALSE - if error occurred

26 5000957 C

MDCS_CONNECTION_GetDetailsFromString

```
BOOL MDCS_CONNECTION_GetDetailsFromString(
MDCS_ST_UserLogin * pstLoginInfo,
LPCSTR pszConnectionString
);
```

Purpose

To generate the MDCS_ST_UserLogin structure from the ODBC connection string.

Parameters

pszConnectionString - ODBC connection string

Output

pstLoginInfo - pointer to the structure that contains login information

Return

FALSE - if error occurred

MDCS_CONNECTION_BeginTransaction

```
BOOL MDCS_CONNECTION_BeginTransaction(
HDBHANDLE hHandle
);
```

Purpose

To start a transaction.

Parameters

hHandle - database connection handle

Return

FALSE - if error occurred

MDCS_CONNECTION_CommitTransaction

```
BOOL MDCS_CONNECTION_CommitTransaction(
HDBHANDLE hHandle
);
```

Purpose

To commit a transaction.

Parameters

hHandle - database connection handle

Return

FALSE - if error occurred

MDCS_CONNECTION_RollbackTransaction

```
BOOL MDCS_CONNECTION_RollbackTransaction(
HDBHANDLE hHandle
);
```

Purpose

To rollback a transaction.

Parameters

hHandle - database connection handle

Return

FALSE - if error occurred

MDCS_CONNECTION_GetDatabaseType

```
MDCS_E_SQLServerType MDCS_CONNECTION_GetDatabaseType(
HDBHANDLE hHandle
);
```

Purpose

To get the type of the database.

Parameters

hHandle - database connection handle

Return

 $MDCS_E_SQLServerType$ - type of database server (if the server type is supported)

28 5000957 C

MDCS_CONNECTION_GetInfo

```
BOOL MDCS_CONNECTION_GetInfo(
HDBHANDLE hHandle,
MDCS_ST_ConnectionInfo& stConnectionInfo
);
```

Purpose

To get information about a connection such as: Server name, type, version, ODBC name, and driver version.

Parameters

hHandle - database connection handle

Output

stConnectionInfo - structure that contains connection information

Return

False - if error occurred

MDCS_CONNECTION_SetAsyncMode

Purpose

To set the database connection to async mode. This affects all statements.

Parameters

hHandle - database connection handle

Output

None

Return

FALSE - if error occurred

MDCS_CONNECTION_CancelQueryExecution(

Purpose

To cancel the current execution statement.

Parameters

hHandle - database connection handle

Output

None

Return

FALSE - if error occurred

MDCS_CONNECTION_SetSilentMode

Purpose

To set silent mode on the connection – no messages will be displayed.

Parameters

hHandle - database connection handle

Output

None

Return

FALSE - if error occurred

30 5000957 C

MDCS_CONNECTION_CheckIfDead

Purpose

To check if the connection is dead.

Parameters

hHandle - database connection handlebReconnect - reconnect if the connection is dead

Output

None

Return

FALSE - if error occurred

Results Set Functions

This chapter contains descriptions of functions you can use to work with Results sets.

Table 3-1: Results Set Functions

Function Name	Description
MDCS_ImportDataToDB	To import result set data into the database
MDCS_GetResultsetData	To get result set data
MDCS_GetAllResultsets	To get a description of all result sets in the database
MDCS_GetResultsetInfo	To get a description of an individual result set

MDCS_ImportDataToDB

```
BOOL MDCS_ImportDataToDB(
HDBHANDLE hHandle,
const MDCS_ST_ResultSet * pst0bjDesc,
const MDCS_ImportResultSet * pExDataSource,
LPSTR psz0bjectID,
int nSize,
MDCS_ProgressCallback * pCallBack
);
```

Purpose

To import the result set data into the database.

Parameters

```
pstObjDesc - result set description
pExDataSource - pointer to the data source that will supply data
hHandle- database connection handle
pCallBack - progress callback object
nSize - size of pszObjectID
```

Output

pszObjectID - ID of the newly created result set

Return

FALSE - if error occurred

MDCS_GetResultsetData

```
BOOL MDCS_GetResultSetData (
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get result set data.

Parameters

hHandle - database connection handlelAssayID - Assay ID that data will be returned for pResultCallback - pointer to a callback function to get data,

Output

Includes the columns

PRINTED_SPOT_ID - feature ID

SUBSTANCE_NAME - feature name

Return

FALSE - if error occurred



Note: Call MDCS_GetResultsetDataTypes on page 37 to find out about the columns that are returned in a callback

34 5000957 C

MDCS_GetAllResultsets

```
BOOL MDCS_GetAllResult sets(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get a description of all result sets in the database.

Parameters

hHandle - database connection handlepResultCallback - pointer to a callback function to get data

Output

Columns from the ASSAYS table

Return

FALSE - if error occurred

MDCS_GetResultsetInfo

```
BOOL MDCS_GetResultsetInfo(
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get a description of individual result set.

Parameters

```
hHandle - database connection handlelAssayID - Assay IDpResultCallback - pointer to a callback function to get data
```

Output

Columns from the ASSAYS table

Return

FALSE - if error occurred

Result Set Data Types Functions

4

This chapter contains functions you can use to work with result set data types. **Table 4-1:** Result Set Data Type Functions

Function Name	Description
MDCS_GetResultsetDataTypes	To get result set data types
MDCS_GetAllDataTypes	To get all data types for all result sets.
MDCS_GetUniqueDataTypes	To get unique data types.
MDCS_GetAllDataTypesOfAssays	To get all data types of a list of assays.
MDCS_GetDataTypesNotInAssaysByDBName	To get all data types that do not belong to the set of assays.

MDCS_GetResultsetDataTypes

```
BOOL MDCS_GetResultsetDataTypes (
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_GetDBResultsCCallback * pResultCallback
):
```

Purpose

To get result set data types.

Parameters

```
hHandle - database connection handleAssayID - Assay IDpResultCallback - pointer to a call back function to get information from database
```

Output

Columns from the TABLE_COLUMNS table

Return

FALSE - if error occurred

MDCS_GetAllDataTypes

```
BOOL MDCS_GetAllDataTypes(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get all data types for all result sets.

Parameters

 $\it hHandle$ - database connection handle $\it pResultCallback$ - pointer to a call back function to get information from database

Output

Columns from the TABLE_COLUMNS table

Return

FALSE - if error occurred

MDCS_GetUniqueDataTypes

```
BOOL MDCS_GetUniqueDataTypes(
HDBHANDLE hHandle,
MDCS GetDBResultsCCallback * pResultCallback);
```

Purpose

To get unique data types.

Parameters

hHandle - database connection handlepResultCallback - pointer to a call back function to get information from database

Output

Columns from the TABLE_COLUMNS table

Return

FALSE - if error occurred

MDCS_GetAllDataTypesOfAssays

Purpose

To get all data types for all assays

Parameters

```
hHandle - database connection handlearrAssayIDs - array of assays IDsnSize - size of array (arrAssayIDs)pResultCallback - pointer to a callback function to get data
```

Output

Columns from the TABLE_COLUMNS table

Return

FALSE - if error occurred

MDCS_GetDataTypesNotInAssaysByDBName

```
BOOL MDCS_GetAllDataTypesOfAssays (
HDBHANDLE hHandle,

LPCSTR pszDBColName,

LONGLONG* arrMicIDs,

INT_PTR nSize,

const MDCS_E_ColumnType& eColType,

const MDCS_E_AverageType& eAvgType,

MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get all data types that don't belong to the provided measurement sets.

Parameters

hHandle - database connection handle
pszDBColName - DB name of data type
arrMicIDs- array of measurement set IDs
nSize - size of array arrMicIDs
eColType - format type of data type such as String or Float...
eAvgType - an average type of data type such as mean, minimal...
pResultCallback - pointer to a callback function to get data

Output

Columns from the TABLE_COLUMNS table COLUMN_NAME, COLUMN_NAME_EXT, COLUMN_TYPE, AVG_TYPE, COLUMN_ID...

Return

FALSE - if error occurred

Error Handling Functions

This chapter contains functions for working with error handling.

Table 5-1: Error Handling Functions

Function Name	Description
MDCS_GetLastError	To get last error and return its error code.
MDCS_GetLastErrorMsg	To get last error message and return error message
MDCS_GetDefaultErrorMsg	To get default error message for the error code

MDCS_GetLastError

```
int MDCS_GetLastError();
```

Purpose

To get the last error code from database.

Return

An error code (for example, MDCS_ERR_SUCCESS).

MDCS_GetLastErrorMsg

```
int MDCS_GetLastErrorMsg(
LPSTR pszError,
int nErrorSize
int* pnErrorSize = NULL
);
```

Purpose

To get the last error message.

Parameters

```
pszError - pointer to an error
nErrorSize - error size, could be as big as MDCS_MAX_ERR_SIZE
pnErrorSize - size of the error not including the null terminated symbol
```

Output

pszError - last error message

Return

An error code (for example, MDCS_ERR_SUCCESS).

MDCS_GetDefaultErrorMsg

```
BOOL MDCS_GetDefaultErrorMsg(
int nError,
LPSTR pszError,
int nSize
);
```

Purpose

Get default error message for the error code MDCS_ERR_XXXXXX.

Parameters

```
nError - error code
nSize - size of the pszError (can use MDCS_MAX_ERR_SIZE)
```

Output

pszError - error message

Return

FALSE - if error is not found

BLOB Functions 6

Table 6-1: BLOB Functions

Function Name	Description
MDCS_BLOB_GetInfo	To get image information
MDCS_BLOB_GetLocationID	To get/create location ID
MDCS_BLOB_GetLocation	To get BLOB location description
MDCS_BLOB_Get	To get BLOB from the storage (database or file server)
MDCS_BLOB_Save	To save BLOB data in database or file server
MDCS_BLOB_RemoveBlobData	To remove BLOB data from database server and file server
MDCS_GetNewDatabaseID	To generate a new database ID
MDCS_BLOB_UpdateBlobDescAndName	To update BLOB's description and name in database and file server
MDCS_BLOB_UpdateBlobDescAndNameEx	To update BLOB's description and name (i.e. attached file) in database
MDCS_BLOB_SaveBlobEx	To save BLOB into database (i.e. attached file to object)
MDCS_BLOB_GetBLOBInfoEx	To get BLOB info using BLOB ID
MDCS_BLOB_GetInfoByReferenceID	To get BLOB info (i.e. image) using reference id
MDCS_BLOB_GetBlobInfoByReferenceIDEx	To get BLOB info (i.e. attached file) using reference id
MDCS_BLOB_GetNumBlobOfObject	To get number of BLOB of an object (i.e. dataset, measurement sets)
MDCS_BLOB_UpdatePlateObjectImageID	To update the image ID in the Plate image table using the BLOB ID
MDCS_BLOB_GetUniqueLocationNamebyFilter	To get a unique location name based on a filter
MDCS_BLOB_GetAllBlobInfoByLocationFilterAnd LocationInfo	To get all image information based on criteria, BLOB type and file location filter
MDCS_BLOB_GetAllBlobInfoByPlateID	To get all image information based on the plate ID
MDCS_BLOB_DeleteUnUsedLocations	To remove unused locations in a File_location table
MDCS_BLOB_DeleteFileLocation	To remove a file location
MDCS_BLOB_LocationIsUsed	To check if any existing BLOB are stored in a location

Table 6-1: BLOB Functions (cont'd)

MDCS_BLOB_RemoveBlobData	To remove BLOB data from datebase server and file server
MDCS_BLOB_RemoveFile	To remove file on file server or UNC path location
MDCS_BLOB_Attach	To attach a BLOB
MDCS_BLOB_UpdateBlobInfo	To update the description and display name of the attached file
MDCS_BLOB_MapLocations	To map BLOB locations
MDCS_BLOB_CanWriteToLocation	To check if the BLOB can be saved to a specified location
MDCS_BLOB_FillOutBlobLocationStruct	To fill out MDCS_ST_BlobLocation structure with results of MDCS_QueryResults

MDCS_BLOB_GetInfo

```
BOOL MDCS_BLOB_GetInfo(
HDBHANDLE hHandle,

MDCS_ST_BlobInfo * pBlobInfo,

LONGLONG lBlobID,

MDCS_E_BlobType eBlobType
);
```

Purpose

To get image information using BLOB ID.

Parameters

hHandle - database connection handle BlobID - database ID of the BLOB eBlobType - BLOB Type

Output

pBlobInfo - BLOB information

Return

FALSE - if error occurred

MDCS_BLOB_GetLocationID

```
BOOL MDCS_BLOB_GetLocationID (
HDBHANDLE hHandle,
const MDCS_ST_BlobLocation * pstLocation,
LONGLONG* lLocationID,
BOOL bCreateNew = FALSE
);
```

Purpose

To get/create a location ID.

Parameters

hHandle - database connection handle
 bCreateNew - if TRUE - function creates a new record
 pstLocation - pointer to the structure that contains location description
 (all fields must be filled except for nLocationID)

Output

nLocationID - location ID, -1 if location is not found

Return

FALSE - if function fails

MDCS_BLOB_GetLocation

```
BOOL MDCS_BLOB_GetLocation (
HDBHANDLE hHandle,
LONGLONG nLocationID,
MDCS_ST_BlobLocation * pstLocation
);
```

Purpose

To get the BLOB location description.

Parameters

hHandle - database connection handlenLocationID - location ID

Output

pstLocation - pointer to the structure that contains location description

Return

FALSE - if error occurred

MDCS_BLOB_Get

```
BOOL MDCS_BLOB_Get(
HDBHANDLE hHandle,
LONGLONG lBlobID,
MDCS_E_BlobType eBlobType,
MDCS_GetBlobCallback * pCallBack
);
```

Purpose

To get BLOB from storage, database or file server.

Parameters

```
    hHandle - database connection handle
    lBlobID - database ID of the BLOB
    eBlobType - BLOB Type (from MDCS_E_BlobType)
    pCallBack - pointer to callback used to transfer data
```

Return

FALSE - if error occurred

MDCS_BLOB_Save

```
BOOL MDCS_BLOB_Save(
HDBHANDLE hHandle,
const MDCS_ST_BlobInfo * pImageInfo,

MDCS_E_BlobType eBlobType,

MDCS_SaveBlobCallback * pCallBack,

LONGLONG* 1BlobIDOut,
const LONGLONG* 1RefObjectID,

MDCS_E_FileStorage eStorage = MDCS_eDatabaseServer,

LPCSTR pszDirectoryExtra = NULL
);
```

Purpose

To save BLOB data in the database or file server.

Parameters

```
    hHandle - database connection handle
    pBlobInfo - pointer to the structure that contains image information
    pCallBack - pointer to callback used to transfer data
    eStorage - type of destination where image will be saved (default is database)
    eBlobType - type of the BLOB
    pszDirectoryExtra - if specified, creates a directory with
    pszDirectoryExtra on FileServer
```

Output

```
IBlobIDOut - database ID of the BLOB
IRefObjectID - object reference ID
```

Return

FALSE - if error occurred

MDCS_BLOB_RemoveBlobData

```
BOOL MDCS_BLOB_RemoveBlobData(
HDBHANDLE hHandle,
LONGLONG lID,
MDCS_E_BlobType eBlobType
);
```

Purpose

To remove BLOB data from database server and file server.

Parameters

```
hHandle - database connection handleIID - ID of the record that needs to be deletedeBlobType - type of the BLOB
```

Return

FALSE - if error occurred

MDCS_GetNewDatabaseID

```
BOOL MDCS_GetNewDatabaseID(
HDBHANDLE hHandle,

LPSTR pszID,
int nIDSize,
int* pnIDSizeOut
);
```

Purpose

To generate a new database ID.

Parameters

```
hHandle - database connection handle
nIDSize - size of ID (use MDCS_MAX_DATABASE_ID)
```

Output

```
pszID - generated ID
pnIDSizeOut - size of the generated ID
```

Return

FALSE - if error occurred

MDCS_BLOB_UpdateBlobDescAndName

```
BOOL MDCS_BLOB_UpdateBlobDescAndName(
HDBHANDLE hHandle,
LONGLONG lBlobID,
Const MDCS_ST_BlobInfo* pBlobInfo,
MDCS_E_BlobType eBlobType
);
```

Purpose

To update a BLOB's description and name.

Parameters

```
    hHandle - database connection handle
    lBlobID - Blob ID in database
    pBlobInfo - Blob info to update in the database (only need description and name fields)
    eBlobType - Blob type
```

Return

FALSE - if error occurred

48 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

MDCS_BLOB_UpdateBlobDescAndNameEx

```
BOOL MDCS_BLOB_UpdateBlobDescAndName(
HDBHANDLE hHandle,
LONGLONG lBlobID,
Const MDCS_ST_BlobInfoEx* pBlobInfo,
MDCS_E_BlobType eBlobType
);
```

Purpose

Function to update a BLOB's description and name.

Parameters

```
    hHandle - database connection handle
    lBlobID - Blob ID in database
    pBlobInfo - Blob info to update in database (only need description and name fields)
    eBlobType - Blob type
```

Return

FALSE - if error occurred

MDCS_BLOB_SaveBlobEx

```
BOOL MDCS_BLOB_SaveBlobEx(
HDBHANDLE hHandle,
Const MDCS_ST_BlobInfoEx* pBlobInfo,
MDCS_E_BlobType eBlobType,
MDCS_SaveBlobCallback* pCallBack,
LONGLONG* lBlobIDOut,
Const LONGLONG* lRefObjectID
);
```

Purpose

To save a BLOB into the database or file server.

Parameters

```
    hHandle - database connection handle
    pBlobInfo - Blob info (ID not required, since this is input)
    eBlobType - Blob type
    pCallBack - callback class to save BLOB
    IRefObjectID - reference object id that the BLOB associated with
```

OUTPUT

IBlobIDOut - output id of saved BLOB

Return

FALSE - if error occurred

MDCS_BLOB_GetBLOBInfoEx

```
BOOL MDCS_BLOB_GetBLOBInfoEx(
HDBHANDLE hHandle,
MDCS_ST_BlobInfoEx* pBlobInfo,
LONGLONG* lBlobID,
MDCS_E_BlobType eBlobType
);
```

Purpose

To get attached file or BLOB information based on provided BLOB ID.

Parameters

```
hHandle - database connection handleeBlobType - BLOB typelBlobID - BLOB ID
```

OUTPUT

pBlobInfo - BLOB info

Return

FALSE - if error occurred

MDCS_BLOB_GetInfoByReferenceID

```
BOOL MDCS_BLOB_ GetInfoByReferenceID (
HDBHANDLE hHandle,
LONGLONG lReferenceID,
MDCS_E_BlobType eBlobType,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get BLOB information based on provided reference id (dataset ID/ assay ID).

Parameters

```
    hHandle - database connection handle
    lReferenceID - reference id such as dataset id or measurement set id
    eBlobType - BLOB type
    pResultCallback - pointer to a callback function to get data
```

OUTPUT

Columns for BLOB Info like ID, SIZE, DESC etc.

Return

FALSE - if error occurred

MDCS_BLOB_GetBlobInfoByReferenceIDEx

```
BOOL MDCS_BLOB_ GetInfoByReferenceIDEx(
HDBHANDLE hHandle,
LONGLONG lReferenceID,
MDCS_E_BlobType eBlobType,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get BLOB information from the database using BLOB id and type.

Parameters

hHandle - database connection handle
 lReferenceID - reference id such as dataset id or measurement set id
 eBlobType - BLOB type
 pResultCallback - pointer to a callback function to get data

Output

Columns for BLOB info like ID, SIZE, DESC etc.

Return

FALSE - if error occurred

MDCS_BLOB_GetNumBlobOfObject

```
BOOL MDCS_BLOB_GetNumBlobOfObject(
HDBHANDLE hHandle,
LONGLONG lReferenceID,
MDCS_E_BlobType eBlobType,
LONGLONG& lNumbBlob
);
```

Purpose

Function to get the number of BLOB files attached to an object (assay, dataset) using the reference ID.

Parameters

```
hHandle - database connection handleIReferenceID - reference id such as dataset id or assay ideBlobType - BLOB type
```

OUTPUT

INumbBlob – number of BLOB files associated with the object (dataset, measurement set)

Return

FALSE - if error occurred

MDCS_BLOB_UpdatePlateObjectImageID

```
BOOL MDCS_BLOB_UpdatePlateObjectImageID(
HDBHANDLE hHandle,
LONGLONG lCurrentImageObjID,
LONGLONG lNewImageObjID,
);
```

Purpose

To update the image ID in the PLATE_IMAGES table using the BLOB ID.

Parameters

hHandle - database connection handlelCurrentImageObjID - BLOB IDlNewImageObjID - Updated value of BLOB ID in PLATE_IMAGES

Return

FALSE - if error occurred

MDCS_BLOB_GetUniqueLocationNamebyFilter

```
BOOL MDCS_BLOB_GetUniqueLocationNamebyFilter (
HDBHANDLE hHandle,
MDCS_CL_BlobLocationCB* pCallBack
MDCS_E_LocationFilter eFilter = MDCS_eServerName
);
```

Purpose

To get unique location information based on a filter.

Parameters

eFilter - filter to get data hHandle - database connection handle pCallBack - callback object

Output

Depending on the filter, the output result set will contain one, some, or all (1-5) of the following:

- 1) server name
- 2) serverRoot
- 3) Directory
- 4) Location type
- 5) port number

Return

FALSE - if function fails

MDCS_BLOB_GetAllBlobInfoByLocationFilterAndLocationInfo

```
BOOL
MDCS_BLOB_GetAllBlobInfoByLocationFilterAndLocationInfo (
HDBHANDLE hHandle,
MDCS_ST_BlobLocation* arrBlobLocation,
INT_PRT NNumLocation,
MDCS_E_BlobType BlobType
MDCS_CL_BlobInfoCB* PCallBack
MDCS_E_LocationFilter eFilter = MDCS_eServerName
);
```

Purpose

To get all image information based on criteria, BLOB type and file location filter (server name, root name, or directory).

Parameters

```
eFilter - filter to get data

hHandle - database connection handle

eBlobType - Type of BLOB

arrBlobLocation - all BLOB Location information

NNumLocation - size of arrBlobLocation array
```

Output

pCallback - call back result

Return

FALSE - if error occurred

MDCS_BLOB_GetAllBlobInfoByPlateID

```
BOOL MDCS_BLOB_GetAllBlobInfoByPlateID (
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_E_BlobType eBlobType
MDCS_CL_BlobInfoCB* PCallBack
);
```

Purpose

To get all image information based on the plate ID.

Parameters

```
hHandle - database connection handleIPlateID - plate IDeBlobType - Type of BLOB
```

Output

pCallback - call back result

Return

FALSE - if error occurred

MDCS_BLOB_DeleteUnUsedLocations

```
BOOL MDCS_BLOB_DeleteUnUsedLocations(
HDBHANDLE hHandle,
MDCS_ProgressCallBack* pProgressCallback
);
```

Purpose

To remove unused locations in FILE_LOCATIONS.

Parameters

hHandle - database connection handle
pProgressCallback - progress callback

Return

FALSE - if error occurred

MDCS_BLOB_DeleteFileLocation

```
BOOL MDCS_BLOB_DeleteFileLocation(
HDBHANDLE hHandle,
const LONGLONG lLocationID
);
```

Purpose

To remove a file location.

Parameters

hHandle - database connection handlelLocationID - Location ID

Return

FALSE - if error occurred

MDCS_BLOB_LocationIsUsed

```
BOOL MDCS_BLOB_LocationIsUsed (
HDBHANDLE hHandle,
LONGLONG lLocationID,
BOOL& bIsUsed
);
```

Purpose

To check if any existing BLOBs are stored in the location.

Parameters

hHandle - database connection handle
lLocationID - location ID

Output

bIsUsed - TRUE if the Location ID is used by a BLOB, otherwise FALSE

Return

FALSE - if error occurred

MDCS_BLOB_RemoveBlobData

```
BOOL MDCS_BLOB_RemoveBlobData (
HDBHANDLE hHandle,
LONGLONG lID,
MDCS_E_BlobType eBlobType
BOOL bNotRemoveRecord = FALSE
);
```

Purpose

To remove BLOB data from a database, file server or UNC path location.

Parameters

```
hHandle - database connection handleIID - ID of the record that needs to be deletedeBlobType - Type of BLOBbNotRemoveRecord - Whether to delete the image record or not
```

Return

FALSE - if error occurred

MDCS_BLOB_RemoveFile

```
BOOL MDCS_BLOB_RemoveFile (
HDBHANDLE hHandle,
const MDCS_ST_BlobLocation& stBlobLocation
LPCSTR szFileName
);
```

Purpose

To remove a file from the image storage location (database, file server or UNC path).

Parameters

```
hHandle - database connection handle
stBlobLocation - structure describing the location of the file to be
deleted
eBlobType - Type of BLOB
```

Return

FALSE - if error occurred

MDCS_BLOB_Attach

```
BOOL MDCS_BLOB_Attach (
MDCS_E_BlobType eBlobType
const MDCS_ST_BlobinfoEX& stBlobInfo,
LONGLONG* plIDOut = NULL
);
```

Purpose

To attach a BLOB (change BLOB information).

Parameters

```
hHandle - database connection handlestBlobInfo - structure that describes BLOB.eBlobType - Type of BLOB
```

Output

plIDOut - created image ID

Return

FALSE - if error occurred

MDCS_BLOB_UpdateBlobInfo

```
BOOL MDCS_BLOB_UpdateBlobInfo (
HDBHANDLE hHandle,
LONGLONG lBlobID
const MDCS_ST_BlobInfo* pBlobInfo,
MDCS_E_BlobType eBlobType
);
```

Purpose

To update the description and display name of the BLOB.

Parameters

```
hHandle - database connection handlelBlobID - BLOB IDpBlobInfo - information to update in databaseeBlobType - Type of BLOB
```

Return

FALSE - if error occurred

MDCS_BLOB_MapLocations

```
BOOL MDCS_BLOB_MapLocations (
HDBHANDLE hHandle,
const MDCS_E_BlobLocation& stSource
const MDCS_ST_BlobLocation& stDest
);
```

Purpose

To map BLOB locations.

Parameters

hHandle - database connection handle stSource - source BLOB location stDest - destination BLOB location

Return

FALSE - if error occurred

MDCS_BLOB_CanWriteToLocation

```
BOOL MDCS_BLOB_CanWriteToLocation(
HDBHANDLE hHandle,
const MDCS_STBlobLocation& stBlobLocation
LPCSTR pszDirectoryExtra = NULL
);
```

Purpose

To check if the BLOB can be saved to a specified location.

Parameters

hHandle - database connection handlestBlobLocation - structure that describes the locationpszDirectoryExtra - if specified, will create a directory with pszDirectoryExtra on File server

Return

FALSE - if error occurred

MDCS_BLOB_FillOutBlobLocationStruct

```
BOOL MDCS_BLOB_FillOutBlobLocationStruct(
const MDCS_QueryResults* pResults,

MDCS_ST_BlobLocation& stBlobLoc,

BOOL bvNativeFormat = TRUE

);
```

Purpose

To fill out MDCS_ST_BlobLocation structure with results of MDCS_QueryResults.

Parameters

pResults - results of the query bvNativeFormat - format type of query results, if not native then all results are treated as strings

Output

stBlobLoc - location structure (if stBlobLoc.lLocationID == 0 - the BLOB is in the database)

Return

FALSE - if error occurred

Measurement Sets and Results Functions

This chapter contains functions that you can use to work with measurement sets and results of measurement sets.

Table 7-1: Measurement Sets and Results of Measurement Sets Functions

Function Name	Description
MDCS_ASSAY_AppendMeasurementSet	To import (append) measurement set data into the database
MDCS_ASSAY_ImportMeasurementSet	To remove shape data from a measurement set
MDCS_ASSAY_ImportMeasurementSet	To import measurement set data into the database
MDCS_ASSAY_InsertDValue	To insert value of double type to the table that contains results of measurement set
MDCS_ASSAY_InsertNumericValue	To insert value of double type to the table that contains results of measurement set
MDCS_ASSAY_InsertSValue	To insert value of string to the table that contains results of measurement set
MDCS_ASSAY_InsertStringValue	To insert value of string type to the table that contains results of measurement set
MDCS_ASSAY_UpdateShapeLines	To insert or update shapes
MDCS_ASSAY_GetShapeLineBlobDesc	To get shape lines blob description
MDCS_ASSAY_GetShapeLineBlobDescBySite	To get shape lines blob description per site
MDCS_ASSAY_GetShapeLineBlobDescBySiteAndSeries	To get shape lines blob description per site and series
MDCS_ASSAY_GetShapeLineBlobDescBySiteSeriesAnd Instance	To get shape lines blob description per site, series and instance
MDCS_ASSAY_GetShapeLineBlobDescByPlate	To get shape lines blob description by plate
MDCS_ASSAY_GetShapeLineBlobPerAssay	To get all the shape lines blob descriptions per measurement set
MDCS_ASSAY_GetShapeLines	To get shape lines of a shape
MDCS_ASSAY_GetMeasurementAttributes	To get a measurement description based on measurement set ID and measurement ID
MDCS_ASSAY_GetAssayMeasurementRecord	To get a measurement record

Table 7-1: Measurement Sets and Results of Measurement Sets Functions (cont'd)

MDCS_ASSAY_GetAssayMeasurementsByPlateAndMeasurement	To get measurement set measurement by plate ID and measurement ID
MDCS_ASSAY_GetMeasurementsBySiteID	To get measurements by Site ID
MDCS_ASSAY_GetMeasurementInfoByAssayAndColumnName	To get a measurement description based on measurement set ID and database column name
MDCS_ASSAY_GetMeasurementByFunctionAndParame terName	To get a measurement attribute description based on measurement and functions names
MDCS_ASSAY_GetMeasurementRecord	To get a data type record by its ID
MDCS_ASSAY_Create	To create a new measurement set in a specific folder
MDCS_ASSAY_GetAssaySiblingFolders	To get measurement set sibling folders
MDCS_ASSAY_CreateMeasurement	To create a new measurement set measurement
MDCS_ASSAY_AddMeasurement	To add a new measurement to a measurement set
MDCS_ASSAY_GetMeasurementByID	To get the description of a measurement set measurement
MDCS_ASSAY_GetMeasurementByName	To get the description of a measurement set measurement
MDCS_ASSAY_Delete	To delete a measurement set
MDCS_ASSAY_DeleteAllForPlate	To delete all measurement sets for a plate
MDCS_ASSAY_DeleteDataForPlate	To delete measurement set results for a plate
MDCS_ASSAY_DeleteFolder	To delete a measurement set folder
MDCS_ASSAY_CreateOutlinesTable	To create a measurement set outlines table
MDCS_ASSAY_ManageFolderSecurity	Call a dialog to manage security access to the measurement set folders
MDCS_ASSAY_CanModifyFolder	To check if a user can modify a folder
MDCS_ASSAY_CreateFolder	Create a new measurement set folder in the database
MDCS_ASSAY_ModifyFolder	Modify a measurement set folder in the database
MDCS_ASSAY_DoesSubFolderExist	To check if the sub folder exists

Table 7-1: Measurement Sets and Results of Measurement Sets Functions (cont'd)

MDCS_ASSAY_GetAllInFolder	To get measurement sets with a basic description in folder
MDCS_ASSAY_GetFolderPath	To get the path of a folder
MDCS_ASSAY_Reindex	To reindex measurement set measurements
MDCS_ASSAY_ReindexShapes	To reindex measurement set cell shapes
MDCS_ASSAY_OptimizeAll	To optimize all measurement set tables
MDCS_ASSAY_Copy	To copy a measurement set
MDCS_ASSAY_Merge	To merge a measurement set
MDCS_ASSAY_GetSiteCount	To count site appearance in measurement set table
MDCS_ASSAY_GetSiteInfoImageByID	To get site image information by provided ID
MDCS_ASSAY_GetSiteInfoImageBySiteAndSeriesID	To get site image info by site and series ID
MDCS_ASSAY_MeasurementGetSiteCount	To count site appearance in a measurement set
MDCS_ASSAY_DeleteMeasurement	To delete a measurement
MDCS_ASSAY_CellOutlinesGetSiteCount	To count site appearance in a assay cell outlines
MDCS_ASSAY_CreateNewName	Create a new measurement set name in the destination folder.
MDCS_ASSAY_CreateRun	To create a new measurement set run record
MDCS_ASSAY_GetLatestAssayRunID	To get the ID of the latest run on a measurement set
MDCS_ASSAY_GetByID	Get measurement set data
MDCS_ASSAY_CreateProfile	Create a new measurement set profile record
MDCS_ASSAY_AssociateWithPlate	To associate a measurement set with a plate
MDCS_ASSAY_GetSpotID	To get the spot ID from instance site and series
MDCS_ASSAY_GetAllForPlate	To get all measurement sets that are associated with a plate
MDCS_ASSAY_GetProfiles	To get all available profiles
MDCS_ASSAY_GetProfile	To get a profile record
MDCS_ASSAY_GetProfileInfo	To get profile information
MDCS_ASSAY_UpdateProfile	To update a profile record
MDCS_ASSAY_DeleteProfile	To delete a measurement set profile

Table 7-1: Measurement Sets and Results of Measurement Sets Functions (cont'd)

MDCS ASSAY GetRecord	To get a measurement set
	record
MDCS_ASSAY_GetAssayByRunID	To get a measurement set record using RUN ID
MDCS_ASSAY_UsedInDatasets	To check if a measurement set is used in datasets
MDCS_ASSAY_HavePermissionsToModify	To check if the current user can modify a measurement set
MDCS_ASSAY_GetImageSourceRecords	To get all records from the IMAGE_SOURCE table
MDCS_ASSAY_GetStatisticalValuesForPlateAndAssay	To get statistical values for a selected plate and assay
MDCS_ASSAY_GetStatisticalValuesForPlateAndMeasur ement	To get statistical values for a selected plate and measurement
MDCS_ASSAY_GetValuesForPlateAndAssay	To get values for the selected plate and assay
MDCS_ASSAY_GetValuesForPlate	To get values for the selected plate
MDCS_ASSAY_GetAllSiteMeasurements	To get all values for the selected site
MDCS_ASSAY_GetMeasurementInfoByAssayAndPlate	To get data types records by plate and assay
MDCS_ASSAY_GetAssayInfoByNamePlateAndSettings	To get assay information by plate name and plate setting
MDCS_ASSAY_GetAssayIDsOfPlate	To get assay IDs that are available for a plate
MDCS_ASSAY_GetAssaySettingsOfPlate	To get unique assay settings of the plate
MDCS_ASSAY_GetMeasurementStatistic	To get statistic use of a measurement (number of assays and datasets that use a measurement set)
MDCS_ASSAY_GetShapeLinesBySite	To get shapeline blob description per site.
MDCS_ASSAY_GetShapeLinesBySiteAndSeries	To get shapeline blob description per site and series.
MDCS_ASSAY_GetMarkedAssaysWithCallback	To get assays marked for deletion
MDCS_ASSAY_Restore	To restore the deleted assay
MDCS_ASSAY_CreateForPlate	To create a new assay for a specific plate.
MDCS_ASSAY_GetMeasurementByName	To get the description of an assay measurement using display name
MDCS_ASSAY_GetMeasurementByDBName	To get the description of an assay measurement using DB column name

64 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Table 7-1: Measurement Sets and Results of Measurement Sets Functions (cont'd)

MDCS_ASSAY_GetScopeAttributeByID To get the description of a scope attribute using DB column name MDCS_ASSAY_DeleteMeasurementInAssays To delete a measurement from all assays that use it but not the default MDCS_ASSAY_UpdateDataType To update a measurement property MDCS_ASSAY_UpdateDataTypeByAssay To modify a measurement that is not a default MDCS_ASSAY_CanModifyAssay To check if the current user can modify an assay id MDCS_ASSAY_CalculateStatisticResults To get statistic results on an array of assay ids MDCS_ASSAY_CalculateStatisticEx To get statistic results on an array of assay ids MDCS_ASSAY_GetZPrime To calculate Z' on assay and measurement type MDCS_ASSAY_GetZPrimeScopeAttribute To calculate Z' on assay and scope attribute MDCS_ASSAY_GetUniqueMeasurementValues To get unique values of measurement type for array of mensuration sets MDCS_ASSAY_FindMeasurementValues To find values from a provided array in a measurement set MDCS_ASSAY_GetAllMSetUniqueAnnotation To get unique annotation of all assays MDCS_ASSAY_GetAttributeInfoByDisplayName To get attribute infor using its display name. MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name. MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay	MDCS_ASSAY_GetScopeAttributeByName	To get the description of a scope attribute using attribute display name
from all assays that use it but not the default MDCS_ASSAY_UpdateDataType MDCS_ASSAY_UpdateDataTypeByAssay To modify a measurement that is not a default measurement for an assay based on DB column name and assay id MDCS_ASSAY_CanModifyAssay To check if the current user can modify an assay with the current user can modify an assay id MDCS_ASSAY_CalculateStatisticResults To get statistic results on an array of assay ids MDCS_ASSAY_CalculateStatisticEx To get statistic results on an array of assay ids MDCS_ASSAY_GetZPrime To calculate Z' on assay and measurement type MDCS_ASSAY_GetZPrimeScopeAttribute To calculate Z' on assay and scope attribute MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_FindMeasurementValues To get unique values from a provided array in a measurement set MDCS_ASSAY_UpdateMeasurementData To update measurement data based on transformation criteria MDCS_ASSAY_GetAltributeUniqueAnnotation To get unique annotation of all assays MDCS_ASSAY_GetAltributeInfoByDisplayName To get attribute inforusing its display name. MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_GetAttributeValueByDisplayName To assign an Assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueCong To assign an Assay attribute value (type string)	MDCS_ASSAY_GetScopeAttributeByID	scope attribute using DB
mDCS_ASSAY_UpdateDataTypeByAssay To modify a measurement that is not a default measurement for an assay based on DB column name and assay id mDCS_ASSAY_CanModifyAssay To check if the current user can modify an assay with the can modify an assay id modify an assay id modify an assay id modify an assay ids mDCS_ASSAY_CalculateStatisticEx To get statistic results on an array of assay ids mDCS_ASSAY_GetZPrime To calculate 2' on assay and measurement type mDCS_ASSAY_GetZPrimeScopeAttribute mDCS_ASSAY_UpdateMeasurementValues mDCS_ASSAY_UpdateMeasurementValues mDCS_ASSAY_UpdateMeasurementData mDCS_ASSAY_UpdateMeasurementData mDCS_ASSAY_GetAllMSetUniqueAnnotation mDCS_ASSAY_GetAllMSetUniqueAnnotation mDCS_ASSAY_GetAllMSetUniqueAnnotation mDCS_ASSAY_GetAltributeInfoByDisplayName mDCS_ASSAY_GetAttributeInfoByDisplayName mDCS_ASSAY_GetAttributeInfoByDBName mDCS_ASSAY_GetAttributeValueByDBName mDCS_ASSAY_GetAttributeValueByDBName mDCS_ASSAY_GetAttributeValueByDisplayName mDCS_ASSAY_AssignAttributeValueString mDCS_ASSAY_AssignAttributeValueString mDCS_ASSAY_AssignAttributeValueLong mDCS_ASSAY_AssignAttrib	MDCS_ASSAY_DeleteMeasurementInAssays	from all assays that use it but
that is not a default measurement for an assay based on DB column name and assay id MDCS_ASSAY_CanModifyAssay To check if the current user can modify an assay MDCS_ASSAY_CalculateStatisticResults To get statistic results on an array of assay ids MDCS_ASSAY_CalculateStatisticEx To get statistic results on an array of assay ids MDCS_ASSAY_GetZPrime To calculate Z' on assay and measurement type MDCS_ASSAY_GetZPrimeScopeAttribute MDCS_ASSAY_GetZPrimeScopeAttribute MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_UpdateMeasurementData To find values from a provided array in a measurement set MDCS_ASSAY_UpdateMeasurementData To update measurement data based on transformation criteria MDCS_ASSAY_GetAllMSetUniqueAnnotation MDCS_ASSAY_GetAllMSetUniqueAnnotation To get unique annotation of all assays MDCS_ASSAY_GetAttribute MDCS_ASSAY_GetAttributeInfoByDisplayName MDCS_ASSAY_GetAttributeInfoByDBName To get attribute information using its database name MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its database column name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute value (type string)	MDCS_ASSAY_UpdateDataType	To update a measurement property
can modify an assay MDCS_ASSAY_CalculateStatisticResults To get statistic results on an array of assay ids MDCS_ASSAY_CalculateStatisticEx To get statistic results on an array of assay ids MDCS_ASSAY_GetZPrime To calculate Z' on assay and measurement type MDCS_ASSAY_GetZPrimeScopeAttribute To get unique values of measurement type for array of mensuration sets MDCS_ASSAY_FindMeasurementValues To find values from a provided array in a measurement set MDCS_ASSAY_UpdateMeasurementData To update measurement data based on transformation criteria MDCS_ASSAY_GetAllMSetUniqueAnnotation MDCS_ASSAY_GetAllMSetUniqueAnnotation To get unique annotation of all assays MDCS_ASSAY_GetAttributeInfoByDisplayName To get attribute info using its display name. MDCS_ASSAY_GetAttributeInfoByDBName MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name. MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name. MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name. MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string)	MDCS_ASSAY_UpdateDataTypeByAssay	that is not a default measurement for an assay based on DB column name
mDCS_ASSAY_GetZPrime MDCS_ASSAY_GetZPrime MDCS_ASSAY_GetZPrimeScopeAttribute MDCS_ASSAY_GetZPrimeScopeAttribute MDCS_ASSAY_GetZPrimeScopeAttribute MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_UpdateMeasurementData MDCS_ASSAY_UpdateMeasurementData MDCS_ASSAY_GetAllMSetUniqueAnnotation MDCS_ASSAY_GetAllMSetUniqueAnnotation MDCS_ASSAY_GetAltribute MDCS_ASSAY_GetAttributeInfoByDisplayName MDCS_ASSAY_GetAttributeInfoByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database name To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its database column name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute value (type string) To assign an Assay attribute value (type string)	MDCS_ASSAY_CanModifyAssay	
array of assay ids MDCS_ASSAY_GetZPrime To calculate Z' on assay and measurement type MDCS_ASSAY_GetZPrimeScopeAttribute To calculate Z' on assay and scope attribute MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_GetUniqueMeasurementValues To get unique values of measurement type for array of mensuration sets MDCS_ASSAY_FindMeasurementValues To find values from a provided array in a measurement set MDCS_ASSAY_UpdateMeasurementData To update measurement data based on transformation criteria MDCS_ASSAY_GetAllMSetUniqueAnnotation To get unique annotation of all assays MDCS_ASSAY_CreateAttribute MDCS_ASSAY_GetAttributeInfoByDisplayName To get attribute info using its display name. MDCS_ASSAY_GetAttributeInfoByDBName MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name To get the assay attribute value using its display name To get the assay attribute value using its display name To get the assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string)	MDCS_ASSAY_CalculateStatisticResults	To get statistic results on an array of assay ids
MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_UpdateMeasurementData MDCS_ASSAY_UpdateMeasurementData MDCS_ASSAY_GetAllMSetUniqueAnnotation MDCS_ASSAY_GetAllMSetUniqueAnnotation MDCS_ASSAY_CreateAttribute MDCS_ASSAY_GetAttributeInfoByDisplayName MDCS_ASSAY_GetAttributeInfoByDBName MDCS_ASSAY_GetAttributeInfoByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_AssignAttributeValueString MDCS_ASSAY_AssignAttributeValueString MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute value (type string) To assign an Assay attribute	MDCS_ASSAY_CalculateStatisticEx	To get statistic results on an array of assay ids
MDCS_ASSAY_GetUniqueMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_FindMeasurementValues MDCS_ASSAY_UpdateMeasurementData MDCS_ASSAY_UpdateMeasurementData MDCS_ASSAY_GetAllMSetUniqueAnnotation MDCS_ASSAY_GetAllMSetUniqueAnnotation MDCS_ASSAY_CreateAttribute MDCS_ASSAY_GetAttributeInfoByDisplayName MDCS_ASSAY_GetAttributeInfoByDBName MDCS_ASSAY_GetAttributeInfoByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_AssignAttributeValueString MDCS_ASSAY_AssignAttributeValueLong MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_GetZPrime	
measurement type for array of mensuration sets MDCS_ASSAY_FindMeasurementValues To find values from a provided array in a measurement set MDCS_ASSAY_UpdateMeasurementData To update measurement data based on transformation criteria MDCS_ASSAY_GetAllMSetUniqueAnnotation To get unique annotation of all assays MDCS_ASSAY_CreateAttribute MDCS_ASSAY_GetAttributeInfoByDisplayName To get attribute info using its display name. MDCS_ASSAY_GetAttributeInfoByDBName To get attribute information using its database name MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_GetZPrimeScopeAttribute	To calculate Z' on assay and scope attribute
mDCS_ASSAY_UpdateMeasurementData To update measurement data based on transformation criteria MDCS_ASSAY_GetAllMSetUniqueAnnotation To get unique annotation of all assays MDCS_ASSAY_CreateAttribute MDCS_ASSAY_GetAttributeInfoByDisplayName To get attribute info using its display name. MDCS_ASSAY_GetAttributeInfoByDBName To get attribute information using its database name MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_GetUniqueMeasurementValues	measurement type for array
based on transformation criteria MDCS_ASSAY_GetAllMSetUniqueAnnotation To get unique annotation of all assays MDCS_ASSAY_CreateAttribute To create an Assay attribute MDCS_ASSAY_GetAttributeInfoByDisplayName To get attribute info using its display name. MDCS_ASSAY_GetAttributeInfoByDBName To get attribute information using its database name MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_FindMeasurementValues	To find values from a provided array in a measurement set
mDCS_ASSAY_CreateAttribute MDCS_ASSAY_GetAttributeInfoByDisplayName MDCS_ASSAY_GetAttributeInfoByDBName To get attribute information using its database name MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_UpdateMeasurementData	based on transformation
MDCS_ASSAY_GetAttributeInfoByDisplayName To get attribute info using its display name. MDCS_ASSAY_GetAttributeInfoByDBName To get attribute information using its database name MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_GetAllMSetUniqueAnnotation	
display name. MDCS_ASSAY_GetAttributeInfoByDBName To get attribute information using its database name MDCS_ASSAY_GetAttributeValueByDBName To get the assay attribute value using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_CreateAttribute	To create an Assay attribute
 MDCS_ASSAY_GetAttributeValueByDBName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_GetAttributeValueByDisplayName MDCS_ASSAY_AssignAttributeValueString MDCS_ASSAY_AssignAttributeValueLong MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute 	MDCS_ASSAY_GetAttributeInfoByDisplayName	
walue using its database column name MDCS_ASSAY_GetAttributeValueByDisplayName To get the assay attribute value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_GetAttributeInfoByDBName	To get attribute information using its database name
value using its display name MDCS_ASSAY_AssignAttributeValueString To assign an Assay attribute value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_GetAttributeValueByDBName	value using its database
value (type string) MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute	MDCS_ASSAY_GetAttributeValueByDisplayName	To get the assay attribute value using its display name
MDCS_ASSAY_AssignAttributeValueLong To assign an Assay attribute value (type Long)	MDCS_ASSAY_AssignAttributeValueString	To assign an Assay attribute value (type string)
	MDCS_ASSAY_AssignAttributeValueLong	To assign an Assay attribute value (type Long)

Table 7-1: Measurement Sets and Results of Measurement Sets Functions (cont'd)

MDCS_ASSAY_AssignAttributeValueFloat	To assign an Assay attribute value (type Float)
MDCS_ASSAY_RenameAttribute	To rename an Assay attribute
MDCS_ASSAY_DeleteAttribute	To delete an Assay attribute based on its DB name
MDCS_ASSAY_GetHeaderAndFileInfo	To get header info and file import location of a measurement set
MDCS_ASSAY_UpdateMeasurementSetName	To update a measurement set name
MDCS_ASSAY_UpdateMeasurementSetDescription	To update a measurement set description
MDCS_ASSAY_GetUniqueMeasurementValuesByPlate	To get unique values of measurement types for a plate

MDCS_ASSAY_AppendMeasurementSet

```
BOOL MDCS_ASSAY_AppendMeasurementSet(
HDBHANDLE hHandle,
const MDCS_CL_ImportDS * pDatasource,
MDCS_ProgressCallback* pCallback
);
```

Purpose

To import measurement set data into the database.

Parameters

hHandle - database handlepDatasource - pointer to a datasource datapCallBack - progress callback class

Return

FALSE - if error occurred

MDCS_ASSAY_ImportMeasurementSet

```
BOOL MDCS_ASSAY_RemoveShapeData(
HDBHANDLE hHandle,
LONGLONG lAssayID,
);
```

Purpose

To remove shape data from a measurement set.

Parameters

```
hHandle - database handle 
IAssayID - assay ID
```

Return

FALSE - if error occurred

MDCS ASSAY ImportMeasurementSet

```
BOOL MDCS_ASSAY_ImportMeasurementSet(
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_ImportMeasurementSet* pDatasource,
MDCS_ProgressCallback* pCallback
);
```

Purpose

To import measurement set data into the database.

Parameters

```
hHandle - database handle

lAssayID - assay ID

pDatasource - pointer to a datasource data

pCallBack - progress callback class
```

Return

FALSE - if error occurred

MDCS_ASSAY_InsertDValue

```
BOOL MDCS_ASSAY_InsertDValue(
HDBHANDLE hHandle,
const MDCS_ST_ShapeInfo* stShapeInfo,
LONGLONG lAssayID,
LONGLONG lMeasurementID,
double dValue
);
```

Purpose

To insert the value of a double type to the table that contains results of a measurement set.

Parameters

```
hHandle - database connection handle pstShapeInfo - pointer to a structure that describes shape lAssayID - ID of an assay that was run MeasurementID - measurement set ID of an assay dValue - value to insert
```

Return

FALSE - if error occurred

MDCS ASSAY InsertNumericValue

```
BOOL MDCS_ASSAY_InsertNumericValue(
HDBHANDLE hHandle,

LPCSTR pszTableName,

LPCSTR pszColumnName,

const MDCS_ST_ShapeInfo* stShapeInfo,
double dValue
);
```

Purpose

To insert the value of a double type to the table that contains results of a measurement set.

Parameters

hHandle - database connection handle pstShapeInfo - pointer to a structure that describes shape pszTableName - table name where values are inserted pszColumnName - column name of a measurement dValue - value to insert

Return

FALSE - if error occurred

MDCS_ASSAY_InsertSValue

```
BOOL MDCS_ASSAY_InsertsValue(
DBHANDLE hHandle,
const MDCS_ST_ShapeInfo* pstShapeInfo,
LONGLONG lAssayID,
LONGLONG lMeasurementID,
LPCSTR pszValue
);
```

Purpose

To insert the value of a string to the table that contains results of a measurement set.

Parameters

hHandle - database connection handle pstShapeInfo - pointer to a structure that describes shape IAssayID - ID of an assay that was run IMeasurementID - measurement set ID of an assay pszValue - value to insert

Return

FALSE - if error occurred

MDCS_ASSAY_InsertStringValue

```
BOOL MDCS_ASSAY_InsertStringValue(
HDBHANDLE hHandle,

LPCSTR pszTableName,

LPCSTR pszColumnName,

const MDCS_ST_ShapeInfo* stShapeInfo,

LPCSTR pszValue
);
```

Purpose

To insert the value of a string type to the table that contains results of measurement set.

Parameters

hHandle - database connection handle pstShapeInfo - pointer to a structure that describes shape pszTableName - table name where values are inserted pszColumnName - column name of a measurement pszValue - value to insert

Return

FALSE - if error occurred

MDCS ASSAY InsertShapeLines

```
BOOL MDCS_ASSAY_InsertShapeLines(
HDBHANDLE hHandle,

LPCSTR pszTableName,

const MDCS_ST_ShapeInfo* pstShapeInfo,

const MDCS_ST_ShapeLine* pShapeLines,

UINT nNumShapeLines
);
```

Purpose

To insert or update shapes.

Parameters

hHandle - database connection handle
pstShapeInfo - pointer to a structure that describes shape
pszTableName - table name where values are inserted
pShapeLines - pointer to the array of structures that describes shape
lines
nNumShapeLines - number of the shape lines in the array pShapeLines

Return

FALSE - if error occurred

MDCS_ASSAY_UpdateShapeLines

```
BOOL MDCS_ASSAY_UpdateShapeLines(
HDBHANDLE hHandle,

LPCSTR pszTableName,

LONGLONG lShapeID,

const MDCS_ST_ShapeInfo* pstShapeInfo,

const MDCS_ST_ShapeLine* pShapeLines,

UINT nNumShapeLines
);
```

Purpose

To update shape lines.

Parameters

hHandle - database connection handle
IShapeID - shape ID to update
pstShapeInfo - pointer to a structure that describes shape
pszTableName - table name where values are inserted
pShapeLines - pointer to the array of structures that describe shape
lines
nNumShapeLines - number of the shape lines in the array
(pShapelines)

Return

FALSE - if error occurred

MDCS_ASSAY_GetShapeLineBlobDesc

```
BOOL MDCS_ASSAY_GetShapeLineBlobDesc(
HDBHANDLE hHandle,

LPCSTR pszTableName,

const MDCS_ST_ShapeInfo* pstShapeInfo,

MDCS_ST_ShapeLinesBlob *pBlobDesc,

LONGLONG* plRecordsFound
);
```

Purpose

To get shape lines BLOB description.

Parameters

hHandle - database connection handle pstShapeInfo - pointer to a structure that describes shape pszTableName - table name where values are inserted pBlobDesc - pointer to the structure that describes a blob of shape lines

Output

plRecordsFound - number of records found

Return

FALSE - if error occurred

MDCS_ASSAY_GetShapeLineBlobDescBySite

```
BOOL MDCS_ASSAY_GetShapeLineBlobDescBySite(
HDBHANDLE hHandle,

LPCSTR pszTableName,

LONGLONG ISiteID

MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get shape lines BLOB description per site.

Parameters

```
hHandle - database connection handle
pszTableName - table name where values are inserted
|SiteID - site ID to get BLOBs for
pResultCallback - pointer to callback function to get data
```

Output

```
Columns:
SITE_ID
SERIES_ID
INSTANCE
SHAPE_TYPE
```

COLOR

 ${\it LINES_TOTAL} \ - \ {\it total lines} \ counted \ for \ the \ shapes \ per \ combination \ (SITE_ID, SERIES_ID, INSTANCE, SHAPE_TYPE, COLOR)$

 $\ensuremath{\mathsf{RECORD_COUNT}}$ - total records found for the shapes per combination above

BLOB_SIZE - data size to retrieve for the shapes per combination above

Return

FALSE - if error occurred

MDCS_ASSAY_GetShapeLineBlobDescBySiteAndSeries

```
BOOL MDCS_ASSAY_GetShapeLineBlobDescBySiteAndSeries(
HDBHANDLE hHandle,

LPCSTR pszTableName,

LONGLONG lSiteID,

LONGLONG lSeriesID,

MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get shape lines BLOB description by site and series.

Parameters

```
hHandle - database connection handle
pszTableName - table name where values are inserted
ISiteID - site ID where to get BLOBs from
ISeriesID - series ID
pResultCallback - pointer to a callback function to get data
```

Output

Columns:

```
SITE_ID

SERIES_ID

INSTANCE

SHAPE_TYPE

COLOR

LINES_TOTAL - total lines counted for the shapes per combination (SITE_ID, SERIES_ID, INSTANCE, SHAPE_TYPE, and COLOR)

RECORD_COUNT - total records found for the shapes per combination above

BLOB_SIZE - data size to retrieve for the shapes per combination above
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetShapeLineBlobDescBySiteSeriesAndInstance

```
BOOL
MDCS_ASSAY_GetShapeLineBlobDescBySiteSeriesAndInstance(
HDBHANDLE hHandle,
LPCSTR pszTableName,
LONGLONG lSiteID,
LONGLONG lSeriesID,
LONGLONG lInstanceID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get shape lines BLOB description per site, series and instance.

Parameters

```
hHandle - database connection handle
pszTableName - table name where values are inserted
ISiteID - site ID where to get BLOBs from
ISeriesID - series ID
IInstance - instance ID
pResultCallback - pointer to callback function to get data
```

Output

```
Columns:

SITE_ID - site ID

SERIES_ID

INSTANCE

SHAPE_TYPE

COLOR

LINES_TOTAL - total lines counted for the shapes per combination (SITE_ID, SERIES_ID, INSTANCE, SHAPE_TYPE, and COLOR)

RECORD_COUNT - total records found for the shapes per combination above

BLOB_SIZE - data size to retrieve for the shapes per combination above
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetShapeLineBlobDescByPlate

```
BOOL MDCS_ASSAY_GetShapeLineBlobDescByPlate(
HDBHANDLE hHandle,

LPCSTR pszTableName,

LONGLONG lPlate,

MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get shape lines BLOB description per plate.

Parameters

```
hHandle - database connection handlepszTableName - table name where values are insertedIPlateID - plate IDpResultCallback - pointer to callback function to get data
```

Output

```
Columns:

SITE_ID - site ID

SERIES_ID

INSTANCE

SHAPE_TYPE

COLOR

LINES_TOTAL - total lines counted for the shapes per combination

SITE_ID, SERIES_ID, INSTANCE, SHAPE_TYPE, and COLOR

RECORD_COUNT - total records found for the shapes per combination above

BLOB_SIZE - data size to retrieve for the shapes per combination above
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetShapeLineBlobPerAssay

```
BOOL MDCS_ASSAY_GetShapeLineBlobPerAssay(
HDBHANDLE hHandle,

LPCSTR pszShapeName,

MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get all the shape descriptions per assay.

Parameters

hHandle - database connection handlepszShapeName - table name where shapes are storedpResultCallback - pointer to callback function to get data

Output

```
Columns:

SITE_ID - site ID

SERIES_ID

INSTANCE

SHAPE_TYPE

COLOR

LINES_TOTAL - total lines counted for the shapes per combination

SITE_ID, SERIES_ID, INSTANCE, and SHAPE_TYPE

RECORD_COUNT - total records found for the shapes per combination above

BLOB_SIZE - data size to retrieve for the shapes per combination
```

Return

above

FALSE - if error occurred

MDCS_ASSAY_GetShapeLines

```
BOOL MDCS_ASSAY_GetShapeLines(
HDBHANDLE hHandle,

LPCSTR pszTableName,
const MDCS_ST_ShapeInfo* pstShapeInfo,
MDCS_ST_ShapeLine *pShapeLines,

LONGLONG lNumShapeLines
);
```

Purpose

To get shape lines of a shape.

Parameters

hHandle - database connection handle

pstShapeInfo - pointer to a structure that describes shape

pszTableName - table name where values are inserted

pShapeLines - pointer to an array of structure that describes shape
lines

INumShapeLines - number of the shape lines the array in pShapeLines

Return

FALSE - if error occurred

MDCS ASSAY GetMeasurementAttributes

```
BOOL MDCS_ASSAY_GetMeasurementAttributes(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lMeasurementID,
MDCS_ST_COLUMNPROP* pAssMeas,
BOOL bInsert = FALSE
);
```

Purpose

Function to get a measurement description based on measurement set ID and assay ID.

```
    hHandle - database connection handle
    lAssayID - assay ID
    lMeasurementID - measurement ID
    bInsert - if TRUE, create a new record if it does not exist
```

Output

pAssMeas - structure that contains results of the query

Return

FALSE - if error occurred

MDCS_ASSAY_GetAssayMeasurementRecord

```
BOOL MDCS_ASSAY_GetAssayMeasurementRecord(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lMeasurementID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get a measurements description based on measurement set ID and assay ID.

Parameters

```
    hHandle - database connection handle
    lAssayID - assay ID
    lMeasurementID - measurement ID
    pResultCallback - pointer to callback function to get data
```

Output

Columns from the ASSAYS and Table_COLUMNS table

Return

FALSE - if error occurred

MDCS_ASSAY_GetAssayMeasurementsByPlateAndMeasurement

```
BOOL
MDCS_ASSAY_GetAssayMeasurementsByPlateAndMeasurement(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG* plMeasurementID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get measurement set measurement by plate ID and measurement ID.

Parameters

```
hHandle - database connection handle

lPlateID - plate ID

plMeasurementID - measurement ID of global measurement

pResultCallback - pointer to callback function to get data
```

Output

```
Columns:

TABLE_ID - table name

COLUMN_NAME

ASSAY_ID

COLUMN_TYPE

COLUMN_NAME_EXT
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementsBySiteID

```
BOOL MDCS_ASSAY_GetMeasurementsBySiteID(
HDBHANDLE hHandle,
LONGLONG lSiteID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get measurements by site ID.

Parameters

```
hHandle - database connection handlelSiteID - site IDpResultCallback - pointer to callback function to get data
```

Output

```
Columns:

TABLE_ID - table name

COLUMN_NAME

ASSAY_ID

COLUMN_TYPE

COLUMN_NAME_EXT
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementInfoByAssayAndColumnName

```
BOOL MDCS_ASSAY_GetMeasurementInfoByAssayAndColumnName(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPCSTR pszAttrDBColumnName,
MDCS_ST_COLUMNPROP& stAssayMeas
);
```

Purpose

Function to get a measurement description based on measurement set ID and database column name.

hHandle - database connection handle
 lAssayID - assay ID
 pszAttrDBColumnName - structure describes the attribute column in the database
 bInsert - if TRUE, create a new record if it does not exist

Output

stAssayMeas - structure that contains results of the query

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementByFunctionAndParameterName

```
BOOL
MDCS_ASSAY_GetMeasurementByFunctionAndParameterName(
HDBHANDLE hHandle,

LPCSTR pszFunctionName,

LPCSTR pszParameterName,

const MDCS_E_ColumnType& eColType,

MDCS_ST_COLUMNPROP& stAssayMeas
);
```

Purpose

Function to get a measurement attribute description based on parameter and functions names.

Parameters

```
hHandle - database connection handle pszFunctionName - function name pszParameterName - parameter name eColType - data format
```

Output

stAssayMeas - structure that contains results of the query

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementRecord

```
BOOL MDCS_ASSAY_GetMeasurementRecord(
HDBHANDLE hHandle,
LONGLONG lColumnID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get a data type record by its ID.

Parameters

hHandle - database connection handlepResultCallback - pointer to a callback function to get data

Output

```
Returns records with columns:

ID - measurement ID

FUNCTION_NAME - function name

PARAMETER_NAME - parameter name

COLUMN_NAME_EXT - display name of column

COLUMN_NAME - data type name as in DATABASE

COLUMN_DESCRIPTION - data type description

COLUMN_TYPE - column type

TABLE_ID - name of the table where can be found

ENTITY_ID - measurement set ID
```

Return

FALSE - if error occurred

MDCS_ASSAY_Create

```
BOOL MDCS_ASSAY_Create(
HDBHANDLE hHandle,
const MDCS_ST_Assay* pAssayIn,
MDCS_ST_Assay* pAssayOut,
LONGLONG lFolderID
);
```

Purpose

Function to create a new measurement set.

hHandle - database connection handle

pAssayIn - pointer to a measurement set structure that will be used to create new measurement set

IFolderID - ID of the folder where measurement set will be created,
if 0 - will create measurement set in the user's default folder

Output

pAssayOut - created measurement set

Return

FALSE - if error occurred

MDCS_ASSAY_GetAssaySiblingFolders

```
BOOL MDCS_ASSAY_GetAssaySiblingFolders(
HDBHANDLE hHandle,
LONGLONG lParentFolderID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get measurement set sibling folders.

Parameters

hHandle - database connection handleIParentFolderID - ID of parent to retrieve siblings ofpResultCallback - pointer to a callback function to get data

Output result columns will be

FOLDER_NAME - folder name FOLDER ID - folder ID

Return

FALSE - if error occurred

MDCS_ASSAY_CreateMeasurement

```
BOOL MDCS_ASSAY_CreateMeasurement(
HDBHANDLE hHandle,
const MDCS_ST_COLUMNPROP* pstMeasurementIn,
MDCS_ST_COLUMNPROP* pstMeasurementOut
);
```

Purpose

Function to create a new measurement set measurement.

Parameters

hHandle - database connection handle
pstMeasurementIn - measurement description

Output

pstMeasurementOut - measurement description



Note: Either szColumnNameFull or the pair of szFunctionName and szMeasName are required, column type is required.

Return

FALSE - if error occurred

MDCS_ASSAY_AddMeasurement

```
BOOL MDCS_ASSAY_AddMeasurement(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lMeasID,
MDCS_ST_COLUMNPROP* pCreatedMeasurementOut
);
```

Purpose

Function to add a new measurement to a measurement set.

```
hHandle - database connection handlelAssayID - assay ID for which to add measurementlMeasID - measurement ID to add
```

Output

pCreatedMeasurementOut - pointer to structure that contains created measurement description

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementByID

```
BOOL MDCS_ASSAY_GetMeasurementByID(
HDBHANDLE hHandle,

LONGLONG lMeasuremntID,

MDCS_ST_COLUMNPROP* pstMeasurement
);
```

Purpose

Function to get a description of a measurement set measurement.

Parameters

```
hHandle - database connection handle IMeasuremntID - measurement ID
```

Output

pstMeasurement - measurement description

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementByName

```
BOOL MDCS_ASSAY_GetMeasurementByName(
HDBHANDLE hHandle,

LPCSTR pszName,

MDCS_ST_COLUMNPROP* pstMeasurement
);
```

Purpose

To get a description of a measurement set measurement.

Parameters

```
hHandle - database connection handlepszName - measurement name
```

Output

pstMeasurement - measurement description

Return

FALSE - if error occurred

MDCS_ASSAY_Delete

```
BOOL MDCS_ASSAY_Delete(
HDBHANDLE hHandle,
LONGLONG lAssayID,
BOOL bCanRestore = FALSE
);
```

Purpose

Function to delete an assay.

Parameters

hHandle - database connection handlelAssayID - assay that needs to be deletedbCanRestore - flag if TRUE, indicates that assay can be restored from Recycle bin

Return

FALSE - if error occurred

MDCS_ASSAY_DeleteAllForPlate

```
BOOL MDCS_ASSAY_DeleteAllForPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID
);
```

Purpose

Function to delete all assays for a plate.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Return

FALSE - if error occurred

MDCS ASSAY DeleteDataForPlate

```
BOOL MDCS_ASSAY_DeleteDataForPlate(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lPlateID
);
```

Purpose

Function to delete assay results for a plate.

Parameters

```
hHandle - database connection handleIPlateID - plate IDIAssayID - assay ID
```

Return

FALSE - if error occurred

MDCS_ASSAY_DeleteFolder

```
BOOL MDCS_ASSAY_DeleteFolder(
HDBHANDLE hHandle,
LONGLONG lFolderID
);
```

Purpose

Function to delete an assay folder.

Parameters

hHandle - database connection handlelFolderID - folder ID in measurement sets tree

Return

FALSE - if error occurred

MDCS ASSAY CreateOutlinesTable

```
BOOL MDCS_ASSAY_CreateOutlinesTable(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPSTR pszTableName,
int nSize,
int* pnSizeOut
);
```

Purpose

Function to create an assay outlines table.

Parameters

```
hHandle - database connection handlelAssayID - assay IDnSize - size of table name
```

Output

pszTableName - name of the table that contains outlines pnSizeOut - actual size the table name

Return

FALSE - if error occurred

MDCS_ASSAY_ManageFolderSecurity

```
BOOL MDCS_ASSAY_ManageFolderSecurity(
HDBHANDLE hHandle,
LONGLONG lFolderID,
HWND hWnd = NULL,
LPCSTR pszDlgTitle = NULL
);
```

Purpose

Calls a dialog to manage security access to the assay folders.

Parameters

```
hHandle - database connection handleIFolderID - folder IDpszDlgTitle - dialog titlehWnd - application window
```

Return

FALSE - if error occurred

MDCS_ASSAY_CanModifyFolder

```
BOOL MDCS_ASSAY_CanModifyFolder(
HDBHANDLE hHandle,
LONGLONG lFolderID,
BOOL& bCanModify
);
```

Purpose

To check if a user can modify a folder.

Parameters

```
hHandle - database connection handlelFolderID - folder ID
```

Output

bCanModify - if TRUE - user can modify the folder

Return

FALSE - if error occurred

MDCS_ASSAY_CreateFolder

```
BOOL MDCS_ASSAY_CreateFolder(
HDBHANDLE hHandle,
const MDCS_ST_FolderInfo& stInfoIn,
MDCS_ST_FolderInfo & stInfoOut
);
```

Purpose

To create a new assay folder in the database.

Parameters

hHandle - database connection handlestInfoIn - structure describes input folder information.stInfoOut - structure describes output (created) folder information

Return

FALSE - if error occurred

MDCS_ASSAY_ModifyFolder

```
BOOL MDCS_ASSAY_ModifyFolder(
HDBHANDLE hHandle,
const MDCS_ST_FolderInfo& stInfoIn
);
```

Purpose

To modify an assay folder in the database.

Parameters

hHandle - database connection handle *stInfoIn* - structure describes folder information.

Return

FALSE - if error occurred

MDCS_ASSAY_DoesSubFolderExist

```
BOOL MDCS_ASSAY_DoesSubFolderExist(
HDBHANDLE hHandle,
LONGLONG lFolderID,
LPCSTR pzName,
LONGLONG &lSubfolderID
);
```

Purpose

Check to see if an assay folder exists in the database.

Parameters

```
hHandle - database connection handlepzName - folder ID
```

Output

 ${\it ISubfolderID}$ - subfolder ID, if 0 subfolder does not exist, otherwise - folder ID

Return

FALSE - if folder does not exist

MDCS_ASSAY_GetAllInFolder

```
BOOL MDCS_ASSAY_GetAllInFolder(
HDBHANDLE hHandle,
LONGLONG lFolderID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get assays with basic descriptions in a folder.

Parameters

```
hHandle - database connection handlelFolderID - measurement set folder IDpResultCallback - pointer to a callback function to get data
```

Output

Fields that are always present

ASSAY_ID, ASSAY_NAME, SETTINGS_NAME, TIME_CREATED (as timeout value), TABLE_ID (table name where results are stored),

SHAPE_TABLE_NAME (name of the table where shape results are
stored)

MDCS_ASSAY_GetFolderPath

```
BOOL MDCS_ASSAY_GetFolderPath(
HDBHANDLE hHandle,
LONGLONG lFolderID,
AxString& strPath
);
```

Purpose

To get the path of a folder.

Parameters

hHandle - database connection handle IFolderID - measurement set folder ID

Output

strPath - will contain the path

MDCS_ASSAY_Reindex

```
BOOL MDCS_ASSAY_Reindex(
HDBHANDLE hHandle,
LONGLONG lAssayID
);
```

Purpose

To reindex assay measurements.

Parameters

hHandle - database connection handlelAssayID - assay ID

MDCS_ASSAY_ReindexShapes

```
BOOL MDCS_ASSAY_ReindexShapes(
HDBHANDLE hHandle,
LONGLONG lAssayID
);
```

Purpose

To reindex assay cell shapes.

Parameters

Handle - database connection handle IAssayID - assay folder ID

MDCS_ASSAY_OptimizeAll

```
BOOL MDCS_ASSAY_OptimizeAll(
HDBHANDLE hHandle,
MDCS_ProgressCallback * pCallBack
);
```

Purpose

To optimize all assay tables.

Parameters

hHandle - database connection handlepCallBack - progress callback

Return

FALSE - if error occurred

MDCS_ASSAY_Copy

```
BOOL MDCS_ASSAY_Copy(
HDBHANDLE hHandle,
LONGLONG lAssayIDSource,
LONGLONG lAssayIDDestination,
MDCS_ProgressCallback * pCallBack
);
```

Purpose

Function to copy an assay.

Parameters

hHandle - database connection handle lAssayIDSource - source assay ID lAssayIDDestination - destination assay ID pCallBack - progress callback

Return

FALSE - if error occurred

MDCS_ASSAY_Merge

```
BOOL MDCS_ASSAY_Merge(
HDBHANDLE hHandle,
LONGLONG lAssayIDSource,
LONGLONG lAssayIDDestination,
BOOL bCopyObjects,
MDCS_ProgressCallback * pCallBack
);
```

Purpose

Function to merge assays.

Parameters

hHandle - database connection handle
lAssayIDSource - source assay ID
lAssayIDDestination - destination assay ID
bCopyObjects - option indicates that attachments should be copied
pCallBack - progress callback

Return

FALSE - if error occurred

MDCS ASSAY GetSiteCount

```
BOOL MDCS_ASSAY_GetSiteCount(
HDBHANDLE hHandle,

LPCSTR pszTableName,

LONGLONG lSiteID,

LONGLONG& lSiteCount
);
```

Purpose

Function to count site appearance in the assay table.

Parameters

```
hHandle - database connection handleISiteID - site ID to look forpszTableName - name of a table that contains outlines
```

Output

ISiteCount - number of times site appears in the assay table

MDCS_ASSAY_GetSiteInfoImageByID

```
BOOL MDCS_ASSAY_GetSiteInfoImageByID(
HDBHANDLE hHandle,
LONGLONG lImageID,
MDCS_ST_SiteImageInfo& stInfoOut
);
```

Purpose

Function to get information about the image using the image ID.

Parameters

```
hHandle - database connection handlelImageID - image ID
```

Output

StInfoOut - structure of site information

Return

FALSE - if error occurred

MDCS_ASSAY_GetSiteInfoImageBySiteAndSeriesID

```
BOOL MDCS_ASSAY_GetSiteInfoImageBySiteAndSeriesID(
HDBHANDLE hHandle,
LONGLONG lSiteID,
LONGLONG lSeriesID,
MDCS_ST_SiteImageInfo& stInfoOut
);
```

Purpose

Function to get information about the image using the Site and Series ID.

Parameters

```
hHandle - database connection handleISiteID - Site IDISeriesID - series ID
```

Output

StInfoOut - structure of Site information

Return

FALSE - if error occurred

MDCS_ASSAY_MeasurementGetSiteCount

```
BOOL MDCS_ASSAY_MeasurementGetSiteCount(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lSiteID,
LONGLONG& lSiteCount
);
```

Purpose

Function to count site appearance in an assay.

Parameters

```
hHandle - database connection handlelAssayID - assay IDlSiteID - site ID
```

Output

ISiteCount - number of times site appears in the assay table

Return

FALSE - if error occurred

MDCS_ASSAY_DeleteMeasurement

```
BOOL MDCS_ASSAY_DeleteMeasurement(
HDBHANDLE hHandle,
LONGLONG lMeasurementID
);
```

Purpose

Function to delete a measurement.

Parameters

hHandle - database connection handle IMeasurementID - measurement ID

Return

FALSE - if error occurred

MDCS_ASSAY_CellOutlinesGetSiteCount

```
BOOL MDCS_ASSAY_CellOutlinesGetSiteCount(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lSiteID,
LONGLONG& lSiteCount
);
```

Purpose

Function to count site appearance in assay cell outlines.

Parameters

```
hHandle - database connection handlelAssayID - assay IDlSiteID - site ID
```

Output

ISiteCount - number of times site appears in the assay table

Return

FALSE - if error occurred

MDCS_ASSAY_CreateNewName

Purpose

```
BOOL MDCS_ASSAY_CreateNewName(
HDBHANDLE hHandle,
LONGLONG lFolderID,
LPCSTR pszOrigName,
AxString& strNewName,
BOOL bCreateAsCopy
);
```

Function to create a new assay name.

hHandle - database connection handle

IFolderID - destination folder ID

pszOrigName - to check availability or create a new assay name in destination folder

strNewName - to check availability or create a new assay name in destination folder

bCreateAsCopy -if TRUE - the name is created using the prefix "Copy of <orig name>", if the name exists

Return

FALSE - if error occurred

MDCS_ASSAY_CreateRun

```
BOOL MDCS_ASSAY_CreateRun(
HDBHANDLE hHandle,
const MDCS_ST_AssayRun& stAssayRun,
LONGLONG* plAssayRunID
);
```

Purpose

To create a new assay run record.

Parameters

hHandle - database connection handlestAssayRun - assay run description

Output

plAssayRunID - generated assay runID

Return

FALSE - if error occurred

MDCS_ASSAY_GetLatestAssayRunID

```
BOOL MDCS_ASSAY_GetLatestAssayRunID(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG* plRunID
);
```

Purpose

To get ID of the latest run on an assay.

```
hHandle - database connection handlelAssayID - assay ID
```

Output

```
plRunID - assay run ID
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetByID

```
BOOL MDCS_ASSAY_GetByID(
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_ST_Assay* pstAssay
);
```

Purpose

Function to get assay data by ID.

Parameters

```
hHandle - database connection handle lAssayID - assay ID
```

Output

pstAssay - structure with assay information

Return

FALSE - if error occurred

MDCS_ASSAY_CreateProfile

```
BOOL MDCS_ASSAY_CreateProfile(
HDBHANDLE hHandle,
const MDCS_ST_AssayProfile& stAssayProfile,
LONGLONG* lAssayProfileID
);
```

Purpose

To create a new assay profile record.

```
hHandle - database connection handlestAssayProfile - assay profile description
```

Output

IAssayProfileID - generated assay profile ID

Return

FALSE - if error occurred

MDCS_ASSAY_AssociateWithPlate

```
BOOL MDCS_ASSAY_AssociateWithPlate(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lPlateID,
BOOL* bRecordExisted,
LONGLONG* plPlateAssayId = NULL
);
```

Purpose

To associate an assay with a plate.

Parameters

```
hHandle - database connection handlelAssayID - assay IDlPlateID - plate ID
```

Output

```
bRecordExisted - if TRUE - indicates that record already existed in
database
plPlateAssayId - assay plate ID
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetSpotID

```
BOOL MDCS_ASSAY_GetSpotID(
HDBHANDLE hHandle,
LONGLONG lSpot,
LONGLONG lInstance,
LONGLONG lSeries,
LONGLONG lSite,
LONGLONG& lSpotID
);
```

Purpose

To get the spot ID from the instance, site and series.

Parameters

```
hHandle - database connection handle ISpot - Measurements table number IInstance - instance ID ISite - site ID ISeries - series ID
```

Output

ISpotID - Spot ID

Return

FALSE - if error occurred

MDCS_ASSAY_GetAllForPlate

```
BOOL MDCS_ASSAY_GetAllForPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

Function to get all assays that are associated with a plate.

```
hHandle - database connection handleIPlateID - ID of a platepResultCallback - pointer to a callback function to get data
```

Output

```
Columns:
```

ASSAY_ID - assay ID
ASSAY_NAME - assay name
ASSAY_DESC - assay description
CREATOR_NAME - name of the user who created the assay

Return

FALSE - if error occurred

MDCS_ASSAY_GetProfiles

```
BOOL MDCS_ASSAY_GetProfiles(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

Function to get all available profiles.

Parameters

hHandle - database connection handlepResultCallback - pointer to a callback function to get data

Output

Columns:

ID - unique ID

NAME - record name

DESCRIPTION - record description

SETTINGS_NAME - name of the settings

OPERATOR - name of the user who created settings

DATA_STORAGE_TYPE - type of the data storage

SETTINGS_STORAGE_TYPE - type of a storage where the settings BLOB is stored (DATABASE or File Server)

Return

FALSE - if error occurred

MDCS_ASSAY_GetProfile

```
BOOL MDCS_ASSAY_GetProfile(
HDBHANDLE hHandle,
LONGLONG lProfileID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

Function to get a profile.

Parameters

```
hHandle - database connection handleIProfileID - profile IDpResultCallback - pointer to a callback function to get data
```

Output

```
Columns:

ID - unique ID

NAME - record name

DESCRIPTION - record description

SETTINGS_NAME - name of the settings

OPERATOR - name of user who created settings

DATA_STORAGE_TYPE - type of the data storage

SETTINGS_STORAGE_TYPE - type of a storage where settings BLOB is stored (DATABASE or File Server)
```

Return

FALSE - if error occurred

MDCS ASSAY GetProfileInfo

```
BOOL MDCS_ASSAY_GetProfileInfo(
HDBHANDLE hHandle,
LONGLONG lProfileID,
MDCS_ST_AssayProfile& stProfile
);
```

Purpose

Function to get profile information.

```
hHandle - database connection handleIProfileID - profile ID
```

Output

```
Returns columns: 
stProfile - structure that contains profile information
```

Return

FALSE - if error occurred

MDCS_ASSAY_UpdateProfile

```
BOOL MDCS_ASSAY_UpdateProfile(
HDBHANDLE hHandle,
const MDCS_ST_AssayProfile& stProfileToUpdate
);
```

Purpose

Function to update a profile record.

Parameters

hHandle - database connection handle stProfileToUpdate - profile to update

Return

FALSE - if error occurred

MDCS_ASSAY_DeleteProfile

```
BOOL MDCS_ASSAY_DeleteProfile(
HDBHANDLE hHandle,
LONGLONG lProfileID
);
```

Purpose

To delete an assay profile.

Parameters

```
hHandle - database connection handle
lProfileID - profile ID
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetRecord

```
BOOL MDCS_ASSAY_GetRecord(
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

Function to get an assay record.

Parameters

```
hHandle - database connection handlelAssayID - assay IDpResultCallback - callback class
```

Output

```
Columns:

ID - unique ID

NAME - record name

DESCRIPTION - record description

SETTINGS_NAME - name of the settings

DATA_STORAGE_TYPE - type of the data storage

SETTINGS_STORAGE_TYPE - type of a storage where the settings BLOB is stored (DATABASE or File Server)
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetAssayByRunID

```
BOOL MDCS_ASSAY_GetAssayByRunID(
HDBHANDLE hHandle,
LONGLONG lRunID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

Function to get a measurement set record using the RUN ID.

```
hHandle - database connection handlelRunID - run IDpResultCallback - pointer to callback function to get data
```

Output

Columns:

ID - unique ID

NAME - record name

DESCRIPTION - record description

SETTINGS_NAME - name of the settings

DATA_STORAGE_TYPE - type of the data storage

SETTINGS_STORAGE_TYPE - type of a storage where the settings BLOB is stored (DATABASE or File Server)

Return

FALSE - if error occurred

MDCS_ASSAY_UsedInDatasets

```
BOOL MDCS_ASSAY_UsedInDatasets(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG& lDatasets
);
```

Purpose

Function to check if the assay is used in datasets.

Parameters

```
hHandle - database connection handleAssayID - assay ID
```

Output

IDatasets - number of datasets where the assay is used

Return

FALSE - if dialog cancelled.

MDCS_ASSAY_HavePermissionsToModify

```
BOOL MDCS_ASSAY_HavePermissionsToModify(
HDBHANDLE hHandle,
LONGLONG lAssayID,
BOOL& bCanModify
);
```

Purpose

Function to check if the current user can modify an assay.

Parameters

```
hHandle - database connection handlelAssayID - assay ID
```

Output

bCanModify - if TRUE - user can modify the assay

Return

FALSE - if dialog cancelled

MDCS_ASSAY_GetImageSourceRecords

```
BOOL MDCS_ASSAY_GetImageSourceRecords(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get all records from the IMAGE_SOURCE table.

Parameters

hHandle - database connection handle pResultCallback - pointer to callback function to get data

Output

```
Returns records with columns:

ID - image source ID

ACQ_INSTANCE_ID - acquisition instance ID

ASSAY_INSTANCE_ID - measurement set ID

SOURCE_DESCRIPTION - source description

SOURCE_ILLUMINATION - source illumination
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetStatisticalValuesForPlateAndAssay

```
BOOL MDCS_ASSAY_GetStatisticalValuesForPlateAndAssay(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG lAssayID,
LONGLONG lMeasurementID,
const MDCS_E_Statistic& eStatistic,
BOOL bGroupSites,
BOOL bAllSites,
int nSiteX,
int nSiteY,
LPCSTR pszSeriesInfoIndex,
LPCSTR pszFilter,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get statistical values for the selected plate and assay.

Parameters

```
hHandle - database connection handle

pResultCallback - pointer to a callback function to get data

lPlateID - plate ID

lAssayID - assay ID

lMeasurementID - measurement ID

eStatistic - type of statistic to apply to data

bGroupSites - if TRUE, the data will be grouped by site

bAllSites - if TRUE, the data will be fetched for all series for a site

nSiteX - X position of a site

pszSeriesInfoIndex - series information index

pszFilter - filter that should be applied to the data
```

Output

Contains columns with statistical values

Return

FALSE - if error occurred

MDCS_ASSAY_GetStatisticalValuesForPlateAndMeasurement

```
BOOL
MDCS_ASSAY_GetStatisticalValuesForPlateAndMeasurement(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG lMeasurementID,
const MDCS_E_Statistic& eStatistic,
BOOL bGroupSites,
BOOL bAllSites,
int nSiteX,
int nSiteY,
LPCSTR pszSeriesInfoIndex,
LPCSTR pszFilter,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get statistical values for the selected plate and measurement.

Parameters

```
hHandle - database connection handle

pResultCallback - pointer to a callback function to get data

IPlateID - plate ID

IMeasurementID - measurement ID

eStatistic - type of statistic to apply to data

bGroupSites - if TRUE, the data is grouped by site

bAllSites - if TRUE, the data is fetched for all series for a site

nSiteX - X position of a site

pszSeriesInfoIndex - series information index

pszFilter - filter that should be applied to the data
```

Output

Contains columns with statistical values

Return

FALSE - if error occurred

MDCS_ASSAY_GetValuesForPlateAndAssay

```
BOOL MDCS_ASSAY_GetValuesForPlateAndAssay(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG lAssayID,
LONGLONG lMeasurementID,
BOOL bSingleSite,
int nSiteX,
int nSiteY,
LPCSTR pszFilter,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get values for the selected plate and assay.

Parameters

```
hHandle - database connection handle
pResultCallback - pointer to a callback function to get data
IPlateID - plate ID
IAssayID - assay ID
IMeasurementID - measurement ID
bSingleSites - if TRUE, data will be fetched for single sites.
nSiteX - X position of a site
nSiteY - Y position of a site
pszFilter - filter that should be applied to the data
```

Output

Contains columns with statistical values

Return

FALSE - if error occurred

MDCS_ASSAY_GetValuesForPlate

```
BOOL MDCS_ASSAY_GetValuesForPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG lMeasurementID,
BOOL bSingleSite,
int nSiteX,
int nSiteY,
LPCSTR pszFilter,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get values for the selected plate and measurement.

Parameters

```
hHandle - database connection handle
pResultCallback - pointer to a callback function to get data
IPlateID - plate ID
IMeasurementID - measurement ID
bSingleSites - if TRUE, data will be fetched for single sites
nSiteX - X position of a site
nSiteY - Y position of a site
pszFilter - filter that should be applied to data
```

Output

Includes all values for a plate (including WELL_X, WELL_Y, INSTANCE, SERIES_INFO_ID, SITE_ID)

Return

FALSE - if error occurred

MDCS_ASSAY_GetAllSiteMeasurements

```
BOOL MDCS_ASSAY_GetAllSiteMeasurements(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lSiteID,
BOOL bNumeric,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get all measurements of a selected site.

Parameters

```
IAssayID - Assay ID

hHandle - database connection handle

pResultCallback - pointer to callback function to get data

ISiteID - site ID

bNumeric - if TRUE - fetch all numeric columns
```

Output

Contains site measurement columns from the TABLE_COLUMNS table

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementInfoByAssayAndPlate

```
BOOL MDCS_ASSAY_GetMeasurementInfoByAssayAndPlate(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get data types records by plate and assay.

Parameters

```
hHandle - database connection handlelAssayID - assay IDlPlateID - Plate IDpResultCallback - pointer to a callback function to get data
```

Output

```
Returns columns:

TABLE_ID - (table name)

COLUMN_NAME

ASSAY_ID

COLUMN_TYPE

COLUMN_NAME_EXTN (column type - S for string, N for numeric)

FUNCTION_NAME,

PARAMETER_NAME
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetAssayInfoByNamePlateAndSettings

```
BOOL MDCS_ASSAY_GetAssayInfoByNamePlateAndSettings(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LPCSTR pszAssayName,
LPCSTR pszSettingsName,
MDCS_ST_Assay& stAssayOut
);
```

Purpose

To get assay information by Plate ID, Assay name and settings.

Parameters

```
hHandle - database connection handle
IPlateID - plate ID
pszAssayName - assay name
pszSettingsName - settings name
```

Output

stAssayOut - structure that contains assay information

Return

FALSE - if error occurred

MDCS_ASSAY_GetAssayIDsOfPlate

```
BOOL MDCS_ASSAY_GetAssayIDsOfPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get assay IDs that are available for a plate.

Parameters

```
hHandle - database connection handleIPlateID - plate IDpResultCallback - pointer to callback function to get data
```

Output

Assay IDs for the plate

Return

FALSE - if error occurred

MDCS_ASSAY_GetAssaySettingsOfPlate

```
BOOL MDCS_ASSAY_GetAssaySettingsOfPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get unique measurement set settings of the plate.

Parameters

```
hHandle - database connection handlelPlateID - plate IDpResultCallback - pointer to callback function to get data
```

Output

```
Results set with fields:
ASSAY_ID
ASSAY_NAME
SETTINGS_NAME
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementStatistic

```
BOOL MDCS_ASSAY_GetMeasurementStatistic(
HDBHANDLE hHandle,
LONGLONG lMeasurementID,
MDCS_ST_MeasurementUseStatistic& stMeasurementSta
);
```

Purpose

To get the number of datasets and the number of measurement sets that use a measurement.

Parameters

hHandle - database connection handleMeasurementID - measurement ID

Output

stMeasurementSta – structure that contains statistic information about datasets and measurement sets that use the measurement

Return

FALSE - if error occurred

MDCS_ASSAY_GetShapeLinesBySite

```
BOOL MDCS_ASSAY_GetShapeLinesBySite (
HDBHANDLE hHandle,

LPCSTR pszTableName,

LONGLONG lSiteID,

MDCS_GetBlobColumnResults* pResultCallback,

LONGLONG lMaxSize = -1
);
```

Purpose

To get shape lines BLOB description per site.

Parameters

hHandle - database connection handle
 pszTableName - table name
 lSiteID - site id
 lMaxSize - max size
 pResultCallback - pointer to callback function to get data

Output

Result set with fields:

SITE ID

SERIES ID

INSTANCE

SHAPE_TYPE

LINES_TOTAL - total lines count for the shapes per combination

 $\ensuremath{\mathsf{RECORD_COUNT}}$ - total records found for the shapes per combination above

BLOB_SIZE - data size to retrieve for the shapes per combination above

Return

FALSE - if error occurred

MDCS_ASSAY_GetShapeLinesBySiteAndSeries

```
BOOL MDCS_ASSAY_GetShapeLinesBySiteAndSeries(
HDBHANDLE hHandle,

LPCSTR pszTableName,

LONGLONG lSiteID,

LONGLONG lSeriesID,

MDCS_GetBlobColumnResults* pResultCallback,

LONGLONG lMaxSize = -1);
```

Purpose

To get shape lines BLOB description per site by using site and series id

Parameters

hHandle - database connection handle

pszTableName - table name

ISiteID - site id

ISeriesID - series id

IMaxSize - max size

pResultCallback - pointer to callback function to get data

Output

Columns:

SITE ID - site id

SERIES ID

INSTANCE

SHAPE_TYPE

LINES_TOTAL - total lines count for the shapes per combination RECORD_COUNT - total records found for the shapes per combination above

BLOB_SIZE - data size to retrieve for the shapes per combination

Return

FALSE - if error occurred

MDCS_ASSAY_GetMarkedAssaysWithCallback

```
BOOL MDCS_ASSAY_GetMarkedAssaysWithCallback (
HDBHANDLE hHandle,

MDCS_GetDBResultsCCallback* pResultCallback,

BOOL bMarkPerm = TRUE,

BOOL bWithSecurity = FALSE);
```

Purpose

To get deleted assays using callback

Parameters

hHandle - database connection handle

bMarkPerm – flag to indicate to get permanently deleted assay or restorable assay

bWithSecurity – flag to indicate whether to include plate security when getting assay info

pResultCallback - pointer to callback function to get data

Output

Output structures that contains results of the query

Return

FALSE - if error occurred

MDCS_ASSAY_Restore

```
BOOL MDCS_ASSAY_GetMarkedAssaysWithCallback (
HDBHANDLE hHandle,
LONGLONG lAssayID);
```

Purpose

To restore a deleted assay from the recycle bin.

Parameters

```
hHandle - database connection handlelAssayID - assay ID
```

Output

None

Return

FALSE - if error occurred

MDCS_ASSAY_CreateForPlate

Purpose

To create an assay for a plate.

Parameters

```
    hHandle - database connection handle
    pAssayIn - a structure that describes assay informations (excluding the assay id)
    IPlateID - plate id associated with the assay
```

Output

pAssayOut – a structure that describe assay information about the newly created assay (including the assay id)

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementByName

```
BOOL MDCS_ASSAY_GetMeasurementByName(
HDBHANDLE hHandle,

LPCSTR pszName,

MDCS_ST_COLUMNPROP* pstMeasurement,

LPCSTR pszScopeName = MDCS_cszMDCSCellMeasurementScope);
```

Purpose

To get a description of an assay measurement.

Parameters

```
hHandle - database connection handlepszName - a name of the measurementpszScopeName - Scope name of the measurement
```

Output

pstMeasurement – output structure that describes the measurement

Return

FALSE - if error occurred

MDCS_ASSAY_GetMeasurementByDBName

```
BOOL MDCS_ASSAY_GetMeasurementByDBName (
HDBHANDLE hHandle,

LPCSTR pszColumnDBName,

MDCS_ST_COLUMNPROP& pstMeasurement,

LPCSTR pszScopeName = MDCS cszMDCSCellMeasurementScope);
```

Purpose

To get a description of an assay measurement using the database name of the measurement.

Parameters

```
hHandle - database connection handlepszColumnDBName - database name of the measurementpszScopeName - Scope name of the measurement
```

Output

pstMeasurement – output structure that describes the measurement

Return

FALSE - if error occurred

MDCS_ASSAY_GetScopeAttributeByName

```
BOOL MDCS_ASSAY_GetScopeAttributeByName(
HDBHANDLE hHandle,

LPCSTR pszName,

MDCS_ST_ScopeAttribute& stAttr,

LPCSTR pszScopeName = MDCS cszMDCSCellMeasurementScope);
```

Purpose

To get a description of a scope attribute using the display name.

Parameters

```
hHandle - database connection handlepszName - display name of the attributepszScopeName - scope name of the attribute
```

Output

stAttr - output structure that describe the attribute

Return

FALSE - if error occurred

MDCS_ASSAY_GetScopeAttributeByID

```
BOOL MDCS_ASSAY_GetScopeAttributeByID (
HDBHANDLE hHandle,

LPCSTR pszColumnID,

MDCS_ST_ScopeAttribute& stAttr,

LPCSTR pszScopeName = MDCS_cszMDCSCellMeasurementScope);
```

Purpose

To get a description of a scope attribute using the database name.

Parameters

```
hHandle - database connection handlepszColumnID - a database name of the attributepszScopeName - scope name of the attribute
```

Output

stAttr – output structure that describes the attribute

Return

FALSE - if error occurred.

MDCS_ASSAY_DeleteMeasurementInAssays

```
BOOL MDCS_ASSAY_DeleteMeasurementInAssays (
HDBHANDLE hHandle,

LPCSTR pszDBColumnName,

LONGLONG* arrAssayIDs,

INT_PTR nSize);
```

Purpose

To delete a measurement from all assays that use it but not the default. The default measurement is the measurement assay id = 0.

Parameters

```
hHandle - database connection handlepszDBColumnName - the database name of the attributearrAssayIDs - array of assay IDsnSize - size of arrAssayIDs
```

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_UpdateDataType

```
BOOL MDCS_ASSAY_UpdateDataType (
HDBHANDLE hHandle,
LONGLONG lMeasurementID,
Const MDCS_ST_COLUMNPROP& stDataTypeNew);
```

Purpose

To update a measurement property.

Parameters

```
hHandle - database connection handleIMeasurementID -- a measurement IDstDataTypeNew - a structure of the database column used to replace the original
```

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_UpdateDataTypeByAssay

```
BOOL MDCS_ASSAY_UpdateDataTypeByAssay
HDBHANDLE hHandle,
LPCSTR stOrigDBCol,
MDCS_ST_COLUMNPROP& stDataTypeNew,
LONGLONG lAssayID);
```

Purpose

To modify a measurement for an assay based on the database column name and assay ID, but not the default measurement (the measurement where assay ID = 0).

Parameters

```
hHandle - database connection handle
stOrigDBCol - database column name
stDataTypeNew - structure of a measurement used to replace the
original
lAssayID - assay ID
```

Output

NONE

Return

FALSE - if error occurred.

MDCS_ASSAY_CanModifyAssay

```
BOOL MDCS_ASSAY_CanModifyAssay
HDBHANDLE hHandle,
LONGLONG lAssayID,
BOOL& bCanModify);
```

Purpose

To check if the current user can modify an assay.

Parameters

```
hHandle - database connection handle
lAssayID - assay ID
```

Output

bCanModify - if TRUE user can modify assay

Return

FALSE - if error occurred

MDCS_ASSAY_CalculateStatisticResults

```
BOOL MDCS_ASSAY_CalculateStatisticResults(
HDBHANDLE hHandle,
const LONGLONG* plarrAssayIDs,
int nNumElementArray,
const AxStringArray& axarrFilters,
MDCS_E_Statistic* arrApplyStatistic,
int nStatCount,
const MDCS_ST_ScopeAttribute* arrStatisticColumns,
int nStatColumnCount,
const MDCS_ST_ScopeAttribute* axarrGroupBy,
int nGroupByCount,
const MDCS_ST_ScopeAttribute* axarrGroupBy,
int nGroupByCount,
MDCS_GetDBResultsCCallback *pResultCallback,
BOOL bGroupByAssay = FALSE);
```

Purpose

To get statistic results on array of assay ids.

Parameters

hHandle - database connection handle
plarrAssayIDs - array of assay IDs
arrApplyStatistic- statistic that will be applied to columns in
arrStatisticColumns
arrStatisticColumns - array of columns that will be used to calculate
statistic
axarrGroupBy - array of group by columns
axarrOrderBy - array of order by columns
axarrFilters - array of filters
bGroupByAssay - if TRUE - identifier of assay will be added
pResultCallback - pointer to a callback function to get data
nNumElementArray - number of elements in plarrAssayIDs
nStatCount - number of elements in arrApplyStatistic
nStatColumnCount - number of elements in arrStatisticColumns
nGroupByCount - number of elements in axarrGroupBy
nOrderByCount - number of elements in axarrGroupBy

Output

Recordset with fields: Statistic data for arrStatisticColumns columns as *STATS_<column name>* and columns specified in axarrGroupBy and axarrOrderBy. Columns:

in axarrOrderBy should match columns specified in axarrGroupBy

Return

FALSE - if error occurred.

MDCS_ASSAY_CalculateStatisticEx

```
BOOL MDCS ASSAY CalculateStatisticEx(
HDBHANDLE hHandle,
const LONGLONG* plarrAssayIDs,
int nNumElementArray,
const AxStringArray& axarrFilters,
const MDCS E Statistic* arrApplyStatistic,
int nStatCount,
const MDCS ST ScopeAttribute* arrStatisticColumns,
int nStatColumns,
const MDCS_ST_ScopeAttribute* axarrGroupBy,
int nGroupByColumns,
const MDCS ST ScopeAttribute* axarrOrderBy,
int nOrderByColumns,
const MDCS ST ScopeAttribute* axarrInnerGroupBy,
int nInnerGroupBy,
const MDCS_E_Statistic* arrInnerStatistic,
int nInnerStatCount,
MDCS GetDBResultsCCallback *pResultCallback,
BOOL bGroupByAssay = FALSE);
```

Purpose

To get statistic results on array of assay ids.

Parameters

hHandle - database connection handle plarrAssayIDs - array of assay IDs aarApplyStatistic - statistic that will be applied to columns in arrStatisticColumns arrStatisticColumns - array of columns that will be used to calculate statistic axarrGroupBy - array of group by columns axarrOrderBy - array of order by columns axarrFilters - array of filters bGroupByAssay - if TRUE - identifier of assay will be added arrInnerStatistic - statistic operation that will applied to internal query nNumElementArray - number of elements in plarrAssayIDs nStatCount - number of elements in arrApplyStatistic nStatColumns - number of elements in arrStatisticColumns nGroupByColumns - number of elements in axarrGroupBy nOrderByColumns - number of elements in axarrOrderBy *nInnerGroupBy* - number of elements in axarrInnerGroupBy pResultCallback - pointer to a callback function to get data nInnerStatCount - number of elements in arrInnerStatistic axarrInnerGroupBy - array of group by columns for internal query

Output

Recordset with fields:

Statistic data for arr Statistic
Columns columns as $\it STATS_< column name > columns$ specified in axarr Group
By

Columns:

in axarrOrderBy should match columns specified in axarrGroupBy

Return

FALSE - if error occurred.

MDCS_ASSAY_GetZPrime

```
BOOL MDCS_ASSAY_GetZPrime (
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPCSTR pszFilter,
LPCSTR pszMeasurementType,
MDCS_GetDBResultsCCallback *pResultCallback
);
```

Purpose

To calculate Z' on assay and measurement type.

Parameters

```
hHandle - database connection handle
lAssayID - assay ID
pszFilter - filter
pszMeasurementType - measurement type to use in calculation
pResultCallback - pointer to a callback function to get data
```

Output

```
Recordset with fields:
Z_PRIME - calculated Z'
PLATE_ID - plate ID
```

Return

FALSE - if error occurred

MDCS_ASSAY_GetZPrimeScopeAttribute

```
BOOL MDCS_ASSAY_GetZPrimeScopeAttribute (
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPCSTR pszFilter,
MDCS_ST_ScopeAttribute stMeasurement,
MDCS_GetDBResultsCCallback *pResultCallback
);
```

Purpose

To calculate Z' on the assay and scope attribute.

Parameters

```
hHandle - database connection handle
lAssayID - assay ID
pszFilter - filter
stMeasurement - measurement type to use in calculation
pResultCallback - pointer to a callback function to get data
```

Output

```
Recordset with fields:
Z_PRIME - calculated Z'
PLATE_ID - plate ID
```

Return

FALSE - if error occurred.

MDCS_ASSAY_GetUniqueMeasurementValues

```
BOOL MDCS_ASSAY_GetUniqueMeasurementValues (
HDBHANDLE hHandle,
LONGLONG* larrAssayIDs,
const AxStringArray& sarrFilters,
const MDCS_ST_ScopeAttribute& stMeasurement,
MDCS_GetDBResultsCCallback *pResultCallback
);
```

Purpose

To get unique values of measurement type for array of assays.

Parameters

hHandle - database connection handle
larrAssayIDs - assay IDs
sarrFilters - array of assay filters
stMeasurement - measurement to use in calculation
pResultCallback - pointer to a callback function to get data

Output

Recordset with fields: first column contains results the name is equal to the database column name of the measurement..

Return

FALSE - if error occurred

MDCS_ASSAY_FindMeasurementValues

```
BOOL MDCS_ASSAY_ ASSAY_FindMeasurementValues (
HDBHANDLE hHandle,

LONGLONG lAssayID,

LPCSTR pszFilter,

const MDCS_ST_ScopeAttribute* pstMeasurement,

const AxStringArray& arrValues,

MDCS_ST_ScopeAttribute *paarResultsColumns,

int nArrSize,

MDCS_GetDBResultsCCallback *pResultCallback
);
```

Purpose

To finds values from provided array in a measurement set.

Parameters

```
hHandle - database connection handle

lAssayIDs - measurement set to search on

pszFilter - measurement filter

pstMeasurement - measurement column to search on

arrValues - array of values to be searched

nArrSize - number of elements in array

paarResultsColumns - array of result columns paarResultsColumns

pResultCallback - pointer to a callback function to get data
```

Output

Recordset with fields: first column contains results for the column from paarResultsColumn.

Return

FALSE - if error occurred.

MDCS_ASSAY_UpdateMeasurementData

```
BOOL MDCS_ASSAY_UpdateMeasurementData (
HDBHANDLE hHandle,
LONGLONG lAssayID,
const MDCS_ST_Attribute& stColumnToUpdate,
LPCSTR pszTransformCrit,
LPCSTR pszWhereCrit = NULL
);
```

Purpose

To update measurement data based on transformation criteria.

Parameters

```
hHandle - database connection handle

lAssayID - Assay ID

stColumnToUpdate - column to update

pszTransformCrit - transformation criteria (e.g. ((MDCS15 - 15) * 10))

measurement column to search on

pszWhereCrit - filtering criteria (for example, MDCS15 <> 0)
```

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_GetAllMSetUniqueAnnotation

```
BOOL MDCS_ASSAY_GetAllMSetUniqueAnnotation (
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get unique annotation of all measurement sets.

Parameters

hHandle - database connection handle pResultCallback - pointer to a callback function to get data

Output

Columns:

COLUMN_ID - Database column name COLUMN_NAME - Display name COLUMN_TYPE - column datatype.

Return

FALSE - if error occurred

MDCS_ASSAY_CreateAttribute

```
BOOL MDCS_ASSAY_CreateAttribute (
HDBHANDLE hHandle,
const MDCS_ST_Attribute& stAttributeIn,
MDCS_ST_Attribute& stAttributeOut
);
```

Purpose

To create an Assay attribute.

Parameters

hHandle - database connection handlestAttributeIn - attribute description

Output

stAttributeOut - attribute created

Return

FALSE - if error occurred

MDCS_ASSAY_GetAttributeInfoByDisplayName

```
BOOL MDCS_ASSAY_GetAttributeInfoByDisplayName (
HDBHANDLE hHandle,

LPCSTR pszDisplayName,

MDCS_ST_Attribute& stAttributeOut
);
```

Purpose

To get attribute information using its display name.

Parameters

hHandle - database connection handle
pszDisplayName - attribute display name

Output

stAttributeOut - structure that contains attribute information

Return

FALSE - if error occurred

MDCS_ASSAY_GetAttributeInfoByDBName

```
BOOL MDCS_ASSAY_GetAttributeInfoByDBName (
HDBHANDLE hHandle,

LPCSTR pszDBName,

MDCS_ST_Attribute& stAttributeOut
);
```

Purpose

To get attribute information using its internal name.

Parameters

hHandle - database connection handlepszDBName - internal name of the attribute

Output

stAttributeOut - structure contains attribute info

Return

FALSE - if error occurred

MDCS_ASSAY_GetAttributeValueByDBName

```
BOOL MDCS_ASSAY_GetAttributeValueByDBName (
HDBHANDLE hHandle,

LPCSTR pszDBName,

LONGLONG lAssayID,

AxString& strValue
);
```

Purpose

To get an attribute value using its internal name.

Parameters

```
hHandle - database connection handlepszDBName - internal name of the attributelAssayID - assay id
```

Output

strValue - attribute value

Return

FALSE - if error occurred

MDCS_ASSAY_GetAttributeValueByDisplayName

```
BOOL MDCS_ASSAY_GetAttributeValueByDisplayName (
HDBHANDLE hHandle,
LPCSTR pszDisplayName,
LONGLONG lAssayID,
AxString& strValue
);
```

Purpose

To get an attribute value using its display name.

Parameters

```
hHandle - database connection handlepszDisplayName - display name of the attributelAssayID - assay ID
```

Output

strValue - attribute value

RETURN

FALSE - if error occurred

MDCS_ASSAY_AssignAttributeValueString

```
BOOL MDCS_ASSAY_AssignAttributeValueString (
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPCSTR pszDBName,
LPCSTR pszValue
);
```

Purpose

To assign an assay attribute value using its internal name.

Parameters

hHandle - database connection handle
 pszDBName - internal name of attribute
 lAssayID - Assay ID of the assay to change the attribute value for
 pszValue - value to assign to the attribute

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_AssignAttributeValueLong

```
BOOL MDCS_ASSAY_AssignAttributeValueLong (
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPCSTR pszDBName,
LONGLONG* plValue
);
```

Purpose

To assign an assay attribute value using its internal name.

Parameters

hHandle - database connection handle
pszDBName - internal name of the attribute
lAssayID - Assay ID of the assay to change the attribute value for
plValue - long value to assign to the attribute

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_AssignAttributeValueFloat

```
BOOL MDCS_ASSAY_AssignAttributeValueFloat (
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPCSTR pszDBName,
float* pfValue
);
```

Purpose

To an assign assay attribute value using its internal name.

Parameters

```
    hHandle - database connection handle
    pszDBName - internal name of the attribute
    lAssayID - Assay ID of the assay to change the attribute value for
    pfValue - float value to assign to the attribute
```

OutPUt

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_RenameAttribute

```
BOOL MDCS_ASSAY_RenameAttribute (
HDBHANDLE hHandle,
const MDCS_ST_Attribute& stAttrribute,
LPCSTR pszNewName
);
```

Purpose

To rename an assay attribute.

Parameters

```
hHandle - database connection handlestAttrribute - attribute to modifypszNewName - new name of the attribute
```

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_DeleteAttribute

```
BOOL MDCS_ASSAY_DeleteAttribute (
HDBHANDLE hHandle,

LPCSTR pszDBName
);
```

Purpose

To delete an assay attribute using its internal name.

Parameters

hHandle - database connection handlepszDBName - internal name of the attribute to be deleted

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_GetHeaderAndFileInfo

```
BOOL MDCS_ASSAY_GetHeaderAndFileInfo (
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_GetDBResultsCCallback* pCallback
);
```

Purpose

To get header information and file import location of a measurement set.

Parameters

hHandle - database connection handlelAssayID - Assay ID for the assay to get header and file information for pResultCallback - pointer to a callback function to get data

Output

```
Return columns:

FILE_PATH - file path

FILE_COMPUTER_NAME - name of the computer where the file was imported from.

FILE_NAME - file name

HEADER_INFO - header file information
```

Return

FALSE - if error occurred

MDCS_ASSAY_UpdateMeasurementSetName

```
BOOL MDCS_ASSAY_UpdateMeasurementSetName (
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPCSTR pszNewName
);
```

Purpose

To update the name of a measurement set.

Parameters

hHandle - database connection handlelAssayID - Assay ID of the assay for which to update the namepszNewName - new name

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_UpdateMeasurementSetDescription

```
BOOL MDCS_ASSAY_UpdateMeasurementSetDescription(
HDBHANDLE hHandle,
LONGLONG lAssayID,
LPCSTR pszNewDesc
);
```

Purpose

To update the description of a measurement set.

Parameters

```
hHandle - database connection handlelAssayID - Assay ID of the assay to update the description for pszNewDesc - new description
```

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAY_GetUniqueMeasurementValuesByPlate

```
BOOL MDCS_ASSAY_GetUniqueMeasurementValuesByPlate (
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_ST_ScopeAttribute* arrMeasurements,
int nArraySize,
MDCS_GetDBResultsCCallback *pResultCallback
);
```

Purpose

To get unique values of measurement types for a plate.

Parameters

```
    hHandle - database connection handle
    IPlateID - plate ID
    arrMeasurements - array of measurement data types
    nArraySize - array size
    pResultCallback - pointer to a callback function to get data
```

Output

Returns a recordset with fields corresponding to the measurement types in the array, ordered by the first column.

Return

FALSE - if error occurred



Note: Works only on cell attributes.

Common Database Functions

This chapter contains common database functions.

Table 8-1: Common Database Functions

Function Name	Description
MDCS_DATABASE_Optimize	To optimize the current database
MDCS_DATABASE_Compact	To compact a database. Available only for SQL server
MDCS_DATABASE_CountActiveConnections	To count active connections
MDCS_DATABASE_GetSize	To get the size of the database (in bytes). Available only for SQL server
MDCS_DATABASE_GetAvailableDatabases	To get available databases.
MDCS_DATABASE_GetVersion	To get the database version
MDCS_UTILS_RemoveMarkedData	To remove all data that was marked as deleted
MDCS_UTILS_RemoveMarkedDataEx	To remove all data that was marked as deleted
MDCS_UTILS_DropTable	To drop a table
MDCS_UTILS_DoesObjectExist	To check if an object exists in the database
MDCS_UTILS_Execute	To execute a statement that does not return any value
MDCS_UTILS_CreateTable	To create a table
MDCS_UTILS_AddColumnToTable	To add a column to a table
MDCS_UTILS_CreateForeignKey	To create a foreign key
MDCS_UTILS_CreateStoredProc	To create a stored procedure
MDCS_UPD_CreateHistoryRecord	To create a new history record
MDCS_UPD_UpdateHistoryRecordStatus	To update the update history status field
MDCS_UPD_GetHistoryRecord	To get a history record
MDCS_UPD_FindFinishedUpdates	To find a finished update
MDCS_UPD_UpdateDatabaseVersion	To update the DB_VERSIONS table

MDCS_DATABASE_Optimize

```
BOOL MDCS_DATABASE_Optimize(
HDBHANDLE hHandle,

MDCS_ProgressCallback * pCallBack
BOOL bQuickOpt = FALSE
);
```

Purpose

Function to optimize current database.

Parameters

hHandle - database connection handlepCallBack - progress callbackbQuickOpt - if TRUE, performs the quick optimization that does not include any dynamic tables

Return

FALSE - if error occurred

MDCS_DATABASE_Compact

```
BOOL MDCS_DATABASE_Compact(
HDBHANDLE hHandle,

LPCSTR pszDatabaseName,

UINT uPrecentLeave,

MDCS_ProgressCallback * pCallBack = NULL
);
```

Purpose

Function to compact a database. Available only for SQL Server.

Parameters

```
    hHandle - database connection handle
    pszDatabaseName - name of the database to compact
    uPrecentLeave - percentage of free space that should be left after compacting
    pCallback - progress callback
```

Return

FALSE - if error occurred

MDCS_DATABASE_CountActiveConnections

```
BOOL MDCS_DATABASE_CountActiveConnections(
HDBHANDLE hHandle,
LONGLONG& lConnectionFound,
LPCSTR pszDatabaseName = NULL
);
```

Purpose

Function to count active connections. Available only for SQL Server.

Parameters

hHandle - database connection handle
pszDatabaseName - name of the database

Output

IConnectionFound - number of connections found

Return

FALSE - if error occurred

MDCS_DATABASE_GetSize

```
BOOL MDCS_DATABASE_GetSize(
HDBHANDLE hHandle,

LPCSTR pszDatabaseName,

LONGLONG& lSize
);
```

Purpose

Function to get the size of the database (in bytes). Available only for SQL Server.

Parameters

```
hHandle - database connection handle
pszDatabaseName - name of the database
```

Output

ISize - size of the database

Return

FALSE - if error occurred

MDCS_DATABASE_GetAvailableDatabases

```
BOOL MDCS_DATABASE_GetAvailableDatabases(
HDBHANDLE hHandle,
AxStringArray& arrDBs,
AxStringArray* parrFilterTables = NULL
);
```

Purpose

Function to get available databases.

Parameters

hHandle - database connection handlepaarFilterTables - tables that should exist in the database

Output

arrDBS - list of databases found on server

Return

FALSE - if error occurred

MDCS_DATABASE_GetVersion

```
BOOL MDCS_DATABASE_GetVersion(
HDBHANDLE hHandle,
MDCS_ST_DBVersion& stDBVersion
);
```

Purpose

Function to get the database version.

Parameters

hHandle - database connection handle

Output

stDBVersion - database version

Return

FALSE - if error occurred

MDCS_UTILS_RemoveMarkedData

```
BOOL MDCS_UTILS_RemoveMarkedData(
HDBHANDLE hHandle,
MDCS_ProgressCallback * pCallBack
);
```

Purpose

Function to remove all data that was marked as deleted. Only members of admin group can run this function.

Parameters

```
hHandle - database connection handlepCallBack - progress callback
```

Return

FALSE - if error occurred

MDCS_UTILS_RemoveMarkedDataEx

```
BOOL MDCS_UTILS_RemoveMarkedDataEX(
HDBHANDLE hHandle,

MDCS_ProgressCallback * pCallBack

BOOL bRemovePlateByPlate = TRUE
);
```

Purpose

Function to remove all data that was marked as deleted. Only members of the admin group can run this function.

Parameters

```
    hHandle - database connection handle
    pCallBack - progress callback
    bRemovePlateByPlate - if TRUE, deleted plate by plate, instead of table by table
```

Return

FALSE - if error occurred

MDCS_UTILS_DropTable

```
BOOL MDCS_UTILS_DropTable(
HDBHANDLE hHandle,

LPCSTR pszTableName,

BOOL bTempSpace = FALSE
);
```

Purpose

Function to drop a table.

Parameters

hHandle - database connection handlepszTableName - table namebTempSpace - if TRUE, the table will be dropped from the temp space

Return

FALSE - if error occurred

MDCS_UTILS_DoesObjectExist

```
BOOL MDCS_UTILS_DoesObjectExist (
HDBHANDLE hHandle,

LPCSTR pszObjName,

BOOL& bExists,

BOOL bObjTemp = FALSE
);
```

Purpose

Function to check if an object exists in database.

Parameters

```
hHandle -database connection handlepszObjName - name of the object (table, view)bObjTemp - if TRUE, the object is temporary
```

Output

```
bExists - TRUE, if object exists
```

Return

FALSE - if error occurred

MDCS_UTILS_Execute

```
BOOL MDCS_UTILS_Execute (
HDBHANDLE hHandle,

LPCSTR pszQuery
);
```

Purpose

Function to execute a statement that does not return any value.

Parameters

```
hHandle -database connection handle
pszQuery - query string
```

Return

FALSE - if error occurred

MDCS_UTILS_CreateTable

```
BOOL MDCS_UTILS_CreateTable (
HDBHANDLE hHandle,

LPCSTR pszTableName,
const MDCS_ST_DBColun* arrColumns,
int nArraySize,

LPCTSTR pszPrefix = Null
);
```

Purpose

Function to create a table.

Parameters

hHandle -database connection handle pszTableName - name of the table nArraySize - size of array of columns arrColumns - array of table columns pszPrefix - table prefix

Return

FALSE - if error occurred

MDCS_UTILS_AddColumnToTable

```
BOOL MDCS_UTILS_AddColumnToTable (
HDBHANDLE hHandle,
BOOL& bExists,

LPCSTR strTableName,

LPCSTR strColumnName,

const MDCS_E_ColumnType& eColumnType,
int nColumnLength = 0,
BOOL bNull = TRUE,
BOOL bTempTable = FALSE
);
```

Purpose

Function to add a column to a table.

Parameters

```
hHandle -database connection handle
strTableName - name of the table
strColumnName - name of the column
nColumnLength - length of the column
bNull - if TRUE - Null allowed
bTempTable - if TRUE, table is a temporary table
eColumnType - Column type
```

Return

FALSE - if error occurred

MDCS_UTILS_CreateForeignKey

```
BOOL MDCS_UTILS_CreateForeignKey (
HDBHANDLE hHandle,

LPCTSTR pszTableName,

LPCTSTR pszColumnName,

LPCTSTR pszIndexName,

LPCTSTR pszForeignKey,

LPCTSTR pszForeignKeytable,

BOOL bTempTable=FALSE);
);
```

Purpose

Function to create a foreign key.

146 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

```
hHandle -database connection handle
pszTableName - table name
pszColumnName - column name
pszIndexName - index name
pszForeignKey - foreign key column
pszForeignKeyTable - Foreign key table (table in which pszForeignKey is a primary key)
bTempTable - if TRUE, table is a temporary table
```

Return

FALSE - if error occurred

MDCS_UTILS_CreateStoredProc

```
BOOL MDCS_UTILS_CreateStoredProc (
HDBHANDLE hHandle,

LPCTSTR pszName,

LPCTSTR pszSQL
);
```

Purpose

Function to create a stored procedure.

Parameters

```
hHandle - database connection handlepszName - stored procedure namepszSQL - SQL string that defines the stored procedure
```

Return

FALSE - if error occurred

MDCS_UPD_CreateHistoryRecord

```
BOOL MDCS_UPD_CreateHistoryRecord (
HDBHANDLE hHandle,
const MDCS_ST_Update& stUpdate,
LONGLONG& lUpdateID
);
```

Purpose

Function to create a new history record.

```
hHandle - database connection handle
stUpdate - History update structure
IUpdateID - update ID
```

Return

FALSE - if error occurred

MDCS_UPD_UpdateHistoryRecordStatus

```
BOOL MDCS_UPD_UpdateHistoryRecordStatus (
HDBHANDLE hHandle,
MDCS_ST_Update::MDCS_E_UpdateStatus eStatus,
LONGLONG lUHID
);
```

Purpose

Function to update the status field for a history record.

Parameters

```
hHandle - database connection handleeStatus - statusIUHID - update history ID
```

Return

FALSE - if error occurred

MDCS_UPD_GetHistoryRecord

```
BOOL MDCS_UPD_GetHistoryRecord (
(HDBHANDLE hHandle,
LONGLONG lUHID,
MDCS_GetDBResultsCCallback* pResultCallback);
);
```

Purpose

Function to get a history record for an update history ID.

hHandle - database connection handleIUHID - update history IDpResultCallback - pointer to callback function to get data

Output

The output history record will contain columns from the UPDATE_HISTORY table.

Return

FALSE - if error occurred

MDCS_UPD_FindFinishedUpdates

```
BOOL MDCS_UPD_FindFinishedUpdates
(HDBHANDLE hHandle,
LONGLONG lUHID,
MDCS_ST_Update::MDCS_E_SystemType eSType,
MDCS GetDBResultsCCallback* pResultCallback)
```

Purpose

Function to find a finished update that matches the input update ID and type.

Parameters

```
hHandle - database connection handleIUHID - update history IDeSType - update typepResultCallback - pointer to callback function to get data
```

Output

The output history record will contain columns from the UPDATE HISTORY table.

Return

FALSE - if error occurred

MDCS_UPD_UpdateDatabaseVersion

```
BOOL MDCS_UPD_UpdateDatabaseVersion(
HDBHANDLE hHandle,
LONGLONG lDBVerID,
LONGLONG lDBVerPara1,
LONGLONG lDBVerPara2,
LONGLONG lDBVerPara3
);
```

Purpose

Function to update the DB_Versions table.

Parameters

hHandle - database connection handle IDBVerID - database version ID IDBVerPara1 - database version leading number IDBVerPara2 - database version middle number IDBVerPara3 - database version trailing number

Return

FALSE - if error occurred

Plate Functions 9

This chapter contains functions that you can use to work with plates.

Table 9-1: Plate Functions

Function Name	Description
MDCS_PLATE_CreatePlate	Create a new plate and set default access permissions.
MDCS_PLATE_GetInfo	To get plate information
MDCS_PLATE_GetAllPropertyAttributes	Will return description of columns for the plate attributes
MDCS_PLATE_GetAllOrderedByAttributes	Return all plates ordered by attributes
MDCS_PLATE_GetInfoBasedOnAssay	To get plate information for a measurement set
MDCS_PLATE_Delete	To delete a plate
MDCS_PLATE_DeleteImages	To delete plate images
MDCS_PLATE_GetUniqueAttributeValues	To get unique values for a plate attribute
MDCS_PLATE_GetAllByAttributes	To get plates based on specified attributes values
MDCS_PLATE_ManageSecurity	To call a dialog to manage security access for a plate
MDCS_PLATE_FolderManageSecurity	To call a dialog to manage security access for plates in a dynamic folder
MDCS_PLATE_ShareSiteImagesByTimeAndZInd ex	To share images by site, timepoint and z index
MDCS_PLATE_CreateSite	To create a new site
MDCS_PLATE_CreateSeries	To create a new series
MDCS_PLATE_CreateImageSource	To create an image source
MDCS_PLATE_CreateImageRecord	To create an image record
MDCS_PLATE_ImportLayoutData	To create (import) a plate layout
MDCS_PLATE_ApplyLayoutToAssay	To apply a layout to an array
MDCS_PLATE_DeleteLayout	To delete a plate layout
MDCS_PLATE_GetAllRecords	To get all plate records
MDCS_PLATE_GetRecord	To get a single plate record
MDCS_PLATE_GetSiteRecord	To get a site record
MDCS_PLATE_GetSitesByPlate	To get all sites for a plate
MDCS_PLATE_GetImageRecord	To get an image record
MDCS_PLATE_GetImageRecordPerPlate	To get all image records for a plate
MDCS_PLATE_UpdateImageObjectID	To assign an image object ID to the plate image

Table 9-1: Plate Functions (cont'd)

MDCS_PLATE_GetImageObjectIDForImage	To get the ID of the object record where the image is stored for the plate image
MDCS_PLATE_GetSeriesRecord	To get a series record by ID
MDCS_PLATE_GetSiteLocationsForPlate	To get all site positions for a plate
MDCS_PLATE_GetImageSourcesOfPlate	To get the image source records of a plate
MDCS_PLATE_GetMaxTimePointForPlate	To get the maximum the timepoint value for a plate
MDCS_PLATE_GetSeriesIDAtZAndT	To retrieve the series ID based on the site, Z index and T index
MDCS_PLATE_GetPlatesByDate	To get plates information by date range
MDCS_PLATE_UpdateAcquisition	To update an acquisition
MDCS_PLATE_CreateAttribute	To create a plate attribute
MDCS_PLATE_GetAttributeInfoByDisplayName	To get attribute information by display name
MDCS_PLATE_GetAttributeInfoByDBName	To get attribute information by database name
MDCS_PLATE_AssignAttributeValueString	To assign a plate attribute value(string)
MDCS_PLATE_AssignAttributeValueLong	To assign a plate attribute value(Long)
MDCS_PLATE_AssignAttributeValueFloat	To assign a plate attribute value(float)
MDCS_PLATE_UpdateInfo	To update plate information
MDCS_PLATE_GetInfoBasedOnAssaySet	To get plate information for a set of assays.
MDCS_PLATE_CreateLayout	To create a plate layout
MDCS_PLATE_ApplyLayoutToPlate	To apply a layout to a plate
MDCS_PLATE_GetTemplate	To get the annotation template using the id
MDCS_PLATE_GetTemplateByName	To get the annotation template using the name
MDCS_PLATE_GetTemplates	To get all plate annotation templates
MDCS_PLATE_UpdatePlateTemplate	To update plate annotation templates
MDCS_PLATE_RenameAttribute	To rename a plate attribute
MDCS_PLATE_CountPlateDatasets	To get the number of datasets where plate assays are used
MDCS_PLATE_GetAcqSiteCount	To get a count of all acquired sites per plate
MDCS_PLATE_GetAcqWellCount	To get a count of all acquired wells per plate
MDCS_PLATE_GetAcqSeriesCount	To get a count of all acquired series per plate

152 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Table 9-1: Plate Functions (cont'd)

MDCS_PLATE_GetCompoundCount	To get a count of compounds per plate
MDCS_PLATE_GetControlsCount	To get a count of controls per plate
MDCS_PLATE_GetControlStatistic	To get the control statistic per plate
MDCS_PLATE_GetCompleteImageInfo	To get complete image information for a plate
MDCS_PLATE_CanModify	To check if the current user can modify a plate
MDCS_PLATE_GetImageIDs	To get the IDs of plate images
MDCS_PLATE_GetThumbImageIDs	To get the IDs of plate thumbnails
MDCS_PLATE_ChangeStatus	To change the plate status

MDCS_PLATE_CreatePlate

```
BOOL MDCS_PLATE_CreatePlate(
HDBHANDLE hHandle,
const MDCS_ST_PlateInfo& stPlateInfoIn,
MDCS_ST_PlateInfo* pstPlateInfoOut,
LONGLONG* plDefaultGroupID = NULL
);
```

Purpose

Function creates a new plate and sets default access permissions. Permissions are granted to administrators and the current user only.

Parameters

hHandle - database connection handlestPlateInfoIn - plate description to insertplDefaultGroupID - ID of a group that also has Lab Head permissions for the plate (If NULL - this will be ignored)

Output

pstPlateInfoOut - plate description as created in database

Return

FALSE - if error occurred

MDCS_PLATE_GetInfo

```
BOOL MDCS_PLATE_GetInfo(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_ST_PlateInfo& stPlateInfoOut
);
```

Purpose

Function to get plate information.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output

pstPlateInfoOut - plate description as created in the database

Return

FALSE - if error occurred

MDCS_PLATE_GetAllPropertyAttributes

```
BOOL MDCS_PLATE_GetAllPropertyAttributes(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function returns description of columns for the plate attributes.

Parameters

hHandle - database connection handle pResultCallback - pointer to a callback function to get data

Output

Result columns are:

 $\ensuremath{\mathsf{ATTR}}\xspace_{\mathsf{NAME}}$ - name of the column that contains data for the attribute in the database, string

ATTR_NAME - display name of the attribute, string ATTR_FORMAT - format of the data, integer cast to

MDCS_E_ColumnType

ATTR_TYPE - type of the property, integer cast to MDCS_E_AttributeType

Return

FALSE - if error occurred

154 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

MDCS_PLATE_GetAllOrderedByAttributes

```
BOOL MDCS_PLATE_GetAllOrderedByAttributes(
HDBHANDLE hHandle,
const AxStringArray* pAXStringArray,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function returns all plates ordered by attributes (how the attributes are ordered in the array).

Parameters

hHandle - database connection handlepAXStringArray - pointer to string array of property column namespResultCallback - pointer to a callback function to get data

Output

```
Result columns are:
PLATE_ID - plate ID
PLATE_NAME -plate name
ACQ_ID - acquisition ID
GLOBAL_ID - plate global ID
CREATOR_ID creator ID
CREATOR_NAME - creator name
BARCODE - Plate barcode
TIME_CREATED - time when plate was created in seconds from 01/01/1970 00:00:00
ACQ_NAME - acquisition name
```

Return

FALSE - if error occurred

MDCS_PLATE_GetInfoBasedOnAssay

```
BOOL MDCS_PLATE_GetInfoBasedOnAssay(
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_ST_PlateInfo& stPlateInfoOut
);
```

Purpose

Function to get plate information for an assay.



Note: Plate information may not exist for a measurement set.

Parameters

hHandle - database connection handle IAssayID - assay ID for which plate information will be retrieved

Output

stPlateInfoOut - plate information

Return

FALSE - if error occurred

MDCS_PLATE_Delete

```
BOOL MDCS_PLATE_Delete(
HDBHANDLE hHandle,
LONGLONG lPlateID
);
```

Purpose

Function to delete a plate.



Note: Will not delete actual data; will only mark the data as deleted.

Parameters

hHandle - database connection handleIPlateID - plate ID for the plate to be deleted

Return

FALSE - if error occurred

MDCS_PLATE_DeleteImages

```
BOOL MDCS_PLATE_DeleteImages(
HDBHANDLE hHandle,
LONGLONG lPlateID
);
```

Purpose

Function to delete plate images.



Note: Will not delete actual images; will only mark the images as deleted.

Parameters

hHandle - database connection handleIPlateID - plate ID for the plate images to be deleted

Return

FALSE - if error occurred

MDCS_PLATE_GetUniqueAttributeValues

```
BOOL MDCS_PLATE_GetUniqueAttributeValues(
HDBHANDLE hHandle,

LPCSTR pszAttrColumnName,

const AxStringArray* pAXAttrColumnNames,

const AxStringArray* pAXAttrValues,

AxStringArray* pAXValuesOut
);
```

Purpose

Function to get unique values for a plate attribute.

Parameters

hHandle - database connection handle pszAttrColumnName - column that will contain unique values pAXAttrColumnNames - attribute column names to query on pAXAttrValues - values that attributes are equal to

Output

pAXValuesOut - array of unique values

Return

FALSE - if error occurred

MDCS_PLATE_GetAllByAttributes

```
BOOL MDCS_PLATE_GetAllByAttributes(
HDBHANDLE hHandle,
const AxStringArray* pAXDisplayColumns,
const AxStringArray* pAXAttrColumn,
const AxStringArray* pAXAttrValues,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get plates based on specified attribute values.

Parameters

```
HHandle - database connection handle

pAXDisplayColumns - columns to retrieve

pAXAttrColumn - columns to query on

pAXAttrValues - values corresponding to the columns to query on
```

Output

```
Result columns will be:
PLATE_ID - plate ID
PLATE_NAME - plate name
ACQ_ID - acquisition ID
GLOBAL_ID - plate global ID)
CREATOR_ID - plate creator ID
CREATOR_NAME - name of a person who created a plate (string)
BARCODE - plate barcode)
TIME_CREATED - time when plate was created in seconds from 01/01/1970 00:00:00
ACQ_NAME - acquisition name that was used to acquire plate and columns specified in pAXDisplayColumns
```

MDCS_PLATE_ManageSecurity

```
BOOL MDCS_PLATE_ManageSecurity(
HDBHANDLE hHandle,
LONGLONG lPlateID,
HWND hWnd = NULL,
LPCSTR pszDlgTitle = NULL
);
```

Purpose

Calls a dialog to manage security access to the plate.

158 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

```
hHandle - database connection handleLPlateID - plate IDpszDlgTitle - dialog titlehWnd - handle to the application window
```

Return

FALSE - if dialog cancelled

MDCS_PLATE_FolderManageSecurity

```
BOOL MDCS_PLATE_FolderManageSecurity(
HDBHANDLE hHandle,
const AxStringArray& arrAttributes,
const AxStringArray& arrValues,
HWND hWnd = NULL,
LPCSTR pszDlgTitle = NULL
);
```

Purpose

Calls a dialog to manage security access for plates in a dynamic folder.

Parameters

```
hHandle - database connection handle arrAttributes - array of attributes arrValues - array of attribute values pszDlgTitle - dialog title hWnd - handle to the application window
```

Return

FALSE - if dialog cancelled

MDCS_PLATE_ShareSiteImagesByTimeAndZIndex

```
BOOL MDCS_PLATE_ShareSiteImagesByTimeAndZIndex(
HDBHANDLE hHandle,

LONGLONG lSiteID,

int nTIndexIn,

int nZIndexIn,

LONGLONG lSiteIDNew,

int nTIndexNew,

int nZIndexNew,

LONGLONG lSelectBySourceID = 0
);
```

Purpose

To share images by site, timepoint and z index.

Parameters

```
hHandle - database connection handle

ISiteID - site ID to share

nTIndexIn - time index that should be copied

nZIndexIn - Z Index that should be copied

ISelectBySourceID - source ID, if needed to be part of selection
(If 0 - ignored)

ISiteIDNew - destination site ID

nTIndexNew - destination time index

nZIndexNew - destination Z Index
```

Output

ICountUpdated - number of record copied

Return

FALSE - if it fails

MDCS_PLATE_CreateSite

```
BOOL MDCS_PLATE_CreateSite(
HDBHANDLE hHandle,
const MDCS_ST_Site& stSite,
LONGLONG* plSiteID
);
```

Purpose

Function to create a new site.

160 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

```
hHandle - database connection handle
stSite - site description
```

Output

```
plSiteID - site ID
```

Return

FALSE - if error occurred

MDCS_PLATE_CreateSeries

```
BOOL MDCS_PLATE_CreateSeries(
HDBHANDLE hHandle,
const MDCS_ST_SeriesInfo& stSeries,
LONGLONG* plSeriesID
);
```

Purpose

Function to create a new series.

Parameters

```
hHandle - database connection handlestSeries - series description
```

Output

plSeriesID - newly created series ID

Return

FALSE - if error occurred

MDCS_PLATE_CreateImageSource

```
BOOL MDCS_PLATE_CreateImageSource(
HDBHANDLE hHandle,
const MDCS_ST_ImageSource& stImageSource,
LONGLONG* lImageSourceID
);
```

Purpose

Function to create image source.

hHandle - database connection handle
stImageSource - image source description

Output

IImageSourceID - newly created Image source ID

Return

FALSE - if error occurred

MDCS_PLATE_CreateImageRecord

```
BOOL MDCS_PLATE_CreateImageRecord (
HDBHANDLE hHandle,
const MDCS_ST_SiteImageInfo& stImageDesc,
LONGLONG* plRecordID
);
```

Purpose

Function to create an image record.

Parameters

hHandle - database connection handle
stImageDesc - image description

Output

plRecordID - newly created Image record ID

Return

FALSE - if error occurred

MDCS_PLATE_ImportLayoutData

```
BOOL MDCS_PLATE_ImportLayoutData(
HDBHANDLE hHandle,
const MDCS_ST_PlateTemplate& stPlateTemplate,
MDCS_ImportPlateLayout* pDatasource,
LONGLONG* plLayoutID,
MDCS_ProgressCallback* pCallback
);
```

Purpose

Function to import a plate layout.

```
hHandle - database connection handle pDatasource - datasource that contains layout data stPlateLayout - structure that describes plate layout pCallback - callback class
```

Output

plLayoutID - newly created plate layout ID

Return

FALSE - if error occurred

MDCS_PLATE_ApplyLayoutToAssay

```
BOOL MDCS_PLATE_ApplyLayoutToAssay(
HDBHANDLE hHandle,
LONGLONG lLayoutID,
LONGLONG lAssayID,
MDCS_ProgressCallback * pCallBack
);
```

Purpose

Function to apply a layout to an assay.

Parameters

```
hHandle - database connection handlelLayoutID - layout IDlAssayID - assay IDpCallback - callback object
```

Return

FALSE - if error occurred

MDCS_PLATE_DeleteLayout

```
BOOL MDCS_PLATE_DeleteLayout(
HDBHANDLE hHandle,
LONGLONG lLayoutID
);
```

Purpose

Function to delete plate layout.

hHandle - database connection handlelLayoutID - layout ID

Return

FALSE - if error occurred

MDCS_PLATE_GetAllRecords

```
BOOL MDCS_PLATE_GetAllRecords(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get all plate records.

Parameters

hHandle - database connection handlepResultCallback - pointer to a callback function to get data

Output

Returns records with columns
ID as plate ID
BATCH_ID as batch ID
X_WELLS, Y_WELLS number of wells
BARCODE - barcode
GLOBAL_ID - plate global ID
STORAGE_LOCATION - storage location
STORAGE_TYPE - storage type
PLATE_NAME - plate name
PLATE_DESCRIPTION - plate description
CREATOR - name of user who created the batch
TIME_CREATED - time when plate was created
NAME - plate name

Return

FALSE - if error occurred



Note: Storage locations always set to "database"; the DLL will decide where the images are located.

MDCS_PLATE_GetRecord

```
BOOL MDCS_PLATE_GetRecord(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get a single plate record.

Parameters

hHandle - database connection handleIPlateID - plate IDpResultCallback - pointer to a callback function to get data

Output

Returns records with columns: ID as plate ID

BATCH_ID as batch ID

X WELLS, Y WELLS number of wells

BARCODE - barcode

GLOBAL_ID - plate global ID

STORAGE_LOCATION - storage location

STORAGE_TYPE - storage type

PLATE_NAME - plate name

PLATE_DESCRIPTION - plate description

CREATOR - name of user who create the batch

TIME_CREATED - time when plate was created

NAME - plate name

Return

FALSE - if error occurred



Note: Storage locations will always set to "database", the DLL will decide where the images are located.

MDCS_PLATE_GetSiteRecord

```
BOOL MDCS_PLATE_GetSiteRecord(
HDBHANDLE hHandle,
LONGLONG lSiteID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get a site record.

Parameters

```
hHandle - database connection handle

*ISiteID - site ID

*pResultCallback - pointer to a callback function to get data*
```

Output

```
Returns records with columns:

ID as plate ID

BATCH_ID as batch ID

X_WELLS, (number of wells), Y_WELLS (number of wells)

X_POSITION, (position within a well), Y_POSITION (position within a well)

TO_DELETE - (delete flag)
```

Return

FALSE - if error occurred

MDCS_PLATE_GetSitesByPlate

```
BOOL MDCS_PLATE_GetSitesByPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get all sites for a plate.

```
hHandle - database connection handle

IPlateID - plate ID

pResultCallback - pointer to a callback function to get data
```

Output

```
Returns records with columns:

ID as plate ID

BATCH_ID as batch ID

X_WELLS, (number of wells), Y_WELLS (number of wells)

X_POSITION, (position within a well), Y_POSITION (position within a well)

TO_DELETE - (delete flag)
```

Return

FALSE - if error occurred

MDCS_PLATE_GetImageRecord

```
BOOL MDCS_PLATE_GetImageRecord(
HDBHANDLE hHandle,
LONGLONG lImageID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get an image record.

```
hHandle - database connection handlelImageID - image IDpResultCallback - pointer to a callback function to get data
```

Output

```
Returns records with columns:

ID as record ID

IMAGE_DATA_ID

SITE_ID

IMAGE_SOURCE_ID

ACQUIRE_DATE

TIMEPOINT

Z_POSITION

Z_STEP

SERIES_INFORMATION_ID

SOURCE_DESCRIPTION

STORAGE_LOCATION_NAME (location ID)

STORAGE_TYPE_NAME (always "DATABASE")
```

Return

FALSE - if error occurred

MDCS_PLATE_GetImageRecordPerPlate

```
BOOL MDCS_PLATE_GetImageRecordPerPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get all image records for a plate.

```
hHandle - database connection handleIPlateID - plate IDpResultCallback - pointer to a callback function to get data
```

Output

```
Returns records with columns:
ID as record ID
IMAGE_DATA_ID
SITE_ID
IMAGE_SOURCE_ID
ACQUIRE_DATE
TIMEPOINT
Z_POSITION
Z_STEP
SERIES_INFORMATION_ID
SOURCE_DESCRIPTION
STORAGE_LOCATION_NAME (location ID)
```

STORAGE_TYPE_NAME ("DATABASE")

Return

FALSE - if error occurred

MDCS_PLATE_UpdateImageObjectID

```
BOOL MDCS_PLATE_UpdateImageObjectID(
HDBHANDLE hHandle,
LONGLONG lPlateImageID,
LONGLONG lImageObjectID
);
```

Purpose

Function to assign image object ID to the plate image.

Parameters

```
hHandle - database connection handleIPlateImageID - plate image IDIImageObjectID - image object ID
```

Return

FALSE - if error occurred

MDCS_PLATE_GetImageObjectIDForImage

```
BOOL MDCS_PLATE_GetImageObjectIDForImage(
HDBHANDLE hHandle,
LONGLONG lPlateImageID,
LONGLONG& lImageObjectID
);
```

Purpose

To get the image object ID (ID of the image record) for a plate image.

Parameters

hHandle - database connection handleIPlateImageID - plate image ID

Output

IImageObjectID - image object ID

Return

FALSE - if error occurred

MDCS_PLATE_GetSeriesRecord

```
BOOL MDCS_PLATE_GetSeriesRecord(
HDBHANDLE hHandle,
LONGLONG lSeriesID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get a series record by ID.

Parameters

```
hHandle - database connection handleISeriesID - series IDpResultCallback - pointer to a callback function to get data
```

Output

```
Returns records with columns:
ID as series ID
all other columns from the SERIES_INFO table
```

Return

FALSE - if error occurred

MDCS_PLATE_GetSiteLocationsForPlate

```
BOOL MDCS_PLATE_GetSiteLocationsForPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

Function to get all site positions for a plate.

Parameters

```
hHandle - database connection handleIPlateID - Plate IDpResultCallback - pointer to a callback function to get data
```

Output

Returns records with columns X_POSITION, Y_POSITION values

Return

FALSE - if error occurred

MDCS_PLATE_GetImageSourcesOfPlate

```
BOOL MDCS_PLATE_GetImageSourcesOfPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get image source records of a plate.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output columns

```
ID as image source ID SOURCE_DESCRIPTION as source description SOURCE_ILLUMINATION as source illumination
```

Return

FALSE - if error occurred

MDCS_PLATE_GetMaxTimePointForPlate

```
BOOL MDCS_PLATE_GetMaxTimePointForPlate(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG& lMaxTimepoint
);
```

Purpose

To get the maximum value for a timepoint for a plate.

Parameters

IPlateID - plate ID

Output

IMaxTimepoint - maximum timepoint value
hHandle - database connection handle

Return

FALSE - if error occurred

MDCS_PLATE_GetSeriesIDAtZAndT

```
MDCS_PLATE_GetSeriesIDAtZAndT(
HDBHANDLE hHandle,
LONGLONG lSiteID,
LONG lZIndex,
LONG lTIndex,
LONGLONG& lSeriesID
);
```

Purpose

To retrieve the series ID based on site, Z index and T index.

Parameters

```
hHandle - database connection handlelSiteID - site IDlZIndex - Z IndexlTIndex - Time Index
```

Output

ISeriesID - found series ID

Return

FALSE - if error occurred

172 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

MDCS_PLATE_GetPlatesByDate

```
MDCS_PLATE_GetPlatesByDate(
HDBHANDLE hHandle,
LPCSTR pszStartDate,
LPCSTR pszEndDate,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get plate information by date range.

Parameters

```
hHandle - database connection handlepszStartDate - start datepszEndDate - end datepResultCallback - pointer to a callback function to get data
```

Output

```
Result columns are:

PLATE_ID - plate ID

PLATE_NAME - plate name

ACQ_ID - acquisition ID

GLOBAL_ID - plate global ID

CREATOR_ID - creator ID

CREATOR_NAME - creator name

BARCODE - plate barcode

TIME_CREATED - time when plate was created in seconds from 01/01/1970 00:00:00
```

Return

FALSE - if error occurred

MDCS_PLATE_UpdateAcquisition

```
BOOL MDCS_PLATE_UpdateAcquisition(
HDBHANDLE hHandle,
const MDCS_ST_Acquisition& stAcq
);
```

Purpose

To update an acquisition.

hHandle - database connection handlestAcq - acquisition description

Return

FALSE - if error occurred

MDCS_PLATE_CreateAttribute

```
BOOL MDCS_PLATE_CreateAttribute(
HDBHANDLE hHandle,
const MDCS_ST_Attribute& stAttributeIn,
MDCS_ST_Attribute& stAttributeOut
);
```

Purpose

Function to create a plate attribute.

Parameters

hHandle - database connection handle stAttributeIn - attribute description

Output

stAttributeOut - attribute created

Return

FALSE - if error occurred

MDCS_PLATE_GetAttributeInfoByDisplayName

```
BOOL MDCS_PLATE_GetAttributeInfoByDisplayName(
HDBHANDLE hHandle,

LPCSTR pszDisplayName,

MDCS_ST_Attribute& stAttributeOut
);
```

Purpose

Function to get attribute information by display name.

hHandle - database connection handle
pszDisplayName - display name of attribute

Output

stAttributeOut - attribute information

Return

FALSE - if error occurred

MDCS_PLATE_GetAttributeInfoByDBName

```
BOOL MDCS_PLATE_GetAttributeInfoByDBName(
HDBHANDLE hHandle,

LPCSTR pszDBName,

MDCS_ST_Attribute& stAttributeOut
);
```

Purpose

Function to get attribute information by internal database name.

Parameters

```
hHandle - database connection handle
pszDBName - internal name of attribute
```

Output

stAttributeOut - attribute information

Return

FALSE - if error occurred

MDCS_PLATE_AssignAttributeValueString

```
BOOL MDCS_PLATE_AssignAttributeValueString(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LPCSTR pszDBName,
LPCSTR pszValue
);
```

Purpose

Function to assign a value to a plate attribute.

```
hHandle - database connection handlelPlateID - plate IDpszDBName - internal name of attributepszValue - value to assign
```

Output

None

Return

FALSE - if error occurred

MDCS_PLATE_AssignAttributeValueLong

```
BOOL MDCS_PLATE_AssignAttributeValueLong(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LPCSTR pszDBName,
LONGLONG* plValue
);
```

Purpose

Function to assign plate attribute value.

Parameters

```
hHandle - database connection handle

lPlateID - ID of a plate that needs to change attribute value

pszDBName - internal name of attribute

(MDCS_ST_Attribute::szDBcolumnName)

plValue - value
```

Output

stAttributeOut - attribute created

Return

FALSE - if error occurred

MDCS_PLATE_AssignAttributeValueFloat

```
BOOL MDCS_PLATE_AssignAttributeValueFloat(
HDBHANDLE hHandle,
LONGLONG lPlateID,
LPCSTR pszDBName,
float* pfValue
);
```

Purpose

Function to assign a value to a plate attribute.

Parameters

```
hHandle - database connection handleIPlateID - plate IDpszDBName - internal name of attributepfValue - value to assign
```

Output

None

Return

FALSE - if error occurred

MDCS_PLATE_UpdateInfo

```
BOOL MDCS_PLATE_UpdateInfo(
HDBHANDLE hHandle,
const MDCS_ST_PlateInfo& stInfoIn
);
```

Purpose

Function to update plate information.

hHandle - database connection handlestInfoIn - new information about the plate (plate ID must be specified)

Output

NONE

Return

FALSE - if error occurred

MDCS_PLATE_GetInfoBasedOnAssaySet

```
BOOL MDCS_PLATE_GetInfoBasedOnAssaySet (
HDBHANDLE hHandle,
const LONGLONG* larrAssayID,
int nElementInArray,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

Function to get plate information for a set of assays.



Note: Plate information may not exist for an assay.

hHandle - database connection handle
 larrAssayID - array of assay IDs
 nElementInArray - the number of elements in the array (larrAssayID)
 pResultCallback - pointer to a callback function to get data

Output

Result columns will be:

ASSAY_ID - assay ID

PLATE_ID - plate ID

PLATE_NAME - plate name

ACQ_ID - acquisition ID

GLOBAL_ID - plate global ID

CREATOR_ID -creator ID

CREATOR_NAME - creator name

BARCODE - plate barcode

TIME_CREATED - time when plate was created in seconds from 01/01/1970 00:00:00)

ACQ_NAME - acquisition name

Return

FALSE - if error occurred

MDCS_PLATE_CreateLayout

```
BOOL MDCS_PLATE_CreateLayout (
HDBHANDLE hHandle,
const MDCS_ST_PlateTemplate& stPlateTemplate,
MDCS_ST_PlateTemplate* pstPlateTemplateOut
);
```

Purpose

Function to create a plate layout.

Parameters

hHandle - database connection handlestPlateLayout - structure that describes plate layout

Output

pstPlateTemplateOut - created layout

Return

FALSE - if error occurred

MDCS_PLATE_ApplyLayoutToPlate

```
BOOL MDCS_PLATE_ApplyLayoutToPlate (
HDBHANDLE hHandle,
LONGLONG lLayoutID,
LONGLONG lPlateID
);
```

Purpose

Function to apply a layout to a plate.

Parameters

```
hHandle - database connection handlelLayoutID - layout IDlPlateID - plate ID
```

Output

NONE

Return

FALSE - if error occurred

MDCS_PLATE_GetTemplate

```
BOOL MDCS_PLATE_GetTemplate (
HDBHANDLE hHandle,
LONGLONG lTemplID,
MDCS_ST_PlateTemplate& stPlateTemplate
);
```

Purpose

Function to get the plate annotation template by template ID.

Parameters

```
hHandle - database connection handlelTemplID - template ID
```

Output

stPlateTemplate - plate template

Return

FALSE - if error occurred

MDCS_PLATE_GetTemplateByName

```
BOOL MDCS_PLATE_GetTemplateByName (
HDBHANDLE hHandle,

LPCSTR pszTemplName,

MDCS_ST_PlateTemplate& stPlateTemplate
);
```

Purpose

Function to get the plate annotation template by providing its name.

Parameters

hHandle - database connection handle
pszTemplName - template name

Output

stPlateTemplate - plate template

Return

FALSE - if error occurred

MDCS_PLATE_GetTemplates

```
BOOL MDCS_PLATE_GetTemplates (
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

Function to get all existing plate annotation templates.

Parameters

hHandle - database connection handle pResultCallback - pointer to a callback function to get data

Output

Callback will return columns described in MDCS_ST_PlateLayout structure

Return

FALSE - if error occurred

MDCS_PLATE_UpdatePlateTemplate

```
BOOL MDCS_PLATE_UpdatePlateTemplate (
HDBHANDLE hHandle,
const MDCS_ST_PlateTemplate& stPlateTemplateIn,
MDCS_ST_PlateTemplate* pstPlateTemplate
);
```

Purpose

Function to update a plate layout template.

Parameters

hHandle - database connection handle stPlateTemplateIn - template input

Output

stPlateTemplate - updated plate template

Return

FALSE - if error occurred

MDCS_PLATE_RenameAttribute

```
BOOL MDCS_PLATE_RenameAttribute (
HDBHANDLE hHandle,
const MDCS_ST_Attribute& stAttrribute,
LPCSTR pszNewName
);
```

Purpose

Function to rename a plate attribute.

Parameters

```
hHandle - database connection handlestAttrribute - attribute to modifypszNewName - new name of the attribute
```

Output

NONE

Return

FALSE - if error occurred

MDCS_PLATE_CountPlateDatasets

```
BOOL MDCS_PLATE_CountPlateDatasets (
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG& lNumDataset
);
```

Purpose

Function to get the number of the datasets where a plate is used.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output

INumDataset - number of the datasets where the plate is used

Return

FALSE - if error occurred

MDCS_PLATE_GetAcqSiteCount

```
BOOL MDCS_PLATE_GetAcqSiteCount (
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG& lCount
);
```

Purpose

Function to get the count of sites per plate..

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output

ICount - number of sites per plate

Return

FALSE - if error occurred

MDCS_PLATE_GetAcqWellCount

```
BOOL MDCS_PLATE_GetAcqWellCount (
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG& lCount
);
```

Purpose

Function to get the count of all acquired wells per plate.

Parameters

```
hHandle - database connection handle
IPlateID - plate ID
```

Output

ICount - number of acquired wells per plate

Return

FALSE - if error occurred

MDCS_PLATE_GetAcqSeriesCount

```
BOOL MDCS_PLATE_GetAcqSeriesCount (
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG& lCount
);
```

Purpose

Function to get the count of all acquired series per plate.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output

ICount - number of acquired series per plate

Return

FALSE - if error occurred

MDCS_PLATE_GetCompoundCount

```
BOOL MDCS_PLATE_GetCompoundCount (
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG& lCount
);
```

Purpose

Function to get the count of all compounds per plate.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output

ICount - number of compounds per plate

Return

FALSE - if error occurred

MDCS_PLATE_GetControlsCount

```
BOOL MDCS_PLATE_GetControlsCount (
HDBHANDLE hHandle,
LONGLONG lPlateID,
LONGLONG& lCount
);
```

Purpose

Function to get the count of controls per plate.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output

ICount - number of controls per plate

Return

FALSE - if error occurred

MDCS_PLATE_GetControlStatistic

```
BOOL MDCS_PLATE_GetControlStatistic (
HDBHANDLE hHandle,
LONGLONG lPlateID,
AxStringArray& arrControlNames,
AxStringArray& arrControlCount
);
```

Purpose

Function to get name and count of all controls in a plate.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output

arrControlNames - array that contains the control names
arrControlCount - array that contains the count of the corresponding
control

Return

FALSE - if error occurred

MDCS_PLATE_GetCompleteImageInfo

```
BOOL MDCS_PLATE_GetCompleteImageInfo(
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get complete image information for a plate.

```
hHandle - database connection handleIPlateID - plate IDpResultCallback - pointer to a callback function to get data
```

Output

```
Will return records with columns:
SITE ID
SERIES_ID
IMAGE_SOURCE_ID,
OBJ ID
IMAGE ID
OBJ_SIZE
OBJ SERVER NAME
OBJ_EXT, THUMB_SIZE
THUMB_SERVER_NAME
LOCATION_ID,
THUMB_LOCATION_ID
WELL_X
WELL Y
X POSITION
Y_POSITION
Z INDEX
Z_POSITION
T INDEX
SOURCE_DESCRIPTION
SOURCE_ILLUMINATION
```

Return

FALSE - if error occurred

MDCS_PLATE_CanModify

```
BOOL MDCS_PLATE_CanModify (
HDBHANDLE hHandle,
LONGLONG lPlateID,
BOOL* pbCan
);
```

Purpose

Function to check whether the current user can modify a plate.

```
hHandle - database connection handleIPlateID - plate ID
```

Output

```
pbCan - TRUE (if the plate can be modified)
```

Return

FALSE - if error occurred

MDCS_PLATE_GetImageIDs

```
BOOL MDCS_PLATE_GetImageIDs (
HDBHANDLE hHandle,
LONGLONG lPlateID,
AxStringArray& arrIDs
);
```

Purpose

Function to get the image IDs for plate images.

Parameters

```
hHandle - database connection handleIPlateID - plate ID
```

Output

arrIDs - array of Image IDs found for the plate

Return

FALSE - if error occurred

MDCS_PLATE_GetThumbImageIDs

```
BOOL MDCS_PLATE_GetThumbImageIDs (
HDBHANDLE hHandle,
LONGLONG lPlateID,
AxStringArray& arrIDs
);
```

Purpose

Function to get the image IDS for plate thumbnails.

188 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

```
hHandle - database connection handleIPlateID - plate ID
```

Output

arrIDs - array of thumbnail image IDs found for the plate

MDCS_PLATE_ChangeStatus

```
BOOL MDCS_PLATE_ChangeStatus (
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_E_MarkStatus eStatus)
);
```

Purpose

Function to change the plate status.

Parameters

hHandle - database connection handleIPlateID - plate IDeStatus - plate status

190 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Acquisitions Functions 10

This chapter contains functions that you can use to work with acquisitions.

Table 10-1: Acquisitions Functions

Function Name	Description
MDCS_ACQUISITION_Delete	To delete an acquisition
MDCS_ACQUISITION_Create	Create a new acquisition
MDCS_ACQUISITION_CreateProfile	To create an acquisition profile
MDCS_ACQUISITION_GetProfileRecords	To get all records from the acquisition profile table
MDCS_ACQUISITION_GetProfile	To get a record from the acquisition profile table
MDCS_ACQUISITION_GetInstanceRecord	To get an acquisition instance record
MDCS_ACQUISITION_GetBatchRecords	To get all batch records
MDCS_ACQUISITION_GetBatchRecord	To get a batch record

MDCS_ACQUISITION_Delete

```
BOOL MDCS_ACQUISITION_Delete(
HDBHANDLE hHandle,
LONGLONG lAcqID
);
```

Purpose

Function to delete an acquisition.

Parameters

```
hHandle - database connection handle lAcqID - acquisition ID
```

Return

FALSE - if error occurred

MDCS_ACQUISITION_Create

```
BOOL MDCS_ACQUISITION_Create(
HDBHANDLE hHandle,
const MDCS_ST_Acquisition& stAcq,
LONGLONG* plAcqID
);
```

Purpose

To create a new acquisition.

Parameters

hHandle - database connection handlestAcq - acquisition description

Output

plAcqID - acquisition ID for newly created acquisition

Return

FALSE - if error occurred

MDCS_ACQUISITION_CreateProfile

```
BOOL MDCS_ACQUISITION_CreateProfile(
HDBHANDLE hHandle,
const MDCS_ST_AcquisitionProfile& stAcqProf,
LONGLONG* plAcqProfID
);
```

Purpose

Function to create an acquisition profile.

Parameters

```
hHandle - database connection handlestAcqProf - profile description
```

Output

plAcqProfID - Acquisition Profile ID for newly created profile

Return

FALSE - if error occurred

MDCS_ACQUISITION_GetProfileRecords

```
BOOL MDCS_ACQUISITION_GetProfileRecords(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pReslutCallback
);
```

Purpose

To get all records from the acquisition profile table.

Parameters

hHandle - database connection handlepResultCallback - pointer to a callback function to get data

Output

Returns records with columns: ID as ACQUISITION ID NAME as acquisition name DESCRIPTION as acquisition description

Return

FALSE - if error occurred

MDCS_ACQUISITION_GetProfile

```
BOOL MDCS_ACQUISITION_GetProfile(
HDBHANDLE hHandle,
LONGLONG lProfileID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get records for an acquisition profile.

hHandle - database connection handleIProfileID - profile IDpResultCallback - pointer to a callback function to get data

Output

Returns records with columns: ID as ACQUISITION ID NAME as acquisition name DESCRIPTION as acquisition description

Return

FALSE - if error occurred

MDCS_ACQUISITION_GetInstanceRecord

```
BOOL MDCS_ACQUISITION_GetInstanceRecord(
HDBHANDLE hHandle,
LONGLONG lInstanceID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get an acquisition instance record.

Parameters

hHandle - database connection handlelInstanceID - Instance IDpResultCallback - pointer to a callback function to get data

Output

Returns records with columns: ID as ACQUISITION ID NAME as acquisition name DESCRIPTION as acquisition description START_DATE - start date END_DATE - end date

Return

FALSE - if error occurred

MDCS_ACQUISITION_GetBatchRecords

```
BOOL MDCS_ACQUISITION_GetBatchRecords(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get all batch records.

Parameters

hHandle - database connection handlepResultCallback - pointer to a callback function to get data

Output

Returns records with columns:

ID as ACQUISITION ID

ACQ_INSTANCE_ID as ACQUISITION ID

NAME as acquisition name

START_DATE - start date

END_DATE - end date

OPERATOR - name of the user who created batch

DESCRIPTION - batch description

Return

FALSE - if error occurred

MDCS_ACQUISITION_GetBatchRecord

```
BOOL MDCS_ACQUISITION_GetBatchRecord(
HDBHANDLE hHandle,
LONGLONG lBatchID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get a batch record.

hHandle - database connection handlelBatchId - batch ID

Output

Returns records with columns:

ID as ACQUISITION ID

ACQ_INSTANCE_ID as ACQUISITION ID

NAME as acquisition name

START_DATE - start date

END_DATE - end date

OPERATOR - name of user who created batch

DESCRIPTION - batch description

Return

FALSE - if error occurred

Security Functions 11

This chapter contains functions that you can use to work with security.

Table 11-1: Security Functions

Function Name	Description
MDCS_SECURITY_ManageUsersDlg	To call a dialog to manage users
MDCS_SECURITY_ManageGroupsDlg	To call a dialog to manage groups
MDCS_SECURITY_GetUserGroups	To get a list of groups that the user has access to
MDCS_SECURITY_GetUserGroupsDlg	Display a selection dialog to choose groups
MDCS_SECURITY_GetUserInfo	To get user information
MDCS_SECURITY_IsUserInAdminGroup	To check if the user is in the Administrators group
MDCS_SECURITY_GetEveryoneGroupIn fo	To get information about the everyone group
MDCS_SECURITY_ChangePassword	To change the password of the current user

MDCS_SECURITY_ManageUsersDlg

```
BOOL MDCS_SECURITY_ManageUsersDlg(
HDBHANDLE hHandle,

HWND hWnd = NULL,

LPCSTR pszDlgTitle = NULL,

MDCS_ST_DlgUsersType eDlgType = MDCS_eDlgUsersSimple,

HICON hIcon = 0

);
```

Purpose

To call a dialog to manage users.

Parameters

```
    hHandle - database connection handle
    pszDlgTitle - dialog title
    hWnd - handle to the application window
    eDlgType - type of the dialog, if MDCS_eDlgUsersWithStorage will show storage assigned to each user
    hIcon - dialog's icon
```

Return

FALSE - if error occurred or the dialog was cancelled

MDCS_SECURITY_ManageGroupsDlg

```
BOOL MDCS_SECURITY_ManageGroupsDlg(
HDBHANDLE hHandle,
HWND hWnd,

LPCSTR pszDlgTitle,
BOOL bHideUserGroupsControl = TRUE,
BOOL* pbHideUserGroupsDefault = NULL,
HICON hIcon = 0
);
```

Purpose

Calls a dialog to manage groups.

Parameters

```
hHandle - database connection handle
pszDlgTitle - dialog title
hWnd - handle to application window
pbHideUserGroupsDefault - if true, will not show user-defined groups
bHideUserGroupsControl - if true, will hide check box for user-defined
groups
hIcon - dialog's icon
```

Return

FALSE - if error occurred

MDCS_SECURITY_GetUserGroups

```
BOOL MDCS_SECURITY_GetUserGroups(
HDBHANDLE hHandle,
LONGLONG lUserID,
AxStringArray& arrGroupIDs,
AxStringArray& arrGroupNames
);
```

Purpose

Returns a list of groups that user has access to.

Parameters

hHandle - database connection handlepResultCallback - pointer to a callback function to get data

Output

```
arrGroupIDs - group IDs
arrGroupNames - group names
```

Return

FALSE - if error occurred

MDCS_SECURITY_GetUserGroupsDlg

```
BOOL MDCS_SECURITY_GetUserGroupsDlg(
HDBHANDLE hHandle,
HWND hWnd,
MDCS_GroupInfoCallback* pCallback
);
```

Purpose

Displays a selection dialog to choose groups.

Parameters

hHandle - database connection handle pResultCallback - pointer to a callback function to group information hWnd - handle to the application window

Return

FALSE - if error occurred

MDCS_SECURITY_GetUserInfo

```
BOOL MDCS_SECURITY_GetUserInfo(
HDBHANDLE hHandle,
LONGLONG lUserID,
MDCS_ST_UserInfo& stUserInfo);
```

Purpose

Returns information about a user.

Parameters

hHandle - database connection handle *IUserID* - user ID, if 0, the function will return all info for current user

Output

stUserInfo - user information

Return

FALSE - if error occurred

MDCS_SECURITY_IsUserInAdminGroup

```
BOOL MDCS_SECURITY_IsUserInAdminGroup(
HDBHANDLE hHandle,
LONGLONG lUserID,
BOOL& bAdmin
);
```

Purpose

To check if the user is in the Administrators group.

Parameters

```
hHandle - database connection handle lUserID - user ID
```

Output

bAdmin - if TRUE, the user is in Administrators group

Return

FALSE - if error occurred

MDCS_SECURITY_GetEveryoneGroupInfo

```
BOOL MDCS_SECURITY_GetEveryoneGroupInfo (
HDBHANDLE hHandle,
MDCS_ST_GroupInfo& stGroupInfo);
);
```

Purpose

To get information about the Everyone group.

Parameters

hHandle - database connection handle

Output

stGroupInfo - group information

Return

FALSE - if error occurred

MDCS_SECURITY_ChangePassword

```
BOOL MDCS_SECURITY_ChangePassword (
HDBHANDLE hHandle,
);
```

Purpose

To change the password, calls a dialog to allow the password change.

Parameters

hHandle - database connection handle

Return

FALSE - if error occurred

202 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Datasets Functions 12

This chapter contains functions that you can use to work with datasets.

Table 12-1: Dataset Functions

Function Name	Description
MDCS_DATASET_InsertAnalysisAttributes	To create/insert new dataset analysis attributes
MDCS_DATASET_GetAnalysisAttributesByData set	To get all analysis attributes for a dataset
MDCS_DATASET_GetAnalysisAttributes	To get dataset analysis attributes per analysis
MDCS_DATASET_GetAnalysisDescriptions	To get all analysis descriptions per dataset
MDCS_DATASET_FindAnalysisAttributesRecord	To find a record by attributes
MDCS_DATASET_DeleteAnalysisAttributes	To delete analysis attributes
MDCS_DATASET_UpdateAnalysisAttributes	To update analysis attributes
MDCS_DATASET_ManageFolderSecurity	To call a dialog to manage security access to the dataset folders
MDCS_DATASET_DeleteFolder	To delete the dataset folder
MDCS_DATASET_Copy	To copy a dataset
MDCS_DATASET_Update	To update a dataset
MDCS_DATASET_Get	To get dataset information
MDCS_DATASET_Create	To create dataset
MDCS_DATASET_AddAssay	To add an assay to a dataset
MDCS_DATASET_HavePermissionsToModify	To check if the current user can modify a dataset
MDCS_DATASET_GetAnalysisDescription	To get dataset analysis configuration description.
MDCS_DATASET_GetAssayAndFilters	To get assays and filters that were used to create a dataset
MDCS_DATASET_GetAllInFolder	To get assays (with a basic description) in a folder
MDCS_DATASET_GetItemFolder	To get an item's folder
MDCS_DATASET_GetAssaySiblingFolders	To get assay sibling folders
MDCS_DATASET_DoesSubFolderExist	To check if a subfolder exists in the database
MDCS_DATASET_CreateFolder	To create a new dataset folder in the database
MDCS_DATASET_ModifyFolder	To modify a dataset folder in the database
MDCS_DATASET_GetSiblingFolders	To get dataset sibling folders
MDCS_DATASET_GetAllMeasurements	To get data types for a dataset

Table 12-1: Dataset Functions (cont'd)

MDCS_DATASET_AddScriptletAssay	To create and add a scriptlet to dataset
MDCS_DATASET_GetScriptletAssays	To get all scriptlets assays in dataset
MDCS_DATASET_GetScriptletAssayIDByName	To find a scriptlet by name in the dataset
MDCS_DATASET_GetAllForPlate	To get all available datasets for a plate
MDCS_DATASET_GetAllForAssay	To get all available datasets for an assay
MDCS_DATASET_GetAllMSetParamValues	To get the values of attributes in a dataset
MDCS_DATASET_GetResultsInfoByConfig	To get the dataset result info for a specified configuration
MDCS_DATASET_GetResultInfo	To get dataset result info
MDCS_DATASET_CallbackToAnalysisInfo	To convert the result of callback to MDCS_ST_DatasetResultInfo structure
MDCS_DATASET_DoesNameExist	To check if dataset with specified name already exists in the folder
MDCS_DATASET_UpdateFolderItem	To update the link to the folder of a dataset
MDCS_DATASET_Delete	To delete a dataset
MDCS_DATASET_GetAnalysisCount	To get a number of analysis per dataset

MDCS_DATASET_InsertAnalysisAttributes

```
BOOL MDCS_DATASET_InsertAnalysisAttributes(
HDBHANDLE hHandle,

const MDCS_ST_DsetAnalysis* pDsetAnalysis,

const MDCS_ST_DsetAnalysisAttr* pAttributes,

UINT uNumberAttr,

LONGLONG* plAnalysisID

);
```

Purpose

To create/insert new dataset analysis attributes.

Parameters

hHandle - database connection handlepDsetAnalysis - description of analysispAttributes - pointer to the attributes arrayuNumberAttr - number of attributes in array

Output

plAnalysisID - new analysis ID

Return

FALSE - if error occurred

MDCS_DATASET_GetAnalysisAttributesByDataset

```
BOOL MDCS_DATASET_GetAnalysisAttributesByDataset(
HDBHANDLE hHandle,
LONGLONG lDatasetID,
MDCS_GetDBResultsCCallback * pCallBack
);
```

Purpose

To get all analysis attributes for a dataset.

Parameters

hHandle - database connection handlelDatasetID - dataset IDpResultCallback - pointer to a callback function to get data

Output

Callback returns the following fields:

ATTR_ID - unique ID

CALC_METHOD - calculation method

ATTR_TYPE - attribute type (pivot, data type)

ANALYSIS_ID - analysis ID (pointer to a description of analysis)

ORDER_NUM - order number

COLLAPSE_COLUMN - column that is used to collapse data

ATTR_COLUMN - name of the attribute column

Return

FALSE - if error occurred

MDCS_DATASET_GetAnalysisAttributes

```
BOOL MDCS_DATASET_GetAnalysisAttributes(
HDBHANDLE hHandle,
LONGLONG lAnalysisID,
MDCS_GetDBResultsCCallback * pCallBack
);
```

Purpose

To get dataset analysis attributes per analysis.

Parameters

hHandle - database connection handlelAnalsysisID - analysis IDpResultCallback - pointer to a callback function to get data

Output

```
Callback returns the following fields:

ATTR_ID - unique ID

CALC_METHOD - calculation method

ATTR_TYPE - attribute type (pivot, data type)

ANALYSIS_ID - analysis ID (pointer to a description of analysis)

ORDER_NUM - order number

COLLAPSE_COLUMN - column that is used to collapse data

ATTR_COLUMN - name of the attribute column
```

Return

FALSE - if error occurred

MDCS_DATASET_GetAnalysisDescriptions

```
BOOL MDCS_DATASET_GetAnalysisDescriptions(
HDBHANDLE hHandle,

LONGLONG lDatasetID,

MDCS_GetDBResultsCCallback * pCallBack
);
```

Purpose

To get all analysis descriptions per dataset.

Parameters

hHandle - database connection handlelDatasetID - dataset IDpResultCallback - pointer to a callback function to get data

Output

Callback returns the following fields:
ANALYSIS_ID - unique ID
DATASET_ID - ID of a dataset where it belongs
NAME - analysis description name
DESCRIPTION - description

Return

FALSE - if error occurred

MDCS_DATASET_FindAnalysisAttributesRecord

```
BOOL MDCS_DATASET_FindAnalysisAttributesRecord(
HDBHANDLE hHandle,

LONGLONG lDatasetID,

const MDCS_ST_DsetAnalysisAttr* pAttributes,

UINT uNumberAttr

MDCS_ST_DsetAnalysis* pDsetAttributes
):
```

Purpose

To find a record by attributes.

Parameters

```
hHandle - database connection handlelDatasetID - dataset IDpAttributes - pointer to the attributes arrayuNumberAttr - number of attributes in array
```

Output

pDsetAttributes - record that corresponds to attributes

Return

FALSE - if error occurred

MDCS_DATASET_DeleteAnalysisAttributes

```
BOOL MDCS_DATASET_DeleteAnalysisAttributes(
HDBHANDLE hHandle,
LONGLONG lAttributesID
);
```

Purpose

To delete analysis attributes.

Parameters

hHandle - database connection handle *|AttributesID* - attributes ID

Return

FALSE - if error occurred

MDCS_DATASET_UpdateAnalysisAttributes

```
BOOL MDCS_DATASET_UpdateAnalysisAttributes(
HDBHANDLE hHandle,
const MDCS_ST_DsetAnalysis* pDsetAnalysis
);
```

Purpose

To update analysis attributes.

Parameters

hHandle - database connection handlepDsetAnalysis - structure that contains attributes to update

Return

FALSE - if error occurred

MDCS_DATASET_ManageFolderSecurity

```
BOOL MDCS_DATASET_ManageFolderSecurity(
HDBHANDLE hHandle,
LONGLONG lFolderID,
HWND hWnd = NULL,
LPCSTR pszDlgTitle = NULL
);
```

Purpose

Calls a dialog to manage security access to the dataset folders.

Parameters

```
hHandle - database connection handlelFolderID - folder IDpszDlgTitle - dialog titlehWnd - handle to the application window
```

Return

FALSE - if error occurred or dialog cancelled

MDCS_DATASET_DeleteFolder

```
BOOL MDCS_DATASET_DeleteFolder(
HDBHANDLE hHandle,
LONGLONG lFolderID
);
```

Purpose

Function to delete dataset folder.

Parameters

hHandle - database connection handlelFolderID - ID of a folder in datasets tree

Return

FALSE - if error occurred

MDCS_DATASET_Copy

```
BOOL MDCS_DATASET_Copy(
HDBHANDLE hHandle,

LONGLONG lDatasetSourceID,

LONGLONG lDestinationFolderID,

MDCS_ST_Dataset& stDataSetOut,

WORD wCopyDependants,

BOOL bCopySpots = TRUE,

BOOL bCopyAssays = TRUE

);
```

Purpose

Function to copy a dataset.

hHandle - database connection handle

IDatasetSourceID - ID of the dataset that should be copied

 ${\it IDestinationFolderID}$ - ID of the folder where the dataset should be copied

wCopyDependants - option to copy dataset dependants such as results, objects, scriptlets (use MDC_E_DatasetObjectType)

bCopySpots - if TRUE, will copy all spots that were used in dataset

bCopyAssays - if TRUE, will copy all measurement sets associated with dataset

Output

stDataSetOut - dataset that was created

Return

FALSE - if error occurred

MDCS_DATASET_Update

```
BOOL MDCS_DATASET_Update(
HDBHANDLE hHandle,
const MDCS_ST_Dataset& stDataSetIn
);
```

Purpose

Function to update a dataset.

Parameters

hHandle - database connection handlestDataSetIn - dataset information that will be updated, should contain valid dataset ID

Return

FALSE - if error occurred

MDCS_DATASET_Get

```
BOOL MDCS_DATASET_Get(
HDBHANDLE hHandle,
LONGLONG lDatasetID,
MDCS_ST_Dataset& stDataSetOut
);
```

Purpose

Function to get dataset information.

212 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

hHandle - database connection handle | *DatasetID* - dataset ID

Output

stDataSetOut - dataset information

Return

FALSE - if error occurred

MDCS_DATASET_Create

```
BOOL MDCS_DATASET_Create(
HDBHANDLE hHandle,

LONGLONG lDestFolderID,

const MDCS_ST_Dataset& stDataSetIn,

MDCS_ST_Dataset& stDataSetOut
);
```

Purpose

Function to create a new dataset.

Parameters

hHandle - database connection handleIDestFolderID - folder ID where dataset should be created stDataSetIn - dataset description

Output

stDataSetOut - dataset information

Return

FALSE - if error occurred

MDCS_DATASET_AddAssay

```
BOOL MDCS_DATASET_AddAssay(
HDBHANDLE hHandle,
LONGLONG lDatasetID,
LONGLONG lAssayID
);
```

Purpose

Function to create a new dataset assay.

```
hHandle - database connection handleIDatasetID - destination dataset IDIAssayID - assay ID
```

Output

stDataSetOut - dataset information

Return

FALSE - if error occurred

MDCS_DATASET_HavePermissionsToModify

```
BOOL MDCS_DATASET_HavePermissionsToModify(
HDBHANDLE hHandle,
LONGLONG lDatasetID,
BOOL& bCanModify
);
```

Purpose

To check if the current user can modify a dataset.

Parameters

```
hHandle - database connection handle IDatasetID - dataset ID
```

Output

bCanModify - if TRUE, the user can modify the dataset

Return

FALSE - if error occurred

MDCS_DATASET_GetAnalysisDescription

```
BOOL MDCS_DATASET_GetAnalysisDescription (
HDBHANDLE hHandle,

LONGLONG lAnalysisID,

MDCS_ST_DsetAnalysis& stDsetAttributes
);
```

Purpose

To get a dataset analysis configuration description.

```
hHandle - database connection handlelAnalyisisID - configuration ID
```

Output

stDsetAttributes - structure that describes the analysis configuration

Return

FALSE - if error occurred

MDCS_DATASET_GetAssayAndFilters

```
BOOL MDCS_DATASET_GetAssayAndFilters (
HDBHANDLE hHandle,
LONGLONG lDatasetID,
AxStringArray& arrAssayIDs,
AxStringArray& arrAssayFilters
);
```

Purpose

To get assays and filters that were used to create a dataset.

Parameters

```
hHandle - database connection handle IDatasetID - dataset ID
```

Output

```
arrAssayIDs - array of assay IDs
arrAssayFilters - array of assay filters
```

Return

FALSE - if error occurred

MDCS_DATASET_GetAllInFolder

```
BOOL MDCS_DATASET_GetAllInFolder (
HDBHANDLE hHandle,
LONGLONG lFolderID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get assays (with a basic description) in a folder.

```
hHandle - database connection handlelFolderID - folder IDpResultCallback - pointer to a callback function to get data
```

Output

```
Fields that are always present:

DATASET_ID, DATASET_NAME, TIME_CREATED (as time_t value)

DATASET_DESCRIPTION - dataset description

DATASET_TYPE - type of the dataset

DATASET_DTYPE

ENGLISH_DESCR - dataset description

TO DELETE
```

Return

FALSE - if error occurred

MDCS_DATASET_GetItemFolder

```
BOOL MDCS_DATASET_GetItemFolder (
HDBHANDLE hHandle,
LONGLONG litemID,
LONGLONG& lFolderID
);
```

Purpose

To get an item's folder id.

Parameters

```
hHandle - database connection handlelItemID - item ID
```

Output

IFolderID - folder ID

Return

FALSE - if error occurred

MDCS_DATASET_GetAssaySiblingFolders

```
BOOL MDCS_DATASET_GetAssaySiblingFolders (
HDBHANDLE hHandle,
LONGLONG lParentFolderID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get assay sibling folders.

Parameters

hHandle - database connection handleParentFolderID - parent ID to retrieve siblings forpResultCallback - pointer to a callback function to get data

Output

```
Result columns will be:
FOLDER_NAME - folder name
FOLDER_ID - folder ID
```

Return

FALSE - if error occurred

MDCS_DATASET_DoesSubFolderExist

```
BOOL MDCS_DATASET_DoesSubFolderExist (
HDBHANDLE hHandle,
LONGLONG lFolderID,
LPCSTR pzName,
LONGLONG &lSubfolderID
);
```

Purpose

To check if the subfolder exists in the database.

Parameters

```
hHandle - database connection handleIFolderID - parent folder ID to retrieve subfolder for pzName - name of the parent folder
```

Output

ISubfolderID - subfolder ID, if 0 subfolder does not exists

Return

FALSE - if error occurred

MDCS_DATASET_CreateFolder

```
BOOL MDCS_DATASET_CreateFolder (
HDBHANDLE hHandle,
const MDCS_ST_FolderInfo& stInfoIn,
MDCS_ST_FolderInfo & stInfoOut
);
```

Purpose

To create a new dataset folder in the database.

Parameters

hHandle - database connection handlestInfoIn - structure describes the folder to be created in the database

Output

stInfoOut - structure that describes the new folder just created

Return

FALSE - if error occurred

MDCS_DATASET_ModifyFolder

```
BOOL MDCS_DATASET_ModifyFolder (
HDBHANDLE hHandle,
const MDCS_ST_FolderInfo& stInfoIn
);
```

Purpose

To modify a dataset folder in the database.

Parameters

hHandle - database connection handlestInfoIn - structure that describes the folder to be modified

Output

NONE

Return

FALSE - if error occurred

MDCS_DATASET_GetSiblingFolders

```
BOOL MDCS_DATASET_GetSiblingFolders (
HDBHANDLE hHandle,
LONGLONG lParentFolderID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get dataset sibling folders.

Parameters

hHandle - database connection handleIParentFolderID - parent ID to retrieve siblings for

Output

```
pResultCallback - data callback
Result columns will be
FOLDER_NAME - folder name
FOLDER_ID - folder ID
```

Return

FALSE - if error occurred

MDCS_DATASET_GetAllMeasurements

```
BOOL MDCS_DATASET_GetAllMeasurements (
HDBHANDLE hHandle,
LONGLONG lDatasetID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get data types for a dataset.

Parameters

hHandle - database connection handlelDatasetID - dataset IDpResultCallback - pointer to a callback function to get data

Output

Result columns will be: TABLE_ID - table name COLUMN_NAME COLUMN_TYPE COLUMN_NAME_EXT FUNCTION_NAME PARAMETER_NAME

Return

FALSE - if error occurred

MDCS_DATASET_AddScriptletAssay

```
BOOL MDCS_DATASET_AddScriptletAssay (
HDBHANDLE hHandle,
LONGLONG lDatasetID,
LONGLONG lScriptletID,
LPCSTR pszAssayName,
LONGLONG& lAssayIDOut
);
```

Purpose

To create and add a scriptlet to a dataset.

Parameters

hHandle - database connection handle IDatasetID - dataset ID IScriptletID - scriptlet assay ID pszAssayName - name of assay

Output

pAssayIDOut - created scriptlet assay ID

Return

FALSE - if error occurred

MDCS_DATASET_GetScriptletAssays

```
BOOL MDCS_DATASET_GetScriptletAssays (
HDBHANDLE hHandle,
LONGLONG lDatasetID,
LONGLONG lScriptletID,
AxStringArray& arrAssayIDs,
AxStringArray& arrAssayNames
);
```

Purpose

To get a scriptlet assay in a dataset.

Parameters

hHandle - database connection handleIDatasetID - dataset IDIScriptletID - assay scriptlet ID

Output

arrAssayIDs - array of scriptlet assay IDs
arrAssayNames - array of scriptlet assay names

Return

FALSE - if error occurred

MDCS_DATASET_GetScriptletAssayIDByName

```
BOOL MDCS_DATASET_GetScriptletAssayIDByName (
HDBHANDLE hHandle,
LONGLONG lDatasetID,
LPCSTR pszName,
LONGLONG& lAssayID
);
```

Purpose

To find a scriptlet assay ID by name.

Parameters

```
hHandle - database connection handleIDatasetID - dataset IDpszName - assay name
```

Output

IAssayID - assay ID

Return

FALSE - if error occurred

MDCS_DATASET_GetAllForPlate

```
BOOL MDCS_DATASET_GetAllForPlate (
HDBHANDLE hHandle,
LONGLONG lPlateID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get available datasets for a plate.

Parameters

```
hHandle - database connection handleIPlateID - plate IDpResultCallback - pointer to a callback function to get data
```

Output

Records include fields from the DATASET table.

Return

FALSE - if error occurred

MDCS_DATASET_GetAllForAssay

```
BOOL MDCS_DATASET_GetAllForAssay (
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_GetDBResultsCCallback * pResultCallback
);
```

Purpose

To get available datasets for an assay.

Parameters

```
hHandle - database connection handlelAssayID - assay IDpResultCallback - pointer to a callback function to get data
```

Output

Records include fields from the DATASET table.

Return

FALSE - if error occurred

MDCS_DATASET_GetAllMSetParamValues

```
BOOL MDCS_DATASET_GetAllMSetParamValues (
HDBHANDLE hHandle,
LONGLONG lDatasetID,
const MDCS_ST_ScopeAttribute* arrAttributes,
INT_PTR nElemCount,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get values of attributes in a dataset.

Parameters

```
hHandle - database connection handle

lDatasetID - dataset ID

arrAttributes - list of attributes to retrieve values

nElemCount - number of attributes

pResultCallback - pointer to a callback function to get data
```

Output

Records include measurement set attribute values.

Return

FALSE - if error occurred

MDCS_DATASET_GetResultsInfoByConfig

```
BOOL MDCS_DATASET_GetResultsInfoByConfig (
HDBHANDLE hHandle,
LONGLONG lDatasetID,
LONGLONG lConfigID,
MDCS_GetDBResultsCCallback * pCallBack
);
```

Purpose

To get description of results on dataset for specified configuration.

Parameters

```
    hHandle - database connection handle
    lDatasetID - dataset ID
    lConfigID - configuration ID
    pResultCallback - pointer to a callback function to get data
```

Output

Records include analysis info.

Return

FALSE - if error occurred



Note: Use MDCS_DATASET_CallbackToAnalysisInfo to process results.

MDCS_DATASET_GetResultInfo

```
BOOL MDCS_DATASET_GetResultInfo (
HDBHANDLE hHandle,
LONGLONG lResultID,
MDCS_GetDBResultsCCallback * pCallBack
);
```

Purpose

To get dataset.results info.

Parameters

```
hHandle - database connection handleIResultID - result IDpResultCallback - pointer to a callback function to get data
```

Output

Records include analysis info

Return

FALSE - if error occurred

MDCS_DATASET_CallbackToAnalysisInfo

```
BOOL MDCS_DATASET_CallbackToAnalysisInfo (
const MDCS_QueryResults* pQueryResults,

MDCS_ST_DatasetResultInfo& stOut,

BOOL bResultsInOrigForm = FALSE
);
```

Purpose

To convert the result of a callback to the MDCS_ST_DatasetResultInfo structure.

Parameters

```
pQueryResults - call back results
bResultsInOrigForm - if TRUE, the results are in the original form
```

Output

stOut - destination structure

Return

FALSE - if error occurred

MDCS_DATASET_DoesNameExist

```
BOOL MDCS_DATASET_DoesNameExist (
HDBHANDLE hHandle,

LPCSTR pszName,

LONGLONG lFolderID,

BOOL& bExists
);
```

Purpose

To check if the dataset exists.

Parameters

```
hHandle - database connection handlepszName - dataset namelFolderID - folder ID
```

Output

bExists - TRUE if exists

Return

FALSE - if error occurred

MDCS_DATASET_UpdateFolderItem

```
BOOL MDCS_DATASET_UpdateFolderItem(
HDBHANDLE hHandle,
LONGLONG lFolderID,
LONGLONG lDatasetID,
);
```

Purpose

Function update the link to the folder of a dataset.

Parameters

```
hHandle - database connection handlelDatasetID - dataset IDlFolderID - folder ID
```

Return

FALSE - if error occurred

MDCS_DATASET_Delete

```
BOOL MDCS_DATASET_Delete(
HDBHANDLE hHandle,
LONGLONG lDatasetID,
);
```

Purpose

Function to delete a dataset.

Parameters

hHandle - database connection handle IDatasetID - dataset ID

Return

FALSE - if error occurred

MDCS_DATASET_GetAnalysisCount

```
BOOL MDCS_DATASET_GetAnalysisCount (
HDBHANDLE hHandle,
LONGLONG lDatasetID,
LONGLONG lCount,
);
```

Purpose

To get the number of analyses per dataset.

Parameters

hHandle - database connection handle *IDatasetID* - dataset ID

Output

ICount - number of analyses

Return

FALSE - if error occurred

This chapter contains functions that you can use to work with assay images and assay normalization.

Table 13-1: Assay Images and Assay Normalization Functions

Function Name	Description
MDCS_ASSAYIMAGES_GetAllPropertyAttributes	To get the attributes that the measurement set images can be queried on
MDCS_ASSAYIMAGES_GetAllByAttributes	Function to get measurement set images based on specified attributes values
MDCS_ASSAYIMAGES_GetUniqueAttributeValues	To get unique values for a measurement set images attribute
MDCS_ASSAYNORM_CreateConfig	To create a new normalization configuration
MDCS_ASSAYNORM_UpdateConfig	To update an existing normalization configuration
MDCS_ASSAYNORM_GetConfig	To get a normalization configuration using its ID
MDCS_ASSAYNORM_GetConfigByName	To get a normalization configuration using its name
MDCS_ASSAYNORM_GetConfigAll	To get all normalization configurations
MDCS_ASSAYNORM_GetConfigForAssay	To get all normalization configurations that exist for an assay

MDCS_ASSAYIMAGES_GetAllPropertyAttributes

```
BOOL MDCS_ASSAYIMAGES_GetAllPropertyAttributes(
HDBHANDLE hHandle,
AxStringArray &arrayColumnIDs,
AxStringArray &arrayColumnNames
);
```

Purpose

To get the attributes that measurement set images can be queried on.

Parameters

hHandle - database connection handle arrayColumnIDs - array to hold column ID arrayColumnNames - array to hold column names

Return

FALSE - if error occurred

MDCS_ASSAYIMAGES_GetAllByAttributes

```
BOOL MDCS_ASSAYIMAGES_GetAllByAttributes(
HDBHANDLE hHandle,
const AxStringArray& arrayAssayIDs,
const AxStringArray* pAXDisplayColumns,
const AxStringArray* pAXAttrColumn,
const AxStringArray* pAXAttrValues,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get measurement set images based on specified attributes values.

Parameters

```
hHandle - database connection handle
arrayAssayIDS - assay IDs
pAXDisplayColumns - columns to retrieve
pAXAttrColumn - columns to query on
pAXAttrValues - values of the column to query on
pResultCallback - pointer to a callback function to get data
```

Output

```
Result columns are:

OBJ_ID - object ID

PLATE_ID - plate ID

IMAGE_ID - plate ID

SITE_ID - site ID

SERIES_ID - series ID

WELL_X - well X

WELL_Y - well Y

X_POSITION - x position

Y_POSITION - y position

Z_INDEX - z-index

Z_POSITION - z- position

T_INDEX

T_POSITION - T position in seconds

T MICROSECONDS
```

Return

FALSE - if error occurred

MDCS_ASSAYIMAGES_GetUniqueAttributeValues

```
BOOL MDCS_ASSAYIMAGES_GetUniqueAttributeValues(
HDBHANDLE hHandle,
const AxStringArray& arrayAssayIDs,
LPCSTR pszAttrColumnName,
const AxStringArray* pAXAttrColumnNames,
const AxStringArray* pAXAttrValues,
AxStringArray* pAXValuesOut
);
```

Purpose

Function to get unique values for a measurement set images attribute.

Parameters

hHandle - database connection handle
arrayAssayIDs - array of assay IDs
pszAttrColumnName - column that will contain unique values
pAXAttrColumnNames - attribute column names to query on
pAXAttrValues - attribute values

Output

pAXValuesOut - array of unique values

Return

FALSE - if error occurred

MDCS_ASSAYNORM_CreateConfig

```
BOOL MDCS_ASSAYNORM_CreateConfig (
HDBHANDLE hHandle,
const MDCS_ST_NormConfig& stConfigIn,
MDCS_ST_NormConfig* pstConfigOut
);
```

Purpose

Function to a create new normalization configuration.

Parameters

hHandle - database connection handle stConfigIn - configuration to create

Output

pstConfigOut - created configuration

Return

FALSE - if error occurred

MDCS_ASSAYNORM_UpdateConfig

```
BOOL MDCS_ASSAYNORM_UpdateConfig (
HDBHANDLE hHandle,
const MDCS_ST_NormConfig& stConfigIn
);
```

Purpose

Function to update a normalization configuration.

Parameters

hHandle - database connection handlestConfigIn - configuration that will be updated (ID must exist)

Output

NONE

Return

FALSE - if error occurred

MDCS_ASSAYNORM_GetConfig

```
BOOL MDCS_ASSAYNORM_GetConfig (
HDBHANDLE hHandle,
LONGLONG lConfigID,
MDCS_ST_NormConfig& stConfigOut
);
```

Purpose

Function to get a normalization configuration using its ID.

Parameters

hHandle - database connection handlelConfigID - configuration ID

Output

stConfigOut - configuration for IConfigID

Return

FALSE - if error occurred

MDCS_ASSAYNORM_GetConfigByName

```
BOOL MDCS_ASSAYNORM_GetConfigByName (
HDBHANDLE hHandle,

LPCSTR pszConfigName,

MDCS_ST_NormConfig& stConfigOut
);
```

Purpose

Function to get a normalization configuration using its name.

Parameters

hHandle - database connection handlepszConfigName - Name of configuration to retrieve

Output

stConfigOut - configuration for pszConfigName

Return

FALSE - if error occurred

MDCS_ASSAYNORM_GetConfigAll

```
BOOL MDCS_ASSAYNORM_GetConfigAll (
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get all existing normalization configurations.

Parameters

hHandle - database connection handle pResultCallback - pointer to a callback function to get data

Output

Returns columns that could be mapped to the MDCS_ST_NormConfig structure.

Return

FALSE - if error occurred

MDCS_ASSAYNORM_GetConfigForAssay

```
BOOL MDCS_ASSAYNORM_GetConfigForAssay (
HDBHANDLE hHandle,
LONGLONG lAssayID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get all existing normalization configurations for an assay.

Parameters

```
hHandle - database connection handlelAssayID - assay IDpResultCallback - pointer to a callback function to get data
```

Output

Return columns that could be mapped to the MDCS_ST_NormConfig structure.

Return

FALSE - if error occurred

236 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Quicklist Functions 14

This chapter contains functions that you can use to work with quicklists.

Table 14-1: Quicklist Functions

Function Name	Purpose
MDCS_QUICKLIST_ManageFolderSecurity	Call a dialog to manage security access to the quicklist folders

MDCS_QUICKLIST_ManageFolderSecurity

```
BOOL MDCS_QUICKLIST_ManageFolderSecurity(
HDBHANDLE hHandle,
LONGLONG lFolderID,
HWND hWnd = NULL,
LPCSTR pszDlgTitle = NULL
):
```

Purpose

Calls a dialog to manage security access to the quicklist folders.

Parameters

```
hHandle - database connection handlelFolderID - folder IDpszTitle - dialog titlehWnd - handle to the application window
```

Return

FALSE - if dialog was cancelled

238 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Application Data Table Functions 15

This chapter contains functions that you can use to work with application data tables.

Table 15-1: Application Data Table Functions

Function Name	Description
MDCS_APP_ManageLocationOptions	To call a dialog that allows the user to edit location options
MDCS_APP_CreateLocationOption	To create a new record that describes the location option
MDCS_APP_CreateJobRecord	To create a new job
MDCS_APP_CreateParameter	To create a new application parameter
MDCS_APP_GetParameterByName	To get free space in a file location
MDCS_APP_GetParameterByName	To get an application parameter using its name
MDCS_APP_GetParameterByID	To get an application parameter using its ID
MDCS_APP_UpdateParameter	To update an application parameter using its ID
MDCS_APP_GetLocationOptions	To get all records from the LOCATION_OPTIONS table
MDCS_APP_GetLocationOptionRecord	To get a single record from the LOCATION_OPTIONS table
MDCS_APP_GetLocationRecordsByLabel	To get all records from location table filtered by Label value
MDCS_APP_DeleteLocationOption	To delete a location option by its ID
MDCS_APP_DeleteUserLocationOptions	To delete a user location options by the user login
MDCS_APP_SetUserLocationOptions	To assign a user location option
MDCS_APP_GetUserLocationOption	To get a default location option
MDCS_APP_GetJobQueueRecords	To get all records from the JOB_QUEUE table
MDCS_APP_GetJobQueueRecord	To get a record from the JOB_QUEUE table.
MDCS_APP_GetJobQueue	To get job queue information using the JOB ID
MDCS_APP_UpdateJobStatus	To update job status
MDCS_APP_UpdateJobProgress	To update job progress
MDCS_APP_UpdateJobQueueRecord	To update job queue information using the job ID
MDCS_APP_CancelJobProgress	To cancel job
MDCS_APP_ClaimJob	To claim a new job
MDCS_APP_ResetJob	To reset job status and progress

Table 15-1: Application Data Table Functions (cont'd)

MDCS_APP_ResetCrashedJobs	To reset job status and progress of all crashed jobs from last run
MDCS_APP_RefreshAllJobs	To refresh the status of all jobs
MDCS_APP_GetLocationOptionByID	To get a single record from the LOCATION_OPTIONS table by ID

MDCS_APP_ManageLocationOptions

```
BOOL MDCS_APP_ManageLocationOptions(
HDBHANDLE hHandle,

MDCS_ST_LocationOption* pstLocationOption,

HWND hWnd = NULL,

const LONGLONG* plDefaultSelectionID = NULL,

LPCSTR pszDlgTitle = NULL,

HICON hIcon = 0

BOOL bSelectionMode = TRUE

);
```

Purpose

Function calls dialog that allows user editing of location options.

Parameters

```
hHandle - database connection handle
hWnd - handle to the parent window
plDefaultSelectionID - optional ID for the record that should be selected by default in dialog
pszDlgTitle - dialog title
hIcon - dialog's icon
bSelectionMode - if TRUE, dialog is called in the selection mode (otherwise just display mode)
```

Output

pstLocationOption - selected location option

Return

FALSE - if error occurred

MDCS_APP_CreateLocationOption

```
BOOL MDCS_APP_CreateLocationOption(
HDBHANDLE hHandle,
const MDCS_ST_LocationOption& stLocation,
LONGLONG* plLocationID
);
```

Purpose

Function to create a new record that describes a location option.

Parameters

```
hHandle - database connection handle
stLocation - location description
```

Output

plLocationID - location ID for the new location option

Return

FALSE - if error occurred

MDCS_APP_CreateJobRecord

```
BOOL MDCS_APP_CreateJobRecord(
HDBHANDLE hHandle,
const MDCS_ST_Job& stJob,
LONGLONG* plJobID
);
```

Purpose

Function to create a new job.

Parameters

```
hHandle - database connection handlestJob - job description
```

Output

plJobID - Job ID for new job record

Return

FALSE - if error occurred

MDCS_APP_CreateParameter

```
BOOL MDCS_APP_CreateParameter(
HDBHANDLE hHandle,
const MDCS_ST_AppParameter& stParam,
LONGLONG* plParamID
);
```

Purpose

Function to create a new application parameter.

Parameters

```
hHandle - database connection handlestParam - application parameter
```

Output

plParamID - Parameter ID for new parameter

Return

FALSE - if error occurred

MDCS_APP_GetParameterByName

```
BOOL MDCS_APP_GetFreeSpaceFromFileLocation(
HDBHANDLE hHandle,
const MDCS_ST_LocationOption& stLocation,
LONGLONG* 1FreeSpace
);
```

Purpose

Function to get free disk space in a file location.

Parameters

```
hHandle - database connection handle
stLocation - location description
```

Output

IFreeSpace - free disk space (in kilobytes)

Return

FALSE - if error occurred

MDCS_APP_GetParameterByName

```
BOOL MDCS_APP_GetParameterByName(
HDBHANDLE hHandle,

LPCSTR pszParamName,

MDCS_ST_AppParameter& stParam
);
```

Purpose

Function to get an application parameters using its name.

Parameters

hHandle - database connection handle pszParamName - application parameter name

Output

stParam - application parameter information structure

Return

FALSE - if error occurred

MDCS_APP_GetParameterByID

```
BOOL MDCS_APP_GetParameterByID(
HDBHANDLE hHandle,
LONGLONG lParamID,
MDCS_ST_AppParameter& stParam
);
```

Purpose

Function to get application parameter using its ID.

Parameters

```
hHandle - database connection handleIParamID - application parameter ID
```

Output

stParam - application parameter information structure

Return

FALSE - if error occurred

MDCS_APP_UpdateParameter

```
BOOL MDCS_APP_UpdateParameter(
HDBHANDLE hHandle,
LONGLONG lParamID,
const MDCS_ST_AppParameter& stParam
);
```

Purpose

Function to update application parameter using its ID.

Parameters

```
hHandle - database connection handleIParamID - application parameter IDstParam - application parameter information structure
```

Return

FALSE - if error occurred.

MDCS_APP_GetLocationOptions

```
BOOL MDCS_APP_GetLocationOptions(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get all records from the LOCATION_OPTIONS table.

Parameters

```
hHandle - database connection handlepResultCallback - callback object
```

Output

Columns from LOCATION_OPTIONS table

Return

FALSE - if error occurred

MDCS_APP_GetLocationOptionRecord

```
BOOL MDCS_APP_GetLocationOptionRecord(
HDBHANDLE hHandle,
LONGLONG lLocationID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get a single record from the LOCATION OPTIONS table.

Parameters

hHandle - database connection handlepResultCallback - pointer to a callback function to get datalLocationID - location option ID

Output

Columns from LOCATION OPTIONS table

Return

FALSE - if error occurred

${\bf MDCS_APP_GetLocationRecordsByLabel}$

```
BOOL MDCS_APP_GetLocationRecordsByLabel(
HDBHANDLE hHandle,

LPCSTR pszLabel,

MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get all records from the LOCATION_OPTIONS table filtered by Label value.

Parameters

```
hHandle - database connection handlepszLabel - labelpResultCallback - pointer to a callback function to get data
```

Output

Columns from LOCATION OPTIONS table

Return

FALSE - if error occurred

MDCS_APP_DeleteLocationOption

```
BOOL MDCS_APP_DeleteLocationOption(
HDBHANDLE hHandle,
LONGLONG lOptionID
);
```

Purpose

Function to delete a location option by its ID.

Parameters

hHandle - database connection handlelOptionID - location option ID

Return

FALSE - if error occurred

MDCS_APP_DeleteUserLocationOptions

```
BOOL MDCS_APP_DeleteUserLocationOptions(
HDBHANDLE hHandle,
LONGLONG lUserID
);
```

Purpose

Function to delete user location options by user login.

Parameters

hHandle - database connection handle *lUserID* - user ID

Return

FALSE - if error occurred

MDCS_APP_SetUserLocationOptions

```
BOOL MDCS_APP_SetUserLocationOptions(
HDBHANDLE hHandle,

LONGLONG lOptionID,

LONGLONG lUserID,

MDCS_E_UserDataLocationOptions eType =
MDCS_eUserDataOptionDefault
);
```

Purpose

Function to assign user location options.

Parameters

```
    hHandle - database connection handle
    lOptionID - option ID
    lUserID - user ID
    eType - type of the option (MDCS_eUserDataOptionDefault - by default)
```

Return

FALSE - if error occurred

MDCS_APP_GetUserLocationOption

```
BOOL MDCS_APP_GetUserLocationOption(
HDBHANDLE hHandle,

LONGLONG lUserID,

MDCS_ST_LocationOption &stLocationOption,

MDCS_E_UserDataLocationOptions eType =

MDCS_eUserDataOptionDefault
);
```

Purpose

Function to get the default location option.

Parameters

```
    hHandle - database connection handle
    lOptionID - option ID
    lUserID - user ID
    eType - type of the option (MDCS_eUserDataOptionDefault - by default).
```

Output

stLocationOption - location option

Return

FALSE - if error occurred

MDCS_APP_GetJobQueueRecords

```
BOOL MDCS_APP_GetJobQueueRecords(
HDBHANDLE hHandle,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get all records from the JOB_QUEUE table.

Parameters

hHandle - database handlepResultCallback - pointer to a callback function to get data

Output

Returns columns:

ID

TIMESTAMP - as time updated All other columns from JOB_QUEUE table.

Return

FALSE - if error occurred

MDCS_APP_GetJobQueueRecord

```
BOOL MDCS_APP_GetJobQueueRecord(
HDBHANDLE hHandle,
LONGLONG lJobID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

To get all records from the JOB_QUEUEtable.

Parameters

```
hHandle - database handlelJobID - job IDpResultCallback - pointer to a callback function to get data
```

Output

```
NULL columns:
ID
TIMESTAMP - as time updated
All other columns from JOB_QUEUE table.
```

Return

FALSE - if error occurred

MDCS_APP_GetJobQueue

```
BOOL MDCS_APP_GetJobQueue(
HDBHANDLE hHandle,
LONGLONG lJobID,
MDCS_ST_Job& stJob
);
```

Purpose

To get information from the JOB_QUEUE table using the job ID.

Parameters

```
hHandle - database connection handlelJobID - job ID
```

Output

stJob - structure that contains job information

Return

FALSE - if error occurred

MDCS_APP_UpdateJobStatus

```
BOOL MDCS_APP_UpdateJobStatus(
HDBHANDLE hHandle,
LONGLONG lJobID,
LPCSTR pszStatus
);
```

Purpose

To update job status.

Parameters

```
hHandle - database connection handlelJobID - job IDpszStatus - status
```

Output

NONE

Return

FALSE - if error occurred

MDCS_APP_UpdateJobProgress

```
BOOL MDCS_APP_UpdateJobProgress(
HDBHANDLE hHandle,
LONGLONG lJobID,
LPCSTR pszProgress
);
```

Purpose

To update job progress.

Parameters

```
hHandle - database connection handlelJobID - job IDpszProgress - Progress
```

Output

NONE

Return

FALSE - if error occurred

MDCS_APP_UpdateJobQueueRecord

```
BOOL MDCS_APP_UpdateJobQueueRecord(
HDBHANDLE hHandle,
LONGLONG lJobID,
const MDCS_ST_Job& stUpdateJob
);
```

Purpose

To update job queue information using the job ID.

Parameters

```
hHandle - database connection handlelJobID - job IDstUpdateJob - information to be updated
```

Output

NONE

Return

FALSE - if error occurred

MDCS_APP_CancelJobProgress

```
BOOL MDCS_APP_CancelJobProgress(
HDBHANDLE hHandle,
LONGLONG lJobID
);
```

Purpose

To cancel job.

Parameters

hHandle - database connection handlelJobID - job ID

Output

NONE

Return

FALSE - if error occurred

MDCS_APP_ClaimJob

```
BOOL MDCS_APP_ClaimJob(
HDBHANDLE hHandle,
LONGLONG lJobID
BOOL dDistributed = FALSE
);
```

Purpose

To claim a new job.

Parameters

hHandle - database connection handle *dDistributed* - if TRUE, analysis environment is distributed

Output

IJobID - job ID

Return

FALSE - if error occurred

MDCS_APP_ResetJob

```
BOOL MDCS_APP_ResetJob(
HDBHANDLE hHandle,
LONGLONG lJobID
);
```

Purpose

To reset job status and progress.

Parameters

hHandle - database connection handlelJobID - job ID

Output

NONE

MDCS_APP_ResetCrashedJobs

```
BOOL MDCS_APP_ResetJobProgress(
HDBHANDLE hHandle,
LPCSTR pszNetWorkID
);
```

Purpose

To reset job status and progress of all the crashed jobs from the last run.

Parameters

```
hHandle - database connection handle
pszNetworkID - network ID
```

Output

NONE

MDCS_APP_RefreshAllJobs

```
BOOL MDCS_APP_RefreshAllJobs(
HDBHANDLE hHandle
);
```

Purpose

To refresh the status of all jobs.

Parameters

hHandle - database connection handle

Output

NONE

Return

FALSE - if error occurred

MDCS_APP_GetLocationOptionByID

```
BOOL MDCS_APP_GetLocationOptionByID (
HDBHANDLE hHandle,
LONGLONG lLocationID,
MDCS_ST_LocationOption& stLocationOption
);
```

Purpose

To get a single record from LOCATION_OPTIONS table using the ID.

Parameters

hHandle - database connection handlelLocationID - location option ID

Output

stLocationOption - location option structure

Return

FALSE - if error occurred

256 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

File Data and Macro Functions 16

This chapter contains functions that you can use to work with file data and macros.

Table 16-1: File Data and Macro Functions

Function Name	Description
MDCS_FILEDATA_GetRecordByID	To find record ID by GUID
MDCS_FILEDATA_GetRecordByGUID	To find file data record ID by GUID
MDCS_FILEDATA_UpdateRecord	To update a file data record
MDCS_MACROS_GetRecord	To get macros information by macros ID
MDCS_MACROS_GetAllByGroup	To get all macros by group
MDCS_MACROS_UpdateRecord	To update a macros record
MDCS_MACROS_DeleteRecord	To delete a macros record
MDCS_MACROS_Create	To create a new macros
MDCS_MACROS_ConvertQueryResultsTo MacrosStructure	To convert the results of a query callback to a macros structure

MDCS_FILEDATA_GetRecordByID

```
BOOL MDCS_FILEDATA_GetRecordByID(
HDBHANDLE hHandle,
LONGLONG lID,
MDCS_ST_FileDataInfo* pstFileData
);
```

Purpose

Function to find the record ID by GUID.

Parameters

```
hHandle - database connection handle IID - record ID
```

Output

pstFileData - file data

Return

FALSE - if error occurred

MDCS_FILEDATA_GetRecordByGUID

```
BOOL MDCS_FILEDATA_GetRecordByGUID(
HDBHANDLE hHandle,

LPCSTR pszGUID,

MDCS_ST_FileDataInfo* pstFileData
);
```

Purpose

Function to find the file data record ID by GUID.

Parameters

hHandle - database connection handle
pszGUID - GUID

Output

pstFileData - file data

Return

FALSE - if error occurred

MDCS_FILEDATA_UpdateRecord

```
BOOL MDCS_FILEDATA_UpdateRecord(
HDBHANDLE hHandle,
LONGLONG lID,
const MDCS_ST_FileDataInfo& stFileData
);
```

Purpose

Function to update a file data record.

Parameters

```
hHandle - database connection handleIID - record IDpstFileData - file data
```

Output

NONE

Return

FALSE - if error occurred

MDCS_MACROS_GetRecord

```
BOOL MDCS_MACROS_GetRecord (
HDBHANDLE hHandle,
LONGLONG lMacrosID,
MDCS_ST_GroupMacros& stMacrosData
);
```

Purpose

Function to get macros information by macros ID.

Parameters

```
hHandle - database connection handle 
IMacrosID - macros ID
```

Output

stMacrosData - structure describes macros data

Return

FALSE - if error occurred

MDCS_MACROS_GetAllByGroup

```
BOOL MDCS_MACROS_GetAllByGroup (
HDBHANDLE hHandle,
LONGLONG lGroupID,
MDCS_GetDBResultsCCallback* pResultCallback
);
```

Purpose

Function to get all macros by group.

Parameters

```
hHandle - database connection handleIGroupID - Group ID
```

Output

pResultCallback - callback with data of all macros

Return

FALSE - if error occurred

MDCS_MACROS_UpdateRecord

```
BOOL MDCS_MACROS_UpdateRecord (
HDBHANDLE hHandle,
LONGLONG lMacrosID,
const MDCS_ST_GroupMacros& stMacrosData
);
```

Purpose

Function to update the macros record.

Parameters

```
    hHandle - database connection handle
    IMacrosID - macros ID
    stMacrosData - updated macros data
    pResultCallback - pointer to a callback function to get data
```

Output

All macro data

Return

FALSE - if error occurred

MDCS_MACROS_DeleteRecord

```
BOOL MDCS_MACROS_DeleteRecord (
HDBHANDLE hHandle,
LONGLONG lMacrosID,
);
```

Purpose

Function to delete the macros record.

Parameters

```
hHandle - database connection handle IMacrosID - macros ID
```

Return

FALSE - if error occurred

MDCS_MACROS_Create

```
BOOL MDCS_MACROS_Create(
HDBHANDLE hHandle,
const MDCS_ST_GroupMacros& stObjIn,
MDCS_ST_GroupMacros& stObjOut
);
```

Purpose

Function to create a new macros.

Parameters

hHandle - database connection handle stObjIn - macros data to insert

Output

stObjOut - created macros

Return

FALSE - if error occurred

MDCS_MACROS_ConvertQueryResultsToMacrosStructure

```
BOOL MDCS_MACROS_ConvertQueryResultsToMacrosStructure (
MDCS_QueryResults* pResults,
MDCS_ST_GroupMacros& stMacros,
BOOL bInNativeFormat = FALSE
);
```

Purpose

Function to convert results of query callback to macros structure.

Parameters

pResults - pointer to query results objectbInNativeFormat - if TRUE, the query results should be treated in their original format otherwise they will be treated as strings

Output

stMacros - converted data

Return

FALSE - if error occurred

262 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Structures, ENUMS, and Definitions

17

This chapter contains the following topics:

- Structures on page 263
- ENUMS Types on page 265
- Definitions on page 266
- Error Code Definitions on page 268

Structures

The MDCStoreUtils API uses many user-defined data types whose names and descriptions are provided in Table 17-1.

Table 17-1: Structures

Name	Description
1. MDCS_ST_COLUMNPROP	Structure of the column from TABLE_COLUMNS table (represents MDCStore data type information)
2. MDCS_ST_BlobLocation	Describes BLOB location on server
3. MDCS_ST_BlobInfo	Describes basic BLOB information
4. MDCS_ST_ShapeInfo	Describes measurement set result value
5. MDCS_ST_UserLogin	Describes user credentials to connect to the database
6. MDCS_ST_UserInfo	Describes a user
7. MDCS_ST_GroupInfo	Describes a group
8. MDCS_ST_FolderGroups	Describes group access to the folder
9. MDCS_ST_ObjGroups	Describes group access to a plate
10. MDCS_ST_Result set	Common structure to describe objects (result set) in the database
11. MDCS_ST_DBVersion	Describes the database version
12. MDCS_ST_Measurements	Describes measurements
13. MDCS_ST_Assay	Describes measurement sets
14. MDCS_ST_PlateImages	Describes a plate image
15. MDCS_ST_ShapeLine	Describes shape lines
16. MDCS_ST_ShapeLinesBlob	Describes the BLOB that contains shape lines
17. MDCS_ST_DsetAnalysis	Describes dataset analysis
18. MDCS_ST_DsetAnalysisAttr	Describes dataset analysis attributes
19. MDCS_ST_FolderInfo	Describes a folder
20. MDCS_ST_PlateInfo	Describes a plate
21. MDCS_ST_SiteImageInfo	Describes the site images
22. MDCS_ST_SeriesInfo	Describes series information

Table 17-1: Structures (cont'd)

23. MDCS_ST_Dataset	Describes dataset information
24. MDCS_ST_Site	Describes the site
25. MDCS_ST_ImageSource	Describes the image source
26. MDCS_ST_Acquisition	Describes an acquisition
27. MDCS_ST_AcquisitionProfile	Describes an acquisition profile
28. MDCS_ST_AssayRun	Describes an assay runs
29. MDCS_ST_AssayProfile	Describes an assay profile
30. MDCS_ST_LocationOption	Describes a location option
31. MDCS_ST_Job	Describes a job
32. MDCS_ST_AppParameter	Describes an application parameter
33. MDCS_ST_PlateTemplate	Describes a plate layout template
34. MDCS_ST_ConnectionInfo	Describes database connection information
35. MDCS_ST_Attribute	Describes attributes of the plate, dataset, measurement set
36. MDCS_ST_FileDataInfo	Describes the file data records
37. MDCS_ST_WellInfo	Describes well information
38. MDCS_ST_CellDescription	Describes a measurement
39. MDCS_ST_BlobInfoEx	Describes extended BLOB info
40. MDCS_ST_ScopeAttribute	Describes scope attribute of a plate, dataset or measurement set
41. MDCS_ST_GroupMacros	Describes group macros
42. MDCS_ST_Warehouse	Describes annotation warehouse
43. MDCS_ST_NormConfig	Describes normalization configuration
44. MDCS_ST_MeasurementUseStatistic	Describes measurement user statistic
45. MDCS_ST_DatasetResultInfo	Describes dataset analysis results
46. MDCS_ST_DlgParams	Describes the parameters needed in a dialog

ENUMS Types

Table 17-2: ENUM Types

Name	Description
1. MDCS_E_SQLServerType	Describes the type of database server
2. MDCS_E_DBObjectType	Describes the type of object in the database
3. MDCS_E_FileStorage	Describes the location where files should be stored
4. MDCS_E_RatioType	Describes ratio types
5. MDCS_E_Waves	Describes waves
6. MDCS_E_NormalizationType	Describes normalization types
7. MDCS_E_AverageType	Describes data average types
8. MDCS_E_ColumnType	Describes column types
9. MDCS_E_ImageType	Describes image types
10. MDCS_E_BlobType	Describes BLOB types
11. MDCS_E_ShapeType	Describes shapes
12. MDCS_E_AttributeType	Describes attribute types
13. MDSC_E_StatisticType	Describes statistical types
14. MDSC_E_CollapseType	Describes collapse types
15. MDCS_E_UserStatus	Describes user status to access data
16. MDSC_E_GroupType	Describes group types
17. MDCS_E_AccessPermissions	Describes access permissions to data
18. MDCS_E_PlateAccessPermissions	Describes access permissions to plates
19. MDCS_E_CompareOperation	Describes available comparison operations
20. MDCS_E_DatasetType	Describes dataset types
21. MDCS_E_Statistic	Describes defined statistical methods
22. MDCS_E_UserDataLocationOptions	Describes user data location options
23. MDCS_ST_DlgUsersType	Describes the dialog types to manage users
24. MDC_E_DatasetDependantType	Describes the dataset dependent type
25. MDC_E_SQLServerVersion	Describes version of the database server
26. MDCS_E_GroupObjects	Describes group objects
27. MDCS_E_GroupObjectsType	Describes group object types
28. MDCS_E_DataType	Describes type of result data
29. MDCS_E_DSFilteringGroup	Describes types of group objects that can be used in datasets
30. MDCS_E_MeasurementGroupBy	Describes measurement group by variations
31. MDCS_E_PlateElements	Describes measurement group by variations
32. MDCS_E_NormConfigType	Describes type of normalization configuration

5000957 C 265

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Table 17-2: ENUM Types (cont'd)

33. MDCS_QueryBldType	Describes type of query builder	
-----------------------	---------------------------------	--

Definitions

Table 17-3: Max Object Size Definitions

Name	Value	Description
1. MDCS_MAX_DATABASE_ID	20	Max size of the unique database ID
2. MDCS_MAX_OBJ_ID	20	Max size of the object ID
3. MDCS_MAX_COL_NAME	20	Max size of the database column name (user defined, i.e. data types)
4. MDCS_MAX_COL_NAMEEXT	128	Max size of the ext field
5. MDCS_MAX_COL_RATIOTYPE	4	Max size of the ratio type
6. MDCS_MAX_COL_WAVEDEF	4	Max type of the wavelength field
7. MDCS_MAX_COL_NORMTYPE	4	Max size of the normalization type field
8. MDCS_MAX_TABLEID	25	Max size of the table name in the database
9. MDCS_MAX_RESULT SETID	18	Max size of the result set ID
10. MDCS_MAX_COL_DESCRCOLUMN	2000	Max size of the description field
11. MDCS_MAX_COL_ORDERNUM	16	Max size of the order num field
12. MDCS_MAX_COL_TYPE	10	Max size of the type field
13. MDCS_MAX_COL_LENGTH	4	Max length of the column in database (numeric)
14. MDCS_MAX_COL_AVGTYPE	4	Max size of the column that contains average type description
15. MDCS_MAX_COL_NORMFACTOR	50	Max size of the column that contains normalization factor
16. MDCS_MAX_COL_DEFAULT	10	Max size of the column that contains defaults
17. MDCS_MAX_DSN	256	Max size of the DSN (ODBC data source name)
18. MDCS_MAX_USERNAME	50	Max size of the user name
19. MDCS_MAX_PASSW	256	Max size of the user password
20. MDCS_MAX_DESCRIPTION	256	Max size of the description field
21. MDCS_MAX_OBJ_NAME	128	Max size of the object name
22. MDCS_MAX_TABLE_NAME	25	Max size of the table name
23. MDCS_MAX_SERVER	256	Max size of the server name
24. MDCS_MAX_DBNAME	256	Max size of the database name
25. MDCS_MAX_ERR_SIZE	1024	Max size of error message
26. MDCS_MAX_CONNECTION_STRING	1000	Max size of connection string
27. MDCS_MAX_BARCODE	256	Max size of a barcode

266 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Table 17-3: Max Object Size Definitions (cont'd)

28. MDCS_MAX_GLOBAL_ID	256	Max size of a global ID
29. MDCS_MAX_WELL_NAME	25	Max size of a well name
30. MDCS_MAX_UNIT_NAME	25	Max size of a unit name
31. MDCS_MAX_COMPOUND_NAME	128	Max size of a compound name
32. MDCS_MAX_CELL_ID	128	Max size of a cell ID
33. MDCS_MAX_DESCRIPTION_LARGE	1024	Max size of a large description
34. MDCS_MAX_STORAGE_TYPE	256	Max size of a storage type
35. MDCS_MAX_LABEL	256	Max size of a label
36. MDCS_MAX_NETWORK_ID	256	Max size of network ID
37. MDCS_MAX_REQUEST	20	Max size of request
38. MDCS_MAX_PARAM_VALUE	255	Max size of parameter value
39. MDCS_MAX_SERVER_VERSION	255	Max size of database server version
40. MDCS_MAX_SERVER_REVISION	255	Max size of database server revision, i.e. service pack
41. MDCS_MAX_SERVER_TYPE	255	Max size of database server type
42. MDCS_MAX_ODBC_DRIVER_NAME	255	Max size of ODBC driver name
43. MDCS_MAX_ODBC_DRIVER_VERSION	255	Max size of ODBC driver version
44. MDCS_MAX_ODBC_VERSION	255	Max size of ODBC version
45. MDCS_MAX_ODBC_SOURCE_NAME	255	Max size of ODBC datasource name
46. MDCS_MAX_GUID	256	Max size of GUID
47. MDCS_MAX_OBJ_NAME_BIG	255	Max size of the object name

Error Code Definitions

Table 17-4: Error Code Definitions

Name	Value	Description
1. MDCS_ERR_CANCEL	2999	Operation was cancelled
2. MDCS_ERR_SUCCESS	3000	Operation completed successfully without error
3. MDCS_ERR_ODBC	3001	Error was reported by ODBC driver
4. MDCS_ERR_DATABASE_INCORRECT_PARAM	3002	Incorrect parameter was supplied to function
5. MDCS_ERR_DATABASE_OBJECT_NOT_EXIST	3003	Database object does not exist
6. MDCS_ERR_USER_PERMISSIONS	3004	User does not have sufficient privileges to access data, permission denied
7. MDCS_ERR_INCORRECT_OBJ_SIZE	3005	Result does not fit into the object
8. MDCS_ERR_APPLICATION	3006	Error happened in application function
9. MDCS_ERR_OBJECT_EXIST	3007	Cannot continue, object exists
10. MDCS_ERR_DELETE_APPLICATION_OBJECT	3008	Cannot delete application object
11. MDCS_ERR_CANNOT_CONVERT_DATA	3009	Cannot perform data conversion
12. MDCS_ERR_CANNOT_SAVE_DATA	3010	Cannot save data
13. MDCS_ERR_EMPTY_QUERY	3011	Cannot execute empty query
14. MDCS_ERR_INCORRECT_QUERY	3012	Query is incorrect

Virtual Callback Classes 18

In order to provide a consistent way to work with data, MDCStoreUtils provides the following virtual callback classes to users so that they can use them as an interface to define their own class definitions when using callbacks:

- Class MDCS_GroupInfoCallback on page 269
- Class MDCS_SaveBlobCallback on page 271
- Class MDCS_GetBlobCallback on page 273
- Class MDCS_ProgressCallback on page 275
- Class MDCS_GetDBResultsCCallback on page 277
- Class MDCS_QueryResults on page 279
- Class MDCS_ImportMeasurementSet on page 282
- Class MDCS_ImportPlateLayout on page 285
- Class MDCS CL ImportDS on page 290
- Class MDCS_GetProgressStatus on page 291
- Class MDCS CL BlobLocationCB
- class MDCS CL BlobInfoCB

MDCStoreUtils also provides the MDCS_DBHandleSmartPtr class to deal with database handles.

Class MDCS_GroupInfoCallback

Overview

To retrieve multiple MDCS_ST_GroupInfo structure results which describes group information from database.

Public Methods

MDCS_GroupInfoCallback	Default constructor
~MDCS_GroupInfoCallback	Virtual destructor
ItemCount	Pure virtual function. To get the total number of groups.
GetNextInfo	Pure virtual function. To get group information from database record.

Private Methods

MDCS_GroupInfoCallback	Copy constructor
operator=	Assignment operator

Function signature of MDCS_GroupInfoCallback

```
MDCS_GroupInfoCallback () {};
```

Default constructor to construct MDCS_GroupInfoCallback object.

```
virtual ~MDCS GroupInfoCallback() {};
```

Virtual destructor destroys the MDCS_GroupInfoCallback object.

```
virtual ItemCount(LONGLONG lGroupCount) = 0;
```

To get a total number of groups.

Parameters

IGroupCount - number of groups counted

Return

Total number of groups.

```
virtual GetNextInfo(const MDCS_ST_GroupInfo& stGroupInfo)
= 0:
```

To get group information.

Parameters

stGroupInfo - A structure describes GroupInfo

Output

stGrouptInfo - A structure of MDCS_ST_GroupInfo contains information
of the group

```
MDCS GroupInfoCallback (MDCS GroupInfoCallback &);
```

The default copy constructors and it is private.

```
const MDCS_GroupInfoCallback
&operator=(MDCS GroupInfoCallback const &);
```

Assignment (=) operator reinitializes an existing MDCS_GroupInfoCallback object with new data.

Class MDCS_SaveBlobCallback

Overview

MDCS_SaveBlobCallback - callback for saving data blobs.

Public Methods

MDCS_SaveBlobCallback	Default constructor
~MDCS_SaveBlobCallback	Virtual destructor
NextChunk	Pure virtual function. To get the next chunk of data.
Done	Pure virtual function. Called when completely done.
GetPacketSize	Pure virtual function. To get size of transfer package.
GetTotalSize	Pure virtual function. To get total size of data transfer.
Error	Pure virtual function. Sent when there is an error.

Private Methods

MDCS_SaveBlobCallback	Copy constructor
operator=	Assignment operator

Function signatures of MDCS_SaveBlobCallback

```
MDCS_SaveBlobCallback () {};
Default constructor to initialize the MDCS_SaveBlobCallback.
virtual ~MDCS_SaveBlobCallback() {};
Virtual destructor
virtual BOOL NextChunk(
BYTE* pChunk,
const UINT& uSizeChunk,
UINT& uSizeRead) = 0;
```

Purpose

To get the next chunk of data.

Parameters

```
pChunk - pointer to chunk of data.uSizeChunk - size that was requesteduSizeRead - actual size that was read
```

Output

```
pChunk - pointer to next chunk of data
virtual void Done() = 0;

Call this function when the saving process is completely done.
virtual UINT GetPacketSize() const = 0;

To get the size of the transfer package.
virtual LONGLONG GetTotalSize() const = 0;

To get the total size of data transferred.
virtual void Error(LPCSTR pErrorText) = 0;

Sent error text when there is an error.
```

Parameter

pErrorText - pointer to error text

Output

```
pErrorText - pointer to the sent error text
MDCS_SaveBlobCallback (MDCS_SaveBlobCallback &);
A default copy constructor.
const MDCS_SaveBlobCallback
&operator=(MDCS SaveBlobCallback const &);
```

Assignment operator reinitializes the MDCS_SaveBlobCallback object to new data.

Class MDCS_GetBlobCallback

Overview

MDCS_GetBlobCallback - callback class for getting data blobs.

Public Methods

MDCS_GetBlobCallback	Default constructor
~MDCS_GetBlobCallback	Virtual destructor
NextChunk	Pure virtual function. To get the next chunk of data received from the database.
NextResult	Pure virtual function. Return size of the result.
GetPacketSize	Pure virtual function. To get the size of transfer package.
GetTotalSize	Pure virtual function. To get total size of data.
Done	Pure virtual function. Called when 100 percent complete.
Error	Pure virtual function. Sent when there is an error.

Private Methods

MDCS_GetBlobCallback	Copy constructor
operator=	Assignment operator

Function signature of MDCS_GetBlobCallback

```
MDCS_GetBlobCallback () {};
Default constructor
virtual ~MDCS_GetBlobCallback () {};
Virtual destructor
virtual BOOL NextResult(LONGLONG lSizeOfResult) = 0;
To get the size of the result.
```

Parameter

ISizeOfResult - Signed 64-bit integer to hold the result size.

Output

```
ISizeOfResult - the result size
```

```
virtual BOOL NextChunk(const BYTE* pChunk, UINT
uChunkSize) = 0;
```

To get the next chunk of data that will be received from database or file server.

Parameters

```
pChunk - pointer to the data
  uChunkSize - size
```

Output

```
pChunk - pointer to the retrieved data.
virtual LONGLONG GetTotalSize() const = 0;
To get total size of data.
virtual UINT GetPacketSize() const = 0;
To get size of the transfer packet.
```

virtual void Done() = 0;

Call this function when getting data completely done.

virtual void Error(LPCSTR pErrorText) = 0;

The error text is sent when there is an error.

Parameter

pErrorText - pointer to error text

Output

```
pErrorText - pointer to the sent error text

MDCS_GetBlobCallback (MDCS_GetBlobCallback const &);
The default copy constructors.
const MDCS_GetBlobCallback
&operator=(MDCS_GetBlobCallback const &);
```

Assignment = operator reinitializes the MDCS_GetBlobCallback object to a new data.

Class MDCS_ProgressCallback

Overview

An abstract class of encapsulated functions to work with progress callback.

Public Methods

MDCS_ProgressCallback	Default constructor
~MDCS_ProgressCallback	Virtual destructor
Progress	Pure virtual function. To get percentage done and tell how many percent done.
Done	Pure virtual function. Called when 100 percent complete.
Error	Pure virtual function. Sent when there is an error.

Private Methods

MDCS_ProgressCallback	Copy constructor
operator=	Assignment operator

Function signature of MDCS_ProgressCallback

```
MDCS ProgressCallback() {};
```

The default constructor to construct the MDCS_ProgressCallback object virtual ~MDCS_ProgressCallback() {};

Virtual destructor to destroy dynamic data members.

```
virtual BOOL Progress(
int nPercentDone,
LPCSTR pProgressText,
BOOL bComplete) = 0;
```

To determine percentage progress.

Parameters

nPercentDone - between 0 and 100 pProgressText - string telling what is being done bComplete - TRUE if 100% is really 100%

Output

```
bProgressText - string telling what is being done
virtual void Done()=0;
```

Called when completely done.

```
virtual BOOL Error(
LPCSTR pErrorText,
BOOL bWithContinue = FALSE) = 0;
```

The error text is sent when there is an error. If *bWithContinue* is TRUE, error is not critical and process can be continued.

Parameters

pErrorText - pointer to error text
bWithContinue - to indicate whether the process can be continued

Output

pErrorText - pointer to error text

Return

```
FALSE - exist
TRUE - continue
MDCS_ProgressCallback(MDCS_ProgressCallback const &);
The default copy constructors.
const MDCS ProgressCallback
```

&operator=(MDCS_ProgressCallback const &);

Assignment = operator reinitializes object MDCS_ProgressCallback to new data.

Class MDCS_GetDBResultsCCallback

Overview

An abstract class of encapsulated functions to work with database result callback.

Public Methods

MDCS_GetDBResultsCCallback	Default constructor
~MDCS_GetDBResultsCCallback	Virtual destructor
GetResultInfo	Pure virtual function. Called to get information about result structure (field names and types), all fields will be empty.
GetNextResult	Pure virtual function. Called to get next row of data.
Done	Pure virtual function. Called when 100 percent complete.
Error	Pure virtual function. Called when there is an error.
ProcessResultsInOriginalFormat	Pure virtual function. Called by function to check if results should be return in their original format.
SetResultNumber	Pure virtual function. To pass back to client number of results to be retrieved. Called to get number of results.

Private Methods

MDCS_GetDBResultsCCallback	Copy constructor
Operator=	Assignment operator

Function signature of MDCS_GetDBResultsCCallback

```
MDCS GetDBResultsCCallback() {};
```

Default constructor to construct the MDCS_GetDBResultsCCallback.

```
virtual ~MDCS GetDBResultsCCallback() {};
```

Virtual destructor

```
virtual BOOL GetNextResult(MDCS_QueryResults* pQueryRes)
= 0;
```

To get next result of data.

Parameter

pQueryRes - pointer to object MDCS_QueryResults

Output

```
pQueryRes - pointer to object MDCS_QueryResults
virtual void SetResultNumber(LONGLONG nCount) = 0;
To get number of results.
```

Parameter

nCount - Signed 64-bit integer to hold the number of result.

Output

```
nCount - the number of result
```

```
virtual void Error(LPCSTR pErrorText) = 0;
```

To send the error text when an error occurred.

Parameter

pErrorText - pointer to error text

Output

```
pErrorText - pointer to error text
```

```
virtual\ void\ Done() = 0;
```

Called when completely done.

```
virtual BOOL ProcessResultsInOriginalFormat() const = 0;
```

Called by function to check if results should be returned in their original format.

```
virtual BOOL GetResultInfo(MDCS QueryResults* pQueryRes);
```

Called to get information about results structure (field names and types).

All fields will be empty.

Parameter

pQueryRes - pointer to MDCS_QueryResults

Output

pQueryRes - pointer to MDCS_QueryResults

Return

TRUE always

```
MDCS_GetDBResultsCCallback (MDCS_GetDBResultsCCallback
const &);
```

The default copy constructor.

```
const MDCS_GetDBResultsCCallback
&operator=(MDCS GetDBResultsCCallback const &);
```

Assignment operator reinitializes the object MDCS_GetDBResultsCCallback to new data.

Class MDCS_QueryResults

Overview

An abstract class of encapsulated functions to work with QueryResults.

Public Methods

MDCS_QueryResults	Default constructor
~MDCS_QueryResults	Virtual destructor
GetStringValue	Pure virtual and Overloaded function. Function to get string result value of the field.
GetFloatValue	Pure virtual and Overloaded function. Function to get float result value of the field.
GetLongValue	Pure virtual and Overloaded function. Function to get integer result value of the field.
GetFieldType	Pure virtual and Overloaded function. Return data type of column.
GetColumnName	Pure virtual function. Returns name of the column.
GetColumnsCount	Pure virtual function. Returns number of columns.

Private Methods

MDCS_QueryResults	Copy constructor
operator=	Assignment operator

Function signature of MDCS_QueryResults

```
MDCS QueryResults() {};
```

Default constructor to initialize MDCS_QueryResults object.

```
virtual ~MDCS QueryResults(){};
```

Overridable destructor.

virtual LPCSTR GetStringValue(int nCol) const = 0;

Overloaded function to get string result value of the field.

Parameter

nCol - column number

Output

A pointer to a string result.

virtual LPCSTR GetStringValue(LPCSTR strColumnName) const
= 0;

Overloaded function to get string result value of the field.

Parameter

strColumnName - a pointer to column name of the field of result

Return

A pointer to string value result.

```
virtual float GetFloatValue(int nCol) const = 0;
```

Overloaded function to get float result value of the field.

Parameter

nCol - column number or field number of result

Return

A float value result.

```
virtual float GetFloatValue(LPCSTR strColumnName) const =
0:
```

Overloaded function to get float result value of the field.

Parameter

strColumnName - pointer to column name or field name.

Return

A float result value of the field.

```
virtual LONGLONG GetLongValue(int nCol) const = 0;
```

Overloaded function to get integer value of the field.

Parameter

nCol - column number, field number of result

Return

A signed integer value result of the field.

```
virtual LONGLONG GetLongValue(LPCSTR strColumnName) const
= 0;
```

Overloaded function to get integer result value of the field.

Parameter

strColumnName - pointer to a column name or field name

Return

A signed integer value that is result of the field.

```
virtual MDCS_E_ColumnType GetFieldType(int nCol) const =
0;
```

Overloaded function to get type of column.

Parameter

nCol - column number or field number

Return

Type of the column.

```
virtual MDCS_E_ColumnType GetFieldType(LPCSTR
strColumnName) const = 0;
```

Overloaded function to get type of column.

Parameter

strColumnName - pointer to column name

Return

The enum field type of column.

```
virtual LPCSTR GetColumnName(int nCol) const = 0;
```

To get name of the column.

Parameter

nCol - column number or field number

Return

Name of the column.

```
virtual int GetColumnsCount() const = 0;
```

To get number of columns.

Class MDCS_ImportMeasurementSet

Overview

A data source to import measurement set data into the database.

Public Methods

MDCS_ImportMeasurementSet	Default constructor
~MDCS_ImportMeasurementSet	Virtual destructor
GetMeasurementAttributeNumber	Pure virtual function. To get the number of measurement attributes that applied to a cell, not including ones that are in MDCS_ST_CellDescription.
GetCellCount	Pure virtual function. To get the total number of rows to be inserted.
GetMeasurementAttributeInfo	Pure virtual function. To get measurement column information.
GetCellDescription	Virtual function. To get the cell description for each row that we insert.
GetMeasurementValue	Overloaded and pure virtual function. To get measurement set value to insert.
GetNumberRowsInStep	Pure virtual function. To get max number of rows to be inserted in one step.
StatusWhenDone	Virtual function. To set assay status when import done.
StatusWhileProccessing	Virtual function. To set a status of assay while importing.

Private Methods

MDCS_ImportMeasurementSet	Copy constructor
operator=	Assignment operator

Function signatures of MDCS_ImportMeasurementSet

```
MDCS ImportMeasurementSet() {};
```

Default constructor to initialize MDCS_ImportMeasurementSet object.

```
virtual ~MDCS ImportMeasurementSet() {};
```

Virtual destructor.

virtual UINT GetMeasurementAttributeNumber() const = 0;

To get number of measurement attributes that applied to a cell, not including ones that are in MDCS_ST_CellDescription.

```
virtual BOOL GetMeasurementAttributeInfo(
UINT uColumnNum,
LONGLONG& nAttributeID) = 0;
```

To get measurement column information.

282 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Parameters

uColumnNum - Column numbernAttributeID - an ID of a global measurement attribute

Output

```
Measurement column information.
```

```
virtual UINT GetCellCount() const = 0;
To get total number of rows to be inserted.
virtual BOOL GetCellDescription(
UINT uCellIndex,
MDCS_ST_CellDescription& stCellDescription) = 0;
```

To get cell description for each row that we insert.

Parameters

```
nCellIndex - index of the cell stCelDescription - structure of Cell description
```

Output

A structure of cell description.

```
virtual BOOL GetMeasurementValue(
UINT uCellIndex,
UINT uAttrIndex,
LPSTR pszData,
UINT uSize) const = 0;
```

To get measurement set value to insert.

Parameters

```
uCellIndex - index of cell
uAttrIndex - attribute index
pszData - pointer to measurement data
uSize - size of the data, use MDCS_MAX_STRING_VALUE for it.
virtual BOOL GetMeasurementValue(
UINT uCellIndex,
UINT uAttrIndex,
float** pfData) const = 0;
```

5000957 C 283

To get float measurement set value to insert.

Parameter

```
pfData - pointer to float value that will be inserted
virtual BOOL GetMeasurementValue(
```

```
UINT uCellIndex,
UINT uAttrIndex,
int** pnData) const = 0;
```

To get integer measurement set value to insert

Parameter

pnData - pointer to int value that will be inserted

```
virtual UINT GetNumberRowsInStep() const = 0;
```

To get max number of rows to be inserted in one step.

```
MDCS_ImportMeasurementSet (MDCS_ImportMeasurementSet const &);
```

Copy constructor.

```
const MDCS_ImportMeasurementSet
&operator=(MDCS_ImportMeasurementSet const &);
```

Assignment operator reinitializes MDCS_ImportMeasurementSet object to new data.

```
virtual MDCS_E_AssayMarkStatus StatusWhenDone(
) const { return eMarkNormal;};
To set assay status when import is done.
virtual MDCS_E_AssayMarkStatus StatusWhileProccessing(
) const { return eMarkRemove;};
To set a status of assay while importing.
```

Class MDCS_ImportPlateLayout

Overview

A data source to import plate layout data to the database.

Public Methods

MDCS_ImportPlateLayout	Default constructor
~MDCS_ImportPlateLayout	Virtual destructor
GetNumberWells	Pure virtual function. To get the number of columns in the data source.
GetNumberGroups	Pure virtual function. To get the number of groups on a plate.
GetNumberCompounds	Pure virtual function. To get the number of compounds on a plate.
GetWellInfo	Pure virtual function. To get the well description.
GetCompound	Pure virtual function. To get the compound description.
GetConcentration	Pure virtual function. To get the concentration of compound.
GetUnit	Pure virtual function. To get the unit data.
GetWellGroup	Pure virtual function. To get the well group.

Private Methods

MDCS_ImportPlateLayout	Copy constructor.
operator=	Assignment operator

Function signature of MDCS_ImportPlateLayout

```
MDCS_ImportPlateLayout() {};
Default constructor
virtual ~MDCS_ImportPlateLayout() {};
Virtual destructor
virtual UINT GetNumberWells() const = 0;
To get the number of columns in the datasource.
virtual UINT GetNumberGroups() const = 0;
To get the number of groups on a plate.
virtual UINT GetNumberCompounds() const = 0;
To get number of compounds on a plate.
virtual BOOL GetWellInfo(
UINT uWellIndex,
MDCS_ST_WellInfo& stWellInfo) const = 0;
To get the Well description.
```

Parameters

uWellIndex - well index
stWellInfo - Well information structure

Output

stWellInfo - well information structure

```
virtual BOOL GetCompound(
UINT uWellIndex,
UINT uCompoundIndex,
LPSTR pszData,
UINT uSize) const = 0;
```

To get the compound description.

Parameters

uWellIndex - well index of compound
 uCompoundIndex - compound index
 pszData - pointer to string of compound description
 uSize - size of the compound description string

Output

```
pszData - pointer to compound description
virtual BOOL GetConcentration(
UINT uWellIndex,
UINT uCompoundIndex,
float** fConc) const = 0;
```

To get concentration of a compound in a well.

Parameters

```
uWellIndex - well indexuCompoundIndex - compound indexfConc - pointer to float value to be inserted
```

Output

Concentration of the compound.

```
virtual BOOL GetUnit(
UINT uWellIndex,
UINT uCompoundIndex,
LPSTR pszData,
UINT uSize) const = 0;
To get unit of data.
```

Parameters

```
uWellIndex - well indexuCompoundIndex - compound indexpszData - pointer to datauSize - size of data string
```

Output

Unit of data.

```
virtual BOOL GetWellGroup(
UINT uWellIndex,
UINT uGroupIndex,
LPSTR pszGroupName,
UINT uSize) const = 0;
```

To get the Well group description.

Parameters

uWellIndex - well index
uGroupIndex - group index
pszGroupName - group name
usize - size of group name string

Output

pszGroupName - group name description

MDCS ImportPlateLayout (MDCS ImportPlateLayout const &);

Default copy constructor.

const MDCS_ImportPlateLayout
&operator=(MDCS ImportPlateLayout const &);

Assignment operator reinitializes MDCS_ImportPlateLayout object to new data.

Class MDCS_DBHandleSmartPtr

Overview

MDCS_DBHandleSmartPtr is not a pure virtual callback class. It is a smart pointer that is a wrapper around HDBHANDLE to deal with HDBHandles. It initializes and destroys dbhandle, so you do not need to destroy a database handle after use. However, it can be used only on the stack.

Data Member

m_hdbHandle	A HDBHANDLE database handle.

Public Methods

MDCS_DBHandleSmartPtr	Constructor to initialize user login information.
~MDCS_DBHandleSmartPtr	Virtual Destructor
GetHandle	Return the HDBHANDLE database handle.

Private Methods

MDCS_DBHandleSmartPtr	Default constructor
operator=	Overridden = assignment operator.

288 5000957 C

The MDCStore™ SDK and associated documentation are copyrighted and property of Molecular Devices, LLC. Use of the MDCStore™ API (Application Programming Interface) and acceptance of the associated documentation is prohibited without the express written consent of Molecular Devices, LLC, which consent shall limit the scope of use of the MDCStore™ technology, shall be personal to the licensed user, and shall be non-transferable.

Function signature of MDCS_DBHandleSmartPtr

```
MDCS_DBHandleSmartPtr() {};
Default constructor to construct the MDCS_DBHandleSmartPtr object.
MDCS_DBHandleSmartPtr& operator=(
const MDCS_DBHandleSmartPtr &stConnection
):
```

Assignment operator reinitializes the MDCS_DBHanldeSmartPtr to a new data.

```
MDCS_DBHandleSmartPtr(const MDCS_ST_UserLogin&
stUserLogin,
BOOL bSilentMode = FALSE, BOOL bConnectAsAppUser = TRUE);
MDCS_DBHandleSmartPtr(LPCTSTR pszConnectionString,
BOOL bSilentMode = FALSE);
```

Constructor to initialize the user login information.

Parameters

bSilentMode - Indicates if connection should be in silent mode. *bConnectAsAppUser* - Indicates if the user is an application user.

```
virtual ~MDCS_DBHandleSmartPtr();
Virtual destructor.
HDBHANDLE GetHandle();
```

Return the database handle of this object.

Class MDCS_CL_ImportDS

Public Methods

MDCS_CL_ImportDS	Constructor to initialize measurement column database name.
GetAttributeDBName	Virtual function. To get measurement column database name.
GetTableName	Virtual function. To get measurement table name.

Function signature of MDCS_CL_ImportDS

To get Attribute database name

virtual LPCTSTR GetAttributeDBName(UINT uColumnNum) const
= 0;

Parameters

uColumnNum - column number

Output

Attribute Name

To get attribute table name

virtual LPCTSTR GetTableName() const = 0;

Parameters

None

Returns

Table name

} ;

Class MDCS_GetProgressStatus

Public Methods

MDCS_GetProgressStatus	Default constructor
~MDCS_GetProgressStatus	Virtual destructor
GetProgress	Virtual function to get current progress.
GetProgressText	Virtual function to get current progress text.
GetTitleText	Virtual function to get current progress title text.
GetStatus	Virtual function to get a status.
GetError	Virtual function to get an error description.

Private Methods

MDCS_GetProgressStatu	To hide the default copy constructors.
S	

Function signature of MDCS_GetProgressStatus

```
MDCS GetProgressStatus() {};
Default constructor to initialize MDCS_GetProgressStatus object.
virtual ~MDCS GetProgressStatus() {};
To get current progress.
virtual UINT PTR GetProgress() const = 0;
To get current progress.
virtual LPCSTR GetProgressText() const = 0;
To get current progress.
virtual LPCSTR GetTitleText() const = 0;
To get current progress.
virtual E_ProgressStatus GetStatus() const = 0;
To get a status.
virtual LPCSTR GetError() const = 0;
To get an error description.
MDCS_GetProgressStatus (MDCS_GetProgressStatus const &);
      const MDCS GetProgressStatus
      &operator=(MDCS GetProgressStatus const &);
To hide the default copy constructors
```

Class MDCS_CL_BlobLocationCB

This is a class to process an object count.

Public Methods

MDCS_CL_BlobLocationCB	Constructor to initialize object count.
~MDCS_CL_BlobLocationCBr	Virtual Destructor

Private Methods

MDCS_CL_BlobLocationCB	To hide the default copy constructor.
operator=	Overridden = assignment operator.

Return

virtual void NextResult(const MDCS ST BlobLocation&) = 0;

class MDCS_CL_BlobInfoCB

This is a class to process blob information.

Public Methods

MDCS_CL_BlobInfoCB	Constructor to initialize object count.
~MDCS_CL_BlobLocationCBr	Virtual Destructor

Private Methods

MDCS_CL_BlobInfoCB	To hid the default copy constructor
operator=	Overridden = assignment operator.

Return

virtual void NextResult(const MDCS_ST_BlobInfo&) = 0;

Usage Examples 19

This chapter provides some examples of how to use functions in MDCStoreUtils API to do the following:

- Make a connection to a datasource and get a database handle.
- Get plate, site, and image information.
- Show database error messages.
- Derive from pure virtual callback classes.

How to Connect to the Database and Get a Database Handle

In order to work with the MDCStoreTM database, you need to have a database handle, which is a parameter required in most of the functions in MDCStoreUtils API. You can get a database handle easily by calling the function **MDCS_CONNECTION_GetDBHandle**. In order to call the function, you need to have user login information to pass to this function to create a database connection and return a database handle. When filling the user login information, the following members of MDCS_ST_UserLogin need to be specified:

- szUserName
- szPassword
- szDSN
- szDatabase optional. It is used when an ODBC datasource was set up without indicating database name.

You can fill out the information above manually or by calling function MDCS_CONNECTION_GetDetails.

The function **MDCS_CONNECTION_GetDetails** displays a dialog where you can enter the login information.



Note: Comments within the code examples are shown in a green font.

The following example shows how to get user information and a database handle using these functions:

```
MDCS ST UserLogin stUserLogin ;
//holds user login information
   HDBHANDLE hDBHandle;
//database handler
//Call the following function to gather login information
from user and
//store information in stUserLogin structure.
if (!MDCS CONNECTION GetDetails (&stUserLogin, NULL, NULL, "MD
CS Example", NULL, FALSE) )
  {
//show database error, will discuss in next section
      ShowDBDllError("Function Fail.\nReason:");
//implemented by user, its definition is below
//Another way is to initialize the user name, password and
datasource name in stUserLogin
//manually and use it to connect to database.
//For example,
//strcpy(stUserLogin.szUserName, "sa");
//strcpy(stUserLogin.szPassword, "sa");
//strcpy(stUserLogin.szDSN, "IX MDCStoreDemo");
   cout << "Next! We are getting database handle" << endl;</pre>
 //Call function in MDCStoreUtilApi to get database
handle
   hDBHandle = MDCS CONNECTION GetDBHandle(stUserLogin);
```

Another way to get a database handle is to use the MDCS_DBHandleSmartPtr object. Using this object you don't have to destroy the database handle when you are finished using it because the destructor of MDCS_DBHandleSmartPtr takes care of it. This is demonstrated in the ExecMain function below:

```
void ExecMain(const MDCS ST UserLogin& stUserLogin,
LONGLONG lPlateID)
//get database handle
   MDCS DBHandleSmartPtr ptrDB(stUserLogin);
//check if we got a valid connection
   if (ptrDB.GetHandle())
//get images for a plate that have ID 18
      BOOL bReturn = GetImages(ptrDB.GetHandle(),
lPlateID);
      if (!bReturn )
         ShowDBDllError("Function failed.\nReason:");
   else
      cout << "Execution cancelled." << endl;</pre>
   }
   int ch;
   cputs( "\n\nPress any key to exit." );
   ch = getch();
   _cputs( "\r\n" );
}
```

How to Use Error Handling

There is a function call to **ShowDBDIIError** to display an error message in the first example above. **ShowDBDIIError** is a typical way to use the error-handling functions in MDCStoreUtils API to display the database error whenever a MDCStoreUtils API function fails. **ShowDBDIIError** calls **MDCS_GetLastErrorMsg** to get the last database error and displays the error in a console window. Here is an example of how to call **MDCS_GetLastErrorMsg** to do error handling:

```
BOOL ShowDBDllError(const CString& strErrorExtra )
{
    char szError[MDCS_MAX_ERR_SIZE] = "";
    if(MDCS_GetLastErrorMsg(szError, MDCS_MAX_ERR_SIZE) ==
MDCS_ERR_SUCCESS)
        return FALSE;
    else
    {
        CString strErr;
        strErr.Format("%s\n%s", strErrorExtra, szError);
        cout << strErr << endl;
        return TRUE;
    }
    return FALSE;
}</pre>
```

How to Get Plate Information

Once you have a database handle, you can access any data in the MDCStore database, such as plate information, site information belonging to a specific plate, or an image from the database. Below is a function that calls **MDCS_PLATE_GetInfo** to get plate information:

```
//FUNCTION: GetPlateInfo
//PURPOSE: Get plate information and print out the plate
information
//Parameters
     hDBHandle - database handle
     lPlateID - a plate ID
//RETURN Return false - if fails
//-----
BOOL GetPlateInfo(HDBHANDLE hDBHandle, LONGLONG
lPlateID)
  MDCS ST PlateInfo st PlateInfo;
//to hold the plate information
  if (!MDCS PLATE GetInfo(hDBHandle, lPlateID,
st PlateInfo))
   {
//You can show DATABASE error here by calling
//ShowDBDllError("Function fails."); or let the function,
where GetPlateInfo
//function is called handle the error, or the main
program.
     return FALSE;
   }
  else
   {
// Print out plate information
     cout << "The plate information " << endl;</pre>
     cout << "Plate name:\t" << st PlateInfo.szPlateName</pre>
<< endl;
     cout << "Xwells:\t" << st PlateInfo.nXWells <<</pre>
endl;
     cout << "Ywells:\t" << st PlateInfo.nYWells << endl;</pre>
```

```
cout << "Acquisition:\t" << st_PlateInfo.szAcqName
<< endl;
    cout << "Barcode:\t" << st_PlateInfo.szBarcode <<
endl;
    cout << "Creator:\t" << st_PlateInfo.szCreator <<
endl;
    cout << "Description:\t" << st_PlateInfo.szDesc <<
endl;
    cout << "N\n";
}
return TRUE;
}</pre>
```

In the previous example the function returns FALSE when any function call fails, and we did not display the error here. You can handle the error here by calling **ShowDBDIIError**, but we intend to handle the error in the function that calls **GetPlateInfo**, so we simply return false if there is an error. It is up to you to decide where to handle the error. If the error is handled in this function then you do not need to handle the error in the function that called it.

How to Get All Sites Per Plate

MDCStoreUtils API provides the function MDCS_PLATE_GetSitesByPlate to get all sites per plate. It uses a callback function GetNextResult in the pure virtual class MDCS_GetDBResultsCCallback to process the result. Therefore, first you need to derive the MDCS_GetDBResultsCCallback class and provide the definition to the GetNextResult method to process results. Let's create the CGetSitesByPlate_Callback callback class inherited from MDCS_GetDBResultsCCallback to process and get all site IDs of a plate returned from the database by the MDCS_PLATE_GetSitesByPlate function. Deriving virtual callback classes is discussed in next section.

Deriving the MDCS_GetDBResultCCallback Class

When deriving any pure virtual callback class in MDCStoreUtils, the way that you implement the class depends on what you intend to get from the database. The data returned from the database determines the data member variable to hold the result in the derived class.

The function **MDCS_PLATE_GetSitesByPlate** in MDCStoreUtils API gets all sites per plate. It returns records with columns:

What is needed is only site IDs belong to the plate, so in this case the data member is an array of siteIDs **m_plarrSiteIDs** (see below in GetSitesByPlate_CallBack.h). The function **GetNextResult** (see below in GetSitesByPlate_CallBack.cpp) gets the next row of data. When implementing this function, you just need to call the right function of pointer **pQueryRes** that points to the **MDCS_QueryResults** object to get the data you want and to put it in the data member variable.

CGetSitesByPlate_CallBack.h and **GetSitesByPlate_CallBack.cpp** below show how to implement the derived class in this case.

```
//***************
*****
// Copyright (c) 2005 Molecular Devices
// All rights reserved.
//***************
*****
// Module: GetSitesByPlate CallBack.h
// Purpose: Provide functions to process data and get all
Site IDs belonging to a plate
//
#pragma once
#include "mdcstoreutilsapi.h"
#include <afxtempl.h>
class CGetSitesByPlate CallBack :
  public MDCS GetDBResultsCCallback
```

```
{
public:
   CGetSitesByPlate CallBack(CArray<LONGLONG, LONGLONG>*
arrSiteID);
   virtual ~CGetSitesByPlate_CallBack(void){}
   //called to get next row of data
   virtual BOOL GetNextResult(MDCS QueryResults*
pQueryRes);
   //called to get number of results
   virtual void SetResultNumber(LONGLONG nCount)
   { };
   //called when there is an error
   virtual void Error(LPCSTR pErrorText)
      return;
   // called when completely done
   virtual void Done()
      return;
   //called to check if results should be returned in
their original format
   virtual BOOL ProcessResultsInOriginalFormat() const
```

Using the Callback CGetSitesByPlate_Callback

Once you have the callback class **CGetSitesByPlate_CallBack**, you are ready to call the **MDCS_PLATE_GetSitesByPlate** function to get all siteIDs belonging to a plate. The following function is an example to illustrate how to get all siteIDs belonging to a specific plate:

Deriving the MDCS_GetBlobCallback Class

Similar to getting siteIDs per plate above, getting images requires that you provide a callback class. In this case the callback class is derived from **MDCS_GetBlobCallback**. Below is a typical way to implement the callback class you need. The **CBlobCallbackEx** class allows writing blob data directly to file using a file handle. The most important method in this class is **NextChunk**. It is used to get the next chunk of data from the source (database or file server) and write it to a file directly.

```
//***************
// Copyright (c) 2005 Molecular Devices
// All rights reserved.
//****************
*****
//Module: CBlobCallbackEx.cpp
//PURPOSE: Extension of the class CBlobCallback, it
        allows you to save data right into a file
#include "stdafx.h"
#include ".\blobcallbackex.h"
// Constructor to initialize the object
CBlobCallbackEx::CBlobCallbackEx( HANDLE file, LONGLONG
uTotalSize, LPCSTR pszFileName, UINT uChunkSize /*=
DEFAULT CHUNK SIZE */) :
   m file(file), m uChunkSize(uChunkSize),
m uTotalSize(uTotalSize), m strFileName(pszFileName)
{
  m lCurSize = 0;
// FUNCTION: NextResult
// PURPOSE: Next result is retrieved, returns size of the
result
BOOL CBlobCallbackEx::NextResult(LONGLONG lResultSize)
  return TRUE;
}
```

302 5000957 C

```
//============// FUNCTION: NextChunk
// PURPOSE: Next chunk of data that will be received from database
//
// PARAMETERS
// pChunk - data
// uChunkSize - size of data to be retrieved
// RETURN : FALSE - if error occurs
//==========BOOL CBlobCallbackEx::NextChunk(const BYTE* pChunk, UINT uChunkSize )
```

Using the CBlobCallbackEx Callback Class

Having a callback class CBlobCallbackEx ready, you can call the function **MDCS_BLOB_Get** in MDCStoreUtils to get image data.

```
MDCS_BLOB_Get(HDBHANDLE hHandle,
LONGLONG lBlobID,
MDCS_E_BLobType pBlobType,
MDCS GetBlobCallback * pCallBack)
```

The **GetImages** function below illustrates the use of the **CBlobCallbackEx** callback class and the **MDCS_BLOB_Get** function to write image data to a file on the C drive.

5000957 C 303

```
CString strFileNameTemplate2 = "Image_%d.jpg";
    CString strFileName;
//get each image from database to C drive
   for ( int nCount = 0; nCount < larrImageIDs.GetSize();</pre>
nCount++ )
//get image information and use information about
original file extension to
//create files on disk
      MDCS ST BlobInfo stBlobInfo;
      if (!MDCS_BLOB_GetInfo(hDBHandle, &stBlobInfo,
                              larrImageIDs.GetAt (nCount),
MDCS eBlobSiteImage))
         return FALSE;
      if (!stBlobInfo.lBlobID)
         cout << "Image with ID " <<
larrImageIDs.GetAt(nCount) << " is not</pre>
                 found in the database" << endl;
         continue;
      CString strExt;
      if (!strlen(stBlobInfo.szImportExtension))
         strFileName.Format(strFileNameTemplate2,nCount);
      else
         strFileName.Format(strFileNameTemplate, nCount,
stBlobInfo.szImportExtension);
HANDLE hFile;
      //create image file
      hFile = CreateFile(strFileName,
                      FILE WRITE DATA,
                      FILE SHARE WRITE,
                      NULL,
                      CREATE ALWAYS,
                      FILE ATTRIBUTE_NORMAL,
                      NULL);
```

```
if (hFile == INVALID HANDLE VALUE)
      {
         cout << "Could not open file." << endl;</pre>
      CBlobCallbackEx pCallBack(hFile, lTotalSize,
strFileName , 1024*16);
      //1024*16 is default size
      if(!MDCS BLOB Get(hDBHandle, larrImageIDs[nCount],
eBlobType , &pCallBack))
         CloseHandle (hFile);
         return FALSE;
     cout << "Image " << nCount + 1 << " saved into " <<</pre>
strFileName.GetBuffer() <<</pre>
                 endl;
      CloseHandle (hFile);
    return TRUE;
}
```

The above examples demonstrate only some of the MDCStoreUtils API functions used to access the MDCStore database. To see how all the functions and classes discussed above are used in an application, you can examine the application MDCStoreImageExample attached with this documentation.

MDCStoreImageExample is a console application that demonstrates how to use functions in MDCStoreUtils API to get plate, site information, get images of from site, show database error messages, and use callback classes. It accepts the following parameters that are displayed if you call the application with an empty command line:

```
-u <User name>
-p <Password>
-dsn <ODBC datasource name>
-database <database name> - this parameter is optional
-plate <plate ID> - this parameter is optional
```

The application uses most of functions and classes discussed in this section. Always call **MDCStoreUtils_Init** to initialize the MDCStore interface before calling any other functions in MDCStoreUtils API and call **MDCStoreUtils_Finished** to detach the MDCStoreUtils interface when exiting the application.

The following is a list of the main source and header files used in the application:

- BlobCallbackEx.h
- GetImageRecord CallBack.h
- GetSitesByPlate_CallBack.h
- MDCStoreImageExample.h
- BlobCallbackEx.cpp
- GetSitesByPlate_CallBack.cpp
- GetImageRecord_CallBack.cpp
- MDCStoreImageExample.cpp main application

about Everyone 201 acquisition batch records 195 create profile 192 delete 191 instance record 194	attributes 230, 231 dataset 205, 206 delete 209 find 209 image 232 update 210 attributes information 22
profile records 193 profile table 193 update 173 active connections	B batch records 195
count 141 add column to table 146 new measurement 85 AddAssay 213 AddColumnToTable 146 AddMeasurement 85 admin group 200 all values 114	BeginTransaction 27 BLOB getting information about 51 save data 46 update description 48 BLOB data remove 47, 57
analysis 228 descriptions 208	C
annotation template 180 AppendMeasurementSet 66 apply layout 180 ApplyLayoutToPlate 180 assay create 213 description 215 assign	calculate Z prime 126, 127 CalculateStatisticEx 125 CalculateStatisticResults 124 callback 226 CallbackToAnalysisInfo 226 cancel current execution 30
assay attribute value 133, 134, 135 object ID 169 plate attribute 176, 177 AssignAttributeValueFloat 135, 177 AssignAttributeValueLong 134, 176 AssignAttributeValueString 133, 175 associate measurement set 101 AssociateWithPlate 101 Attach 58 attach BLOB 58 attribute information 174, 175	job 252 CancelJobProgress 252, 254 CancelQueryExecution 30 CanModify 187 CanModifyAssay 123 CanModifyFolder 90 CanWriteToLocation 59 CellOutlinesGetSiteCount 98 change password 201 ChangePassword 201 ChangeStatus 189

5000957 C 307

check	dataset 213
admin group 200	dataset analysis attributes 205
BLOB location 56	forgein key 146
BLOB saving 59	history record 147
can user modify dataset 214	image record 162
can user modify folder permissions	image source 161
90	location ID 45
connection 31	macros 261
current user assay modification	measurement set name 98
ability 123	new assay profile record 100
dataset 227	new folder 218
for subfolder 217 is measurement set is used in	new measurement 91 new measurement set 83
datasets 107	new measurement set measurement
object 144	85
plate modification 187	new measurement set run record 99
CheckIfDead 31	new record 241
choose	new series 161
user groups 199	new site 160
claim	outlines table 89
job 253	parameter 242
ClaimJob 253	plate 153
commit a transaction 28	plate attribute 174
	plate layout 179
CommitTransaction 28	stored procedure 147
Compact 140	Create(261
compact	CreateAttribute 131, 174
database 140	CreateConfig 233
configuration 235	CreateConnectionString 24
description 214	CreateFolder 91, 218
existing 235 name 234	CreateForeignKey 146
normalization 233, 234	CreateForPlate 119
ConvertQueryResultsToMacrosStructure	CreateHistoryRecord 147
261	CreateImageRecord 162
Copy 94, 211	CreateImageSource 161
	CreateJobRecord 241
copy dataset 211	CreateLayout 179
measurement set 94	CreateLocationOption 241
count	CreateMeasurement 85
acquired series 184	CreateNewName 98
acquired wells 184	CreateOutlinesTable 89
compounds 185	
controls 185, 186	CreatePlate 153
site appearances 95, 97	CreatePlate 153
count site appearances 98	CreateProfile 100, 192
CountActiveConnections 141	CreateRun 99
CountPlateDatasets 183	CreateSeries 161
Create 83, 192, 213	CreateSite 160
create 179, 192	CreateStoredProc 147
acquisition 192	CreateTable 145
acquisition profile 192	
assay 119, 213	
Assay attribute 131	
database connection 23	

	DeleteDataForPlate 88
D	DeleteFileLocation 56
D	DeleteFolder 89
data	DeleteImages 157
remove 143	DeleteLayout 163
database	DeleteLocationOption 246
available 142	DeleteMeasurement 97
	DeleteMeasurementInAssays 122
compact 140	DeleteProfile 105
get size 141	DeleteRecord 260
get version 142	
optimize 140	DeleteUnUsedLocations 55
Dataset 19	DeleteUserLocationOptions 246
dataset 207	description
assays 215	assays 215
attributes 224 attributres 206	descriptions
	analysis 208
copy 211	destroy
create 213 delete 228	database handle 26
filters 215	DestroyDBHandle 26
information 212	Disconnect 26
plates 223	disconnect from the database 26
record types 220	DoesNameExist 227
resluts 226	DoesObjectExist 144
results 225	DoesSubFolderExist 92, 217
scriptlet 221	DropTable 144
update 212	2.00.000
datasets 224	
default	
location 248	_
Delete 87, 191, 228	E
delete	
acquisition 191	error message
all measurement sets for a plate 88	default 42
analysis attributes 209	Everyone 201
assay attribute 136	Execute 145
dataset 228	execute
dataset folder 211	a statement 145
location 246	
macros 260	
measurement 97, 122	
measurement set 87	F
measurement set results for a plate	F
88	FillOutBlobLocationStruct 60
plate 156	find
plate images 157	record 257
plate layout 163	record by attributes 209
profile 105	FindAnalysisAttributesRecord 209
tree folder 89	FindMeasurementValues 129
user location 246	finds values 129
Delete 156	
DeleteAllForPlate 88	folder 215, 216 create 218
DeleteAnalysisAttributes 209	delete 211
DeleteAttribute 136	UCICLE ZII

measurement set sibling folders 84 measurements description 79 number of datasets 116 object ID 170 one batch record 195 plate attributes 158 ODBC connection 24 ogin structure 27 Get 46, 212 get 166, 167, 180 all available profiles 103 measurement set sibling folders 84 measurement set sibling f
generate one batch record 195 new database ID 48 plate attributes 158 ODBC connection 24 plate information 154, 156 ogin structure 27 profile record 104 Get 46, 212 profile table 193 get 166, 167, 180 record 106 all available profiles 103
new database ID 48 ODBC connection 24 ogin structure 27 Get 46, 212 get 166, 167, 180 all available profiles 103 plate attributes 158 plate information 154, 156 profile record 104 profile table 193 record 106 records 193
ODBC connection 24 plate information 154, 156 ogin structure 27 profile record 104 profile table 193 record 106 all available profiles 103 records 193
ogin structure 27 profile record 104 Get 46, 212 profile table 193 get 166, 167, 180 record 106 all available profiles 103 records 193
Get 46, 212 profile table 193 get 166, 167, 180 record 106 all available profiles 103 records 193
get 166, 167, 180 record 106 all available profiles 103 records 193
all available profiles 103 records 193
all available profiles 105
all data types 38, 39 result set data 34
all image information 54, 55 result set data types 37
all measurement sets 102 shape lines 78
all plate records 164 shape lines BLOB description 73, 74
all records 108 75, 76
all the shape descriptions 77 shape lines blob description 72, 116
all values 113, 114
delibate information 131, 132
attribute value 132, 133 spot 10 102 attributes 206 statistic results 124, 125
available databases 142 statistical values 109, 110
BLOB from storage 46 type of the database 28
BLOB information 51 unique annotation 130
BLOB location description 45 unique data types 38
data type record 81, 83 unique location 53
data types records 113 unique measurement set settings
database handle 25 115 database size 141 unique plate values 157
uatabase 5/26 1/1
444 445
default error message 42 Values 111, 112 deleted assays 118 GetAcqSeriesCount 184
description 120, 121 GetAcqSiteCount 183
description of a measurement set GetAcqWellCount 184
measurement 86, 87 GetAllBlobInfoByPlateID 55
description of all result 35 GetAllByAttributes 158, 231
description of individual result set 35 GetAllByGroup 259
folder description 92 GetAllDataTypes 38
folder path 93 GetAllDataTypesOfAssays 39
neader information 130 CetAllForAssay 224
CotAllEarDate 102 222
ID of the latest run 99 GetAllForPlate 102, 223 image information 44 GetAllInFolder 92, 215
information 29 GetAllMeasurements 220

GetAllMSetParamValues 224 GetAllMSetUniqueAnnotation 130 GetAllOrderedByAttributes 155 GetAllPropertyAttributes 154, 230

GetAllRecords 164 GetAllResult sets 35

GetAllSiteMeasurements 113 GetAnalysisAttributes 207

GetAnalysisAttributesByDataset 206

GetAnalysisCount 228 GetAnalysisDescription 214 GetAnalysisDescriptions 208 GetAssayAndFilters 215 GetAssayByRunID 106 GetAssayIDsOfPlate 115

GetAssayMeasurementRecord 79 GetAssayMeasurementsByPlateAndMea

surement 80

GetAssaySettingsOfPlate 115 GetAssaySiblingFolders 84 GetAttributeInfoByDBNam 175 GetAttributeInfoByDBName 132 GetAttributeInfoByDisplayName 131,

GetAttributes 22

GetAttributeValueByDBName 132 GetAttributeValueByDisplayName 133

GetAvailableDatabases 142

GetBatchRecord 195 GetBatchRecords 195

GetByID 100

GetCompleteImageInfo 186 GetCompoundCount 185

GetConfig 234 GetConfigAll 235 GetConfigByName 234 GetConfigForAssay 235 GetControlsCount 185 GetControlStatistic 186 GetDatabaseType 28 GetDBHandle 25

GetDBHandleFromString 24 GetDefaultErrorMsq 42

GetDetails 22

GetDetailsFromString 27 GetEveryoneGroupInfo 201

GetFolderPath 93

GetFreeSpaceFromFileLocation 242

GetHeaderAndFileInfo 136 GetHistoryRecord 148 GetImageIDs 188

GetImageObjectIDForImage 170

GetImageRecord 167

GetImageRecordPerPlate 168 GetImageSourceRecords 108 GetImageSourcesOfPlate 171

GetInfo 29, 44, 154

GetInfoBasedOnAssay 156 GetInfoBasedOnAssaySet 178 GetInfoByReferenceID 51 GetInfoByReferenceIDEx 51 GetInstanceRecord 194 GetItemFolder 216 GetJobQueue 250

GetJobQueueRecord 250 GetJobOueueRecords 249

GetLastError 41 GetLastErrorMsq 41 GetLatestAssayRunID 99

GetLocation 45 GetLocationID 45

GetLocationOptionByID 255 GetLocationOptionRecord 245

GetLocationOptions 244

GetLocationRecordsByLabel 245

GetMarkedAssaysWithCallback 118, 119

GetMaxTimePointForPlate 172 GetMeasurementAttributes 78 GetMeasurementByDBName 120

GetMeasurementByFunctionAndParame

terName 82

GetMeasurementByID 86

GetMeasurementByName 87, 120 GetMeasurementInfoByAssayAndColum

nName 81

GetMeasurementInfoByAssayAndPlate

113

GetMeasurementRecord 83 GetMeasurementsBySiteID 81 GetMeasurementStatistic 116

GetNumBlobOfObject 52 GetParameterByID 243 GetParameterByName 243 GetPlatesByDate 173 GetProfile 104, 193 GetProfileInfo 104 GetProfileRecords 193

GetNewDatabaseID 48

GetProfiles 103

GetRecord 106, 165, 259 GetRecordByGUID 258 GetRecordByID 257

GetResultInfo 226 GetResultsetData 34 GetResultsetInfo 35 GetResultsInfoByConfig 225 GetScopeAttributeByID 121 GetScopeAttributeByName 121 GetScriptletAssayIDByName 223 GetScriptletAssays 222 GetSeriesIDAtZAndT 172	GetZPrime 126 GetZPrimeScopeAttribute 127 groups 199 choose 199 Everyone 201 manage 198
GetSeriesRecord 170 GetShapeLineBlobDesc 72 GetShapeLineBlobDescByPlate 76 GetShapeLineBlobDescBySite 73 GetShapeLineBlobDescBySiteAndSeries 74 GetShapeLineBlobDescBySiteSeriesAnd Instance 75 GetShapeLineBlobPerAssay(77 GetShapeLinesBySite 116	HavePermissionsToModify 108, 214 history update status 148 history record create 147 get 148
GetShapeLinesBySiteAndSeries 117 GetSiblingFolders 219 GetSiteCount 95 GetSiteInfoImageByID 96 GetSiteInfoImageBySiteAndSeriesID 96 GetSiteLocationsForPlate 171 GetSiteRecord 166 GetSitesByPlate 166 GetSize 141 GetSpotID 102 GetStatisticalValuesForPlateAndAssay 109 GetStatisticalValuesForPlateAndMeasur ement 110 GetTemplate 180 GetTemplate 181 GetTemplates 181 GetThumbImageIDs 188 GetUniqueAttributeValues 157, 232 GetUniqueDataTypes 38 GetUniqueLocationNamebyFilter 53 GetUniqueMeasurementValues 128 GetUniqueMeasurementValues 128 GetUserGroups 199 GetUserGroups 199 GetUserInfo 200 GetUserLocationOption 248 GetValuesForPlateAndAssay 111 GetVersion 142	image attributes 232 create record 162 create source 161 get information 186 image ID 53 image information 96 image record 167 get 168 images share 160 import measurement set data 66, 67 plate layout 162 result set data 33 ImportDataToDB 33 ImportDataToDB 33 ImportMeasurementSet 67 information 201 about users 200 dataset 212 insert shapes 70 InsertAnalysisAttributes 205 InsertDValue 68 InsertNumericValue 68 InsertShapeLines 70 InsertStringValue 70 InsertSValue 69 Instance 19

IsUserInAdminGroup 200 MDCS_BLOB_GetAllBlobInfoByLocation FilterAndLocationInfo 54 Measurement Set 19 MeasurementGetSiteCount 97 Merge 95 J merge measurement sets 95 modify job 251, 252 folder 91, 219 cancel 252 measurement 123 claim 253 ModifyFolder 91, 219 create create job 241 refresh 254 reset 254 N new database handle 25 new series create 161 new site create 160 layout 179 normalization configuration 234 apply 180 location 247 default 248 delete 246 free disk space 242 options 240 location ID creation 45 object ID LocationIsUsed 56 assign 169 login information 22 get 170 Optimize 140 optimize database 140 optimize all measurement set tables 94 М OptimizeAll 94 macros 259, 260, 261 create 261 delete 260 information 259 P manage groups 198 parameter security 237 create 242 security access 210 **ID 243** manage security access 90 name 243 ManageFolderSecurity 90, 210, 237 update 244 ManageGroupsDlg 198 password change 201 ManageLocationOptions 240 Plate 19 ManageSecurity 158 plate ManageUsersDlg 197 annotation template 180 apply layout 163, 180 **BLOB** locations 59 assign attribute 176 MapLocations 59 attributes 154, 158

create 153 create attribute 174 create layout 179 dataset 223 date information 173 delete 156 delete images 157 delete layout 163 get all sites 166 get information 154, 156 get records 164 get single record 165 import layout 162 nformation 178 ordered by attributes 155 rename attribute 182 site positions 171 source records 171 timepoint 172 unique values 157 update information 177 profile 192 table 193	shape data 67 unused locations 55 RemoveBlobData 47, 57 RemoveFile 57 RemoveMarkedData 143 RemoveMarkedDataEX 143 RemoveShapeData 67 rename attribute 182 rename assay attribute 135 RenameAttribute 182 reset job 254 Reset Job 253 Reset job 253 restore assay 119 Result Set 19 result set data types 37 results dataset 225, 226 retrieve series ID 172 rollback a transaction 28 RollbackTransaction 28
quicklist 237	
	_ S
R	Save 46 save
Reconnect 26	BLOB data 46 SaveBlobEx 49
reconnect to the database 26	sciptlet 223
record 245, 258	scriptlet 221, 222
find 257 update 258	security folder 159
record types	manage 158, 237
dataset 220	manage accéss 210
records 249	Series 19
RefreshAllJobs 254	series ID
Reindex 93	retrieve 172 series record
reindex measurement set 93	get by ID 170
measurement set cell shapes 93	set
ReindexShapes 93	database connection mode 29
remove	SetAsyncMode 29
BLOB data 47, 57	setDataTypes 37
data 143	SetSilentMode 30
file 57	SetUserLocationOptions 247

share

file location 56

images 160 ShareSiteImagesByTimeAndZIndex 160 sibling folders 84 Site 19 site record get 166 spot ID 102 start a transaction 27 statement execution 145 statistic results 124, 125 statistical values 109 stored procedure create 147	plate layout template 182 profile 105 record 258 shape lines 71 shapes 70 UpdateAcquisition 173 UpdateAnalysisAttributes 210 UpdateBlobDescAndName 48 UpdateBlobDescAndNameEx 49 UpdateBlobInfo 58 UpdateConfig 233 UpdateDatabaseVersion 150 UpdateDataType 122 UpdateDataTypeByAssay 123
table add column 146 creation 145 DB_versions 150 drop 144 template 181 annotation 181 update 182 TestCredentials 23	UpdateHistoryRecordStatus 148 UpdateImageObjectID 169 UpdateInfo 177 UpdateJobProgress 251 UpdateJobQueueRecord 252 UpdateJobStatus 251 UpdateMeasurementData 129 UpdateMeasurementSetDescription 137 UpdateMeasurementSetName 137 UpdateParameter 244 UpdatePlateObjectImageID 53 UpdatePlateTemplate 182 UpdateProfile 105 UpdateRecord 258, 260
unique values plate 157 unused locations 55	UpdateShapeLines 71 UsedInDatasets 107 user information 200
Update 212 update 53 acquisition 173 attributes 210 BLOB description 48, 49, 58 dataset 212 DB_versions table 150 find finished 149	W well count acquired 184
folder 227 history record status 148 job status 251 measurement 122 measurement data 130 measurement set description 137	Z Z prime 127

5000957 C 315

measurement set name 137 normalization configuration 233